

QUANTITATIVE DETERMINATION OF HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA “NATURAL ANTIBODIES”

NEW PRODUCT

Human Anti-Alpha-Galactosyl IgG ELISA
Human Anti-Alpha-Galactosyl IgM ELISA
Human Anti-Alpha-Galactosyl IgA ELISA

- › Excellent analytical characteristics
- › Validated for human serum and plasma (EDTA, citrate, heparin) samples
- › Preliminary population data

IMMUNOLOGY · IMMUNE RESPONSE
AUTOIMMUNITY · TRANSPLANTATION
ONCOLOGY

HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

Introduction

Natural antibodies circulate in serum of healthy humans naturally, without previous immunization [1]. They can be directed against the individual's own antigens as well as against foreign antigens [2, 3, 4]. The repertoire of natural antibodies is encoded by germ-line genes [5, 4, 6] and is not subjected to recombination. Natural antibodies are polyreactive [2, 3, 6] and react with low affinity but high avidity [7].

Among the natural antibodies directed against saccharide antigens, the best known are:

- 1) isohemagglutinins (anti-A, anti-B) which prevent the transfusion of allogeneic blood [8] and allotransplantation [9];
- 2) xenoagglutinins (80% of the repertoire of human xenoreactive natural antibodies is specific for the terminal α -galactose determinant) [10] which have traditionally been detected by agglutination test with rabbit or porcine erythrocytes [11]. The epitope against which they are directed was identified as Gal α 1-3Gal β 1-4GlcNAc-R [12]. This epitope is present in cell membranes of mammals with the exception of humans and Old World primates [13] in which the gene encoding for α -1,3-galactosyltransferase was inactivated during evolution [14]. This epitope is also present in lipopolysaccharides and capsular polysaccharides of various

gram-negative bacteria colonizing the gastrointestinal tract of humans [15].

In human serum, natural anti- α -galactosyl antibodies (anti- α -Gal) constitute approximately 4 - 8 % of total serum IgM, approx. 1 - 2 % of total serum IgG [16], and they can also be found as the IgA isotope in serum, saliva, milk, colostrum and vaginal washings [17, 15]. About 1 % of circulating B lymphocytes in adults is capable of producing these antibodies [18]. Their production is thought to be induced postnatally after antigenic stimuli by intestinal microflora [19].

The biological roles of natural antibodies are thought to be protection against infectious agents [20], particularly at the mucosal surface [15], effector functions of the immune system, e.g. in removing senescent or otherwise altered cells and physiologically degraded molecules [21], and regulation of immune responses, particularly autoimmune responses [6]. Moreover, the presence of anti- α -galactosyl antibodies in humans represents an important obstacle to xenotransplantations of porcine grafts since these antibodies participate in their rejection [22]. Natural antibodies also play an important role in malignant processes, and they have potential to be used in treatment of tumors that express these (α -Gal) antigens (e.g. breast cancer) [23, 24, 25, 26].

BioVendor Human Anti-Alpha-Galactosyl IgG ELISA (RD199178100R)
BioVendor Human Anti-Alpha-Galactosyl IgM ELISA (RD199178110R)
BioVendor Human Anti-Alpha-Galactosyl IgA ELISA (RD199178120R)

Intended use

The Human Anti-Alpha-Galactosyl IgG (or IgM or IgA, respectively) ELISA is an enzyme immunoassay for the quantitative measurement of human anti-alpha-galactosyl antibodies in the IgG (or IgM or IgA, respectively) class.

- The total assay time is less than 3 hours
- The kit measures IgG (or IgM or IgA, respectively) anti-alpha-galactosyl antibodies in human serum and plasma (EDTA, citrate, heparin)
- Assay format is 96 wells
- Calibrator is human serum based
- Components of the kit are provided ready to use, concentrated or lyophilized

Clinical application

- Immunology, Immune Response
- Autoimmunity
- Infection and Inflammation
- Transplantation
- Oncology


QUANTITATIVE DETERMINATION OF HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

Test principle

In the BioVendor Human Anti-Alpha-Galactosyl IgG (or IgM or IgA, respectively) ELISA, Calibrators and samples are incubated in microtitration wells pre-coated with chemically synthesized molecule displaying terminal beta-disaccharide Gal α 1-3Gal which is recognized by human anti-alpha-galactosyl antibodies. After 60 minutes incubation followed by washing, polyclonal antibody against human IgG (or IgM or IgA, respectively) conjugated with horseradish peroxidase (HRP) is added into the wells and incubated with the captured anti-alpha-galactosyl IgG (or IgM or IgA, respectively) for 60

minutes. The last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration of IgG (or IgM or IgA, respectively) anti-alpha-galactosyl antibodies. A calibration curve is constructed by plotting absorbance values against anti-alpha-galactosyl IgG (or IgM or IgA, respectively) concentrations of calibrators and concentrations of unknown samples are determined using this calibration curve.

HUMAN ANTI-ALPHA-GALACTOSYL IgG ELISA CAT. NO.: RD199178100R	
Standards	3.13 to 100 U/ml
Limit of detection	1.44 U/ml
HUMAN ANTI-ALPHA-GALACTOSYL IgM ELISA CAT. NO.: RD199178110R	
Standards	3.13 to 100 U/ml
Limit of detection	1.35 U/ml
HUMAN ANTI-ALPHA-GALACTOSYL IgA ELISA CAT. NO.: RD199178120R	
Standards	3.13 to 100 U/ml
Limit of detection	0.61 U/ml



Summary of protocol

- Reconstitute Master Calibrator and prepare set of Calibrators
- Dilute samples (100x)
- Add 100 μ l Calibrators and samples
- Incubate at RT for 1 hour
- Wash plate 3 times
- Add 100 μ l Conjugate Solution
- Incubate at RT for 1 hour
- Wash plate 3 times
- Add 100 μ l Substrate Solution
- Incubate at RT for 30 min
- Add 100 μ l Stop Solution
- Read absorbance and calculate results

HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

Precision

Intra-assay (Within-Run) (n=8)

Sample	Coefficient of Variation (%)		
	Human Anti- α -Gal IgG ELISA	Human Anti- α -Gal IgM ELISA	Human Anti- α -Gal IgA ELISA
1	3.3	2.7	1.7
2	3.7	3.1	6.8

Inter-assay (Run-to-Run) (n=5)

Sample	Coefficient of Variation (%)		
	Human Anti- α -Gal IgG ELISA	Human Anti- α -Gal IgM ELISA	Human Anti- α -Gal IgA ELISA
1	5.1	8.4	5.7
2	7.9	4.4	8.0

Spiking recovery

Serum samples were spiked with different amounts of human anti-alpha-galactosyl IgG (or IgM or IgA, respectively) antibody and assayed. Results are expressed as observed/expected ratio (O/E) range.

ELISA	O/E Range (%)
Human Anti- α -Gal IgG ELISA	80.8 - 97.9
Human Anti- α -Gal IgM ELISA	88.4 - 93.6
Human Anti- α -Gal IgA ELISA	83.2 - 104.0


Linearity

Serum samples were serially diluted with Dilution Buffer and assayed. Results are expressed as observed/expected ratio (O/E) range.

ELISA	O/E Range (%)
Human Anti- α -Gal IgG ELISA	83.9 - 96.2
Human Anti- α -Gal IgM ELISA	94.2 - 102.9
Human Anti- α -Gal IgA ELISA	91.5 - 105.7

Effect of sample matrix

EDTA, citrate and heparin plasma samples were compared to respective serum samples from the same 10 individuals. Results obtained for Human Anti-Alpha-Galactosyl IgG ELISA are shown in the graph below for illustration:



The values of coefficient of determination are summarized below:

ELISA	Coefficient of Determination R ²			
	Serum	EDTA plasma	Citrate plasma	Heparin plasma
Human Anti- α -Gal IgG ELISA	-	0.97	0.97	0.98
Human Anti- α -Gal IgM ELISA	-	1.00	0.97	0.99
Human Anti- α -Gal IgA ELISA	-	0.98	0.97	0.99

QUANTITATIVE DETERMINATION OF HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

BioVendor Human Anti-Alpha-Galactosyl IgG ELISA (RD199178100R)

Cross-reactivity


Mammalian serum Sample	Bovine	Cat	Dog	Goat	Hamster	Horse	Monkey	Mouse	Pig	Rabbit	Rat	Sheep
Observed cross-reactivity	no	no	no	no	no	no	yes	no	no	no	no	no

Preliminary Population Data

The following results were obtained when serum samples from 168 unselected donors (91 men + 77 women) 20 - 65 years old were assayed with the BioVendor Human Anti-Alpha-Galactosyl IgG ELISA in our laboratory.

Age and Sex Dependent Distribution of Hu Anti-Alpha-Galactosyl IgG

Sex	Age (years)	n	Mean Anti- α -Gal IgG (U/ml)	Median Anti- α -Gal IgG (U/ml)	SD Anti- α -Gal IgG (U/ml)	Min. Anti- α -Gal IgG (U/ml)	Max. Anti- α -Gal IgG (U/ml)
Male	20-29	18	41.98	34.35	35.13	5.22	126.94
	30-39	30	55.49	46.03	41.32	5.88	148.86
	40-49	33	65.22	32.48	76.40	4.57	305.41
	50-65	10	123.43	80.26	144.55	9.85	534.55
Female	22-29	14	48.45	26.81	42.53	3.06	138.40
	30-39	29	59.85	33.38	69.99	4.00	351.77
	40-49	27	38.75	23.14	42.91	5.30	183.06
	50-61	7	43.30	34.16	25.53	17.98	86.21



A very high degree of variability in concentrations was observed, indicating that the inherent variability of anti-alpha-galactosyl IgG levels appears to be very high in this population.

An appropriate dilution should be assessed by the researcher (due to the high variability of serum anti-alpha-galactosyl IgG level between individuals) in advance to batch measurement. Recommended starting dilution is 100x.

HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

BioVendor Human Anti-Alpha-Galactosyl IgM ELISA (RD199178110R)

Cross-reactivity


Mammalian serum Sample	Bovine	Cat	Dog	Goat	Hamster	Horse	Monkey	Mouse	Pig	Rabbit	Rat	Sheep
Observed cross-reactivity	no	no	no	no	no	no	yes	no	no	no	no	no

Preliminary Population Data

The following results were obtained when serum samples from 168 unselected donors (91 men + 77 women) 20 - 65 years old were assayed with the BioVendor Human Anti-Alpha-Galactosyl IgM ELISA in our laboratory.

Age and Sex Dependent Distribution of Hu Anti-Alpha-Galactosyl IgM

Sex	Age (years)	n	Mean Anti- α -Gal IgM (U/ml)	Median Anti- α -Gal IgM (U/ml)	SD Anti- α -Gal IgM (U/ml)	Min. Anti- α -Gal IgM (U/ml)	Max. Anti- α -Gal IgM (U/ml)
Male	20-29	18	26.54	21.34	19.27	3.60	65.80
	30-39	30	32.05	24.91	25.03	7.16	123.85
	40-49	33	33.04	24.38	26.85	1.61	113.54
	50-65	10	24.64	20.75	10.02	10.45	42.46
Female	22-29	14	45.14	45.15	21.67	16.00	83.75
	30-39	29	52.02	29.35	92.51	7.70	520.32
	40-49	27	44.19	19.82	51.02	5.53	222.03
	50-61	7	39.11	29.95	32.55	3.86	91.16



A very high degree of variability in concentrations was observed, indicating that the inherent variability of anti-alpha-galactosyl IgM levels appears to be very high in this population.

An appropriate dilution should be assessed by the researcher (due to the high variability of serum anti-alpha-galactosyl IgM level between individuals) in advance to batch measurement. Recommended starting dilution is 100x.

QUANTITATIVE DETERMINATION OF HUMAN ANTI-ALPHA-GALACTOSYL IgG, IgM, IgA

BioVendor Human Anti-Alpha-Galactosyl IgA ELISA (RD199178120R)

Cross-reactivity


Mammalian serum Sample	Bovine	Cat	Dog	Goat	Hamster	Horse	Monkey	Mouse	Pig	Rabbit	Rat	Sheep
Observed cross-reactivity	no	yes	yes	no	no	yes	yes	no	no	no	no	no

Preliminary Population Data

The following results were obtained when serum samples from 168 unselected donors (91 men + 77 women) 20 - 65 years old were assayed with the BioVendor Human Anti-Alpha-Galactosyl IgG ELISA in our laboratory.

Age and Sex Dependent Distribution of Hu Anti-Alpha-Galactosyl IgA

Sex	Age (years)	n	Mean Anti- α -Gal IgA (U/ml)	Median Anti- α -Gal IgA (U/ml)	SD Anti- α -Gal IgA (U/ml)	Min. Anti- α -Gal IgA (U/ml)	Max. Anti- α -Gal IgA (U/ml)
Male	20-29	18	20.00	15.01	13.17	2.09	54.69
	30-39	30	33.80	22.39	33.42	0.00	138.59
	40-49	33	63.79	39.78	76.15	6.16	423.00
	50-65	10	40.11	29.06	34.24	2.27	121.73
Female	22-29	14	18.90	12.09	19.02	0.16	67.90
	30-39	29	29.16	19.36	30.64	2.79	148.40
	40-49	27	28.48	16.48	31.17	4.27	151.37
	50-61	7	60.00	33.91	60.94	12.22	195.19



A very high degree of variability in concentrations was observed, indicating that the inherent variability of anti-alpha-galactosyl IgA levels appears to be very high in this population.

An appropriate dilution should be assessed by the researcher (due to the high variability of serum anti-alpha-galactosyl IgA level between individuals) in advance to batch measurement. Recommended starting dilution is 100x.

References

1. Coutinho A, Kazatchkine MD, Avrameas S. Natural autoantibodies. *Curr Opin Immunol*; 7(6):812-8. (1995)
2. Avrameas S, Dighiero G, Lymberi P, Guilbert B. Studies on natural antibodies and autoantibodies. *Ann Immunol (Paris)*; 134D(1):103-13. (1983)
3. Avrameas S, Guilbert B, Mahana W, Matsiota P, Ternynck T. Recognition of self and non-self constituents by polyspecific autoreceptors. *Int Rev Immunol*; 3(1-2):1-15. (1988)
4. Tomer Y, Shoenfeld Y. The significance of natural autoantibodies. *Immunol Invest*; 17(5):389-424. (1988)
5. Dighiero G, Lymberi P, Guilbert B, Ternynck T, Avrameas S. Natural autoantibodies constitute a substantial part of normal circulating immunoglobulins. *Ann NY Acad Sci*; 475:135-45. (1986)
6. Avrameas S. Natural autoantibodies: from 'horror autotoxicus' to 'gnothi seauton'. *Immunol Today*; 12(5):154-9. (1991)
7. Ternynck T, Avrameas S. Murine natural monoclonal autoantibodies: a study of their polyspecificities and their affinities. *Immunol Rev*; 94:99-112. (1986)
8. Landsteiner K. *The Specificity of Serological Reactions*. Harvard University Press, Cambridge, 1945.
9. Starzl TE, Marchioro TL, Holmes JH, Hermann G, Brittain RS, Stonington OH, Talmage DW, Waddell WR. Renal homografts in patients with major donor-recipient blood group incompatibilities. *Surgery*; 55:195-200. (1964)
10. Parker W, Bruno D, Holzknacht ZE, Platt JL. Characterization and affinity isolation of xenoreactive human natural antibodies. *J Immunol*; 153(8):3791-803. (1994)
11. Tönder O, Natvig JB, Matre R. Antibodies in human sera to rabbit erythrocytes. *Immunology*; 12(6):629-37. (1967)
12. Galili U, Macher BA, Buehler J, Shohet SB. Human natural anti-alpha-galactosyl IgG. II. The specific recognition of alpha (1----3)-linked galactose residues. *J Exp Med*; 162(2):573-82. (1985)
13. Galili U, Clark MR, Shohet SB, Buehler J, Macher BA. Evolutionary relationship between the natural anti-Gal antibody and the Gal alpha 1----3Gal epitope in primates. *Proc Natl Acad Sci U S A*; 84(5):1369-73. (1987)
14. Galili U, Swanson K. Gene sequences suggest inactivation of alpha-1,3-galactosyltransferase in catarrhines after the divergence of apes from monkeys. *Proc Natl Acad Sci U S A*; 88(16):7401-4. (1991)
15. Hamadeh RM, Galili U, Zhou P, Griffiss JM. Anti-alpha-galactosyl immunoglobulin A (IgA), IgG, and IgM in human secretions. *Clin Diagn Lab Immunol*; 2(2):125-31. (1995)
16. McMorrow IM, Comrack CA, Sachs DH, DerSimonian H. Heterogeneity of human anti-pig natural antibodies cross-reactive with the Gal(alpha1,3)Galactose epitope. *Transplantation*; 64(3):501-10. (1997)
17. Galili U. Evolution and pathophysiology of the human natural anti-alpha-galactosyl IgG (anti-Gal) antibody. *Springer Semin Immunopathol*; 15(2-3):155-71. (1993)
18. Galili U, Anaraki F, Thall A, Hill-Black C, Radic M. One percent of human circulating B lymphocytes are capable of producing the natural anti-Gal antibody. *Blood*; 82(8):2485-93. (1993)
19. Galili U, Mandrell RE, Hamadeh RM, Shohet SB, Griffiss JM. Interaction between human natural anti-alpha-galactosyl immunoglobulin G and bacteria of the human flora. *Infect Immun*; 56(7):1730-7. (1988)
20. Ochsenbein AF, Zinkernagel RM. Natural antibodies and complement link innate and acquired immunity. *Immunol Today*; 21(12):624-30. (2000)
21. Grabar P. Autoantibodies and the physiological role of immunoglobulins. *Immunol Today*; 4(12):337-340. (1983)
22. Good AH, Cooper DK, Malcolm AJ, Ippolito RM, Koren E, Neethling FA, Ye Y, Zuhdi N, Lamontagne LR. Identification of carbohydrate structures that bind human antiporcine antibodies: implications for discordant xenografting in humans. *Transplant Proc*; 24(2):559-62. (1992)
23. Castronovo V, Colin C, Parent B, Foidart JM, Lambotte R, Mahieu P. Possible role of human natural anti-Gal antibodies in the natural antitumor defense system. *J Natl Cancer Inst*; 81(3):212-6. (1989)
24. Galili U, LaTemple DC. Natural anti-Gal antibody as a universal augmentor of autologous tumor vaccine immunogenicity. *Immunol Today*; 18(6):281-5. (1997)
25. Galili U, Wigglesworth K, Abdel-Motal UM. Intratumoral injection of alpha-gal glycolipids induces xenograft-like destruction and conversion of lesions into endogenous vaccines. *J Immunol*; 178(7):4676-87. (2007)
26. Manches O, Plumas J, Lui G, Chaperot L, Molens JP, Sotto JJ, Bensa JC, Galili U. Anti-Gal-mediated targeting of human B lymphoma cells to antigen-presenting cells: a potential method for immunotherapy using autologous tumor cells. *Haematologica*; 90(5):625-34. (2005)

Contact Information



BioVendor – Laboratorni medicina a.s.
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185, Fax: +420 549 211 460
E-mail: info@biovendor.com

BioVendor GmbH
Otto-Hahn-Straße 16, 34123 Kassel, Germany
Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
Phone: +43 1 890 9025, Fax: +43 1 890 5163
E-mail: infoAustria@biovendor.com

BioVendor, LLC
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
Phone: +1-800-404-7807, Phone: +1-828-575-9250
Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com