

QUANTITATIVE DETERMINATION OF HUMAN GDF-15/MIC-1


IVD CE

Human GDF-15/MIC-1 ELISA

- > High sensitivity (22 pg/ml)
- > Very good analytical characteristics
- > Validated for human serum, plasma (EDTA, citrate, heparin)

CARDIOLOGY
PREGNANCY
ONCOLOGY
HEMATOLOGY

HUMAN GDF-15/MIC-1 ELISA



Introduction

Growth differentiation factor 15 (GDF-15) is a member of the transforming growth factor b (TGF- b) cytokine superfamily and was originally cloned as macrophage-inhibitory cytokine 1 (MIC-1).

Circulating GDF-15 concentrations are increased across a wide spectrum of cardiovascular diseases, including acute and chronic coronary artery disease, congestive heart failure, and ischemic stroke. GDF-15 is also upregulated by other cardiovascular events triggering oxidative stress, including pressure overload, and atherosclerosis.

Moreover, increased circulating GDF-15 concentrations have been linked to an enhanced risk of future adverse cardiovascular events in elderly women and it was describe as a new biomarker of the risk of death in patients with non-ST-elevation acute coronary syndrome.

Serum GDF-15 concentrations increase in maternal serum with advancing gestation in normal pregnancy. Low GDF-15 concentrations reportedly are associated with an increased risk of preterm labor or miscarriage.

Increased GDF-15 expression has been documented in a variety of epithelial cell lines, including breast, pancreas, colorectal, and prostate cancers. Microarray studies have revealed increased expression of GDF-15 in patients with

breast cancer, and serum GDF-15 levels are the best single predictor of the presence of pancreatic carcinoma. In the case of prostate cancer, serum GDF-15 levels increase with progression of disease to metastasis. In colon cancer, increasing GDF-15 expression is associated with the progression of colonic adenomas to invasive cancer and subsequent metastasis, with serum levels at presentation being an independent predictor of subsequent disease-free status and overall survival.

GDF-15 levels in blood plasma have been found to be dramatically elevated in beta-thalassemia patients compared to healthy donors and patients with hereditary hemochromatosis, another form of iron overload disease. In addition, the role of GDF-15 in other disorders characterized by ineffective erythropoiesis, as well as the role of GDF-15 in regulation of iron metabolism is under investigation. There are some hypotheses for treatment of thalassemia by administration of GDF-15 antagonist, and to reduce hepcidin levels by administration of GDF-15, a GDF-15 substitute, or GDF-15 agonist.

QUANTITATIVE DETERMINATION OF HUMAN GDF-15/MIC-1

BioVendor Human GDF-15/MIC-1 ELISA (RD191135200R)

Intended use

The RD191135200R Human GDF-15 ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human GDF-15/MIC-1 (growth differentiation factor 15 / macrophage-inhibitory cytokine 1).

- The total assay time is less than 3.5 hours
- The kit measures human GDF-15 in serum and plasma (EDTA, citrate, heparin)
- Assay format is 96 wells
- Quality Controls are human serum based
- Standard is recombinant protein based
- Components of the kit are provided ready to use, concentrated or lyophilized

Clinical application


- Cardiology
- Pregnancy
- Oncology
- Hematology

HUMAN GDF-15/MIC-1 ELISA CAT. NO.: RD191135200R

Assay format	Sandwich ELISA, Biotin-labelled antibody, 96 wells/kit
Samples	Plasma, Serum
Standards	70 to 4480 pg/ml
Limit of detection	22 pg/ml

Test principle

In the BioVendor Human GDF-15/MIC-1 ELISA, the standards, quality controls and samples are incubated in microtiterate wells pre-coated with polyclonal anti-human GDF-15 antibody. After 60 min incubation and a washing, biotin labelled polyclonal anti-human GDF-15 antibody is added and incubated with captured GDF-15 for 60 min. After another washing, the streptavidin-HRP conjugate is added. After 30 min incubation and the last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration of GDF-15. A standard curve is constructed by plotting absorbance values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.



HUMAN GDF-15/MIC-1 ELISA

Precision

Intra-assay (Within-Run) (n=8)

Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
1	7279.7	456.9	6.3
2	870.9	62.9	7.2

Inter-assay (Run-to-Run) (n=8)

Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
1	1771.0	150.6	8.5
2	397.8	37.6	9.5

Spiking recovery

Serum samples were spiked with different amounts of human GDF-15 and assayed.

Sample	Observed (pg/ml)	Expected (pg/ml)	Recovery O/E (%)
1	1650.2	-	-
	4520.9	4450.2	101.6
	3070.6	3050.2	100.7
	2136.9	2350.2	90.9
2	1051.95	-	-
	3666.6	3851.9	95.2
	2316.7	2451.9	94.5
	1661.6	1751.9	94.8


Linearity

Serum samples were serially diluted with Dilution Buffer and assayed.

Sample	Dilution	Observed (pg/ml)	Expected (pg/ml)	Recovery O/E (%)
1	-	8719.2	-	-
	2x	4053.5	4359.6	93.0
	4x	1949.1	2179.8	89.4
	8x	873.1	1089.9	80.1
2	-	7694.1	-	-
	2x	3888.9	3847.0	101.1
	4x	1897.7	1923.5	98.7
	8x	943.0	961.8	98.0

Effect of sample matrix

Heparin, citrate and EDTA plasmas were compared to respective serum samples from the same 10 individuals. Results are shown below:



Summary of protocol

- Reconstitute QCs and Master Standard and prepare set of Standards
- Dilute samples 5x
- Add 100 µl Standards, QCs and samples
- Incubate at RT for 1 hour/300 rpm
- Prepare Wash Solution
- Wash plate 3 times
- Prepare Biotin Labelled Antibody Solution
- Add 100 µl Biotin Labelled Antibody Solution
- Incubate at RT for 1 hour/300 rpm
- Wash plate 3 times
- Add Streptavidin-HRP Conjugate
- Incubate at RT for 30 min
- Wash plate 3 times
- Add 100 µl Substrate Solution
- Incubate at RT for 30 min
- Add 100 µl Stop Solution
- Read absorbance and calculate results


QUANTITATIVE DETERMINATION OF HUMAN GDF-15/MIC-1

Preliminary Population Data

The following results were obtained when serum samples from 142 unselected donors (87 women + 55 men) 6-86 years old were assayed with the BioVendor Human GDF-15/MIC-1 ELISA in our laboratory:

Age and Sex Dependent Distribution of GDF-15/MIC1

Sex	Age (years)	n	Mean GDF-15 (pg/ml)	Median GDF-15 (pg/ml)	SD GDF-15 (pg/ml)	Min. GDF-15 (pg/ml)	Max. GDF-15 (pg/ml)
Men	14 - 19	2	573.5	573.5	20.3	553.2	593.7
	21 - 49	20	840.8	814.9	272.9	386.6	1344.3
	50 - 85	33	2972.5	2162.0	2437.0	758.0	11340.0
Women	6 - 18	6	852.1	743.3	251.0	581.3	1204.7
	20 - 49	32	980.0	798.9	586.2	402.8	3307.8
	50 - 86	49	2584.4	1986.5	2067.5	658.5	13652.8



Related products

· RBG10164005 GDF-15/MIC-1 Human Cell Culture

References

1. Arslan D, Cihan T, Kose D, Vatanshev H, Cimen D, Koksaly Y, Oran B, Akyurek F: Growth-differentiation factor-15 and tissue doppler imaging in detection of asymptomatic anthracycline cardiomyopathy in childhood cancer survivors. Clin Biochem. 2013;46(13-14):1239-43.
2. Eggers KM, Kempf T, Wallentin L, Wollert KC, Lind L: Change in growth differentiation factor 15 concentrations over time independently predicts mortality in community-dwelling elderly individuals. Clin Chem 2013;59(7):1091-8
3. Kempf T, Zarbock A, Widera C, Butz S, Stadtmann A, Rossaint J, Bolomini-Vittoni M, Korf-Klingebiel M, Napp LC, Hansen B, Kanwischer A, Bavendiek U, Beutel G, Hapke M, Sauer MG, Laudanna C, Hogg N, Vestweber D, Wollert KC: GDF-15 is an inhibitor of leukocyte integrin activation required for survival after myocardial infarction in mice. Nat Med. 2011;17(5):581-8
4. Eggers KM, Kempf T, Allhoff T, Lindahl B, Wallentin L, Wollert KC: Growth -differentiation factor-15 for early risk stratification in patients with acute chest pain. Eur Heart J. 2008;29(19):2327-35
5. Brown DA, Moore J, Johnen H, Smeets TJ, Bauskin AR, Kuffner T, Weedon H, Milliken ST, Tak PP, Smith MD, Breit SN.: Serum macrophage inhibitory cytokine 1 in rheumatoid arthritis: a potent marker of erosive joint destruction. Arthritis Rheum. 2007; 56(3):753-64
6. Huang CY, Beer TM, Higano CS, True LD, Vessella R, Lange PH, Garzotto M, Nelson PS.: Molecular alternations in prostate carcinomas that associate with in vivo exposure to chemtherapy: identification of a cytoprotective mechanism involving growth differentiation factor-15. Clin Cancer Res. 2007;13(19):5825-33
7. Kempf T, Horn-Wichmann R, Brabant G, Peter T, Allhoff T, Klein G, Drexler H, Johnston N, Wallentin L, Wollert KC.: Circulating concentrations of growth-differentiation factor 15 in apparently healthy elderly individuals and patients with chronic heart failure as assessed by a new immunoradiometric sandwich assay. Clin Chem. 2007;53(2):284-91
8. Kempf T, Björklund E, Olofsson S, Lindahl B, Allhoff T, Peter T, Tongers J, Wollert KC, Wallentin L.: Growth-differentiation factor-15 improves risk stratification in ST-segment elevation myocardial infarction. Eur Heart J. 2007;28(23):2858-65
9. Selander KS, Brown DA, Sequeiros GB, Hunter M, Desmond R, Parpala T, Risteli J, Breit SN, Jukkola-Vuorinen A.: Serum macrophage inhibitory cytokine-1 concentrations correlate with the presence of prostate cancer bone metastases. Cancer Epidemiol Biomarkers Prev. 2007;16(3):532-7
10. Tanno T, Bhanu NV, Oneal PA, Goh SH, Staker P, Lee YI, Moroney JW, Reed CH, Luban NL, Wang RH, Eling TE, Childs R, Ganz T, Leitman SF, Fucharoen S, Miller JL.: High levels of GDF15 in thalassemia suppress expression of the iron regulatory protein hepcidin. Nat Med. 2007;13(9):1096-101
11. Wollert KC, Kempf T, Lagerqvist B, Lindahl B, Olofsson S, Allhoff T, Peter T, Siegbahn A, Venge P, Drexler H, Wallentin L.: Growth differentiation factor 15 for risk stratification and selection of an invasive treatment strategy in non ST-elevation acute coronary syndrome. Circulation. 2007;116(14):1540-8
12. Ago T, Sadoshima J.: GDF15, a Cardioprotective TGF- β Superfamily Protein. Circ Res. 2006;98(3):294-7
13. Kempf T, Eden M, Strelau J, Naguib M, Willenbockel C, Tongers J, Heineke J, Kotlarz D, Xu J, Molkenkin JD, Niessen HW, Drexler H, Wollert KC.: The Transforming Growth Factor- β Superfamily Member Growth-Differentiation Factor-15 Protects the Heart From Ischemia/Reperfusion Injury. Circ Res. 2006;98(3):294-7
14. Brown DA, Ward RL, Buckhaults P, Liu T, Romans KE, Hawkins NJ, Bauskin AR, Kinzler KW, Vogelstein B, Breit SN.: MIC-1 serum level and genotype: associations with Progress and Prognosis of Colorectal Carcinoma. Clin Cancer Res. 2003;9(7):2642-50
15. Keelan JA, Wang K, Chaiworapongsa T, Romero R, Mitchell MD, Sato TA, Brown DA, Fairlie WD, Breit SN.: Macrophage inhibitory cytokine 1 in fetal membranes and amniotic fluid from pregnancies with and without preterm labour and premature rupture of membranes. Mol Hum Reprod. 2003;9(9):535-40
16. Liu T, Bauskin AR, Zaunders J, Brown DA, Pankhurst S, Russell PJ, Breit SN.: Macrophage inhibitory cytokine 1 reduce cell adhesion and induces apoptosis in prostate cancer cells. Cancer Res. 2003;63(16):5034-40
17. Moore AG, Brown DA, Fairlie WD, Bauskin AR, Brown PK, Munier ML, Russell PK, Salamonsen LA, Wallace EM, Breit SN.: The transforming growth factor- β Superfamily Cytokine-1 Is Present in High Concentrations in the Serum of Pregnant Women. J Clin Endocrinol Metab. 2000;85(12):4781-8
18. Bootcov MR, Bauskin AR, Valenzuela SM, Moore AG, Bansal M, He XY, Zhang HP, Donnellan M, Mahler S, Pryor K, Walsh BJ, Nicholson RC, Fairlie WD, Por SB, Robbins JM, Breit SN.: MIC-1, a novel macrophage inhibitory cytokine, is a divergent member of the TGF-beta superfamily. Proc Natl Acad Sci U S A. 1997;94(21):11514-9

References to this product:

1. Stejskal D, Karpisek M, Humenanska V, Lacnak B, Svstak M. Macrophage-inhibitory cytokine-1 (mic-1) in differential diagnosis of dyspnea-A pilot study. Clin Biochem. 2009 Mar 31

Contact Information



BioVendor – Laboratorni medicina a.s.
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185, Fax: +420 549 211 460
E-mail: info@biovendor.com

BioVendor GmbH
Otto-Hahn-Straße 16, 34123 Kassel, Germany
Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
E-mail: infoEU@biovendor.com

➤ www.biovendor.com

BioVendor GesmbH
Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
Phone: +43 1 890 9025, Fax: +43 1 890 5163
E-mail: infoAustria@biovendor.com

BioVendor, LLC
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
Phone: +1-800-404-7807, Phone: +1-828-575-9250
Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com