

QUANTITATIVE DETERMINATION OF HUMAN ADIPOCYTE FABP

IVD CE


HUMAN ADIPOCYTE FABP (FABP4) ELISA

- › High sensitivity (0.05 ng/ml)
- › Excellent analytical characteristics
- › Validated for human plasma (Citrate, EDTA, Heparin), serum



**ENERGY METABOLISM
AND BODY WEIGHT REGULATION
DIABETOLOGY**

HUMAN ADIPOCYTE FABP



Introduction

Adipocyte fatty acid binding protein AFABP is a 15 kDa member of the intracellular fatty acid binding protein (FABP) family, which is known for the ability to bind fatty acids and related compounds (bile acids or retinoids) in an internal cavity. AFABP is expressed in a differentiation-dependent fashion in adipocytes and is a critical gene in the regulation of the biological function of these cells. In mice, targeted mutations in FABP4 (mouse gene is also called aP2 and its relevant protein P2 adipocyte protein or 3T3-L1 lipid binding protein) provide significant protection from hyperinsulinemia and insulin resistance in the context of both dietary and genetic obesity. Adipocytes obtained from AFABP-deficient mice also have reduced efficiency of lipolysis in vitro and in vivo, and these mice exhibited moderately improved systemic dyslipidemia. Recent studies also demonstrated AFABP

expression in human macrophages upon differentiation and activation. In these cells, AFABP modulates inflammatory responses and cholesterol ester accumulation, and total or macrophage-specific AFABP deficiency confers dramatic protection against atherosclerosis in the apoE^{-/-} mice. These results indicate a central role for AFABP in the development of major components of the metabolic syndrome through its distinct actions in adipocytes and macrophages.

Besides being active within the cell, AFABP appears to be a secreted protein (for normal levels and correlations with certain metabolic parameters see chapter 15). The extracellular role of secreted AFABP remains to be determined.

QUANTITATIVE DETERMINATION OF HUMAN ADIPOCYTE FABP

BioVendor Human Adipocyte FABP (FABP4) ELISA (RD191036200R)

Intended use

The RD191036200R Human AFABP (FABP4) ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human AFABP.

➤ **It is intended for research use only**

- The total assay time is less than 4 hours
- The kit measures total AFABP in serum and plasma (EDTA, citrate, heparin)
- Assay format is 96 wells
- Quality Controls are human serum based
- Standard is recombinant protein based
- Components of the kit are provided ready to use, concentrated or lyophilized
- Patent Application Number: DE 10 2005 034 788.6


Clinical application

- Energy metabolism and body weight regulation
- Diabetology

HUMAN ADIPOCYTE FABP (FABP4) ELISA CAT. NO.: RD191036200R	
Assay format	Sandwich ELISA, Biotin-labelled antibody, 96 wells/kit
Samples	Plasma (Citrate, EDTA, Heparin), Serum
Standards	0.5 to 25 ng/ml
Limit of detection	0.05 ng/ml

Test principle

In the BioVendor Human AFABP ELISA, standards, quality controls and samples are incubated in microplate wells pre-coated with polyclonal anti-human AFABP antibody. After 60 minutes incubation and washing, biotin labelled polyclonal anti-human AFABP antibody is added and incubated for 60 minutes with captured AFABP. After another washing, streptavidin-HRP conjugate is added. After 30 minutes incubation and the last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration of AFABP. A standard curve is constructed by plotting absorbance values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.



Precision

Intra-assay (Within-Run) (n=8)

Sample	Mean (ng/ml)	SD (ng/ml)	CV (%)
1	9.08	0.21	2.3
2	22.75	0.62	2.7

Inter-assay (Run-to-Run) (n=6)

Sample	Mean (ng/ml)	SD (ng/ml)	CV (%)
1	31.88	1.74	5.5
2	53.13	1.24	2.3

Spiking recovery

Serum samples were spiked with different amounts of human AFABP and assayed.

Sample	Observed (ng/ml)	Expected (ng/ml)	Recovery O/E (%)
1	18.63	-	-
	69.43	68.63	101.2
	45.31	43.63	103.9
	31.87	28.63	111.3
2	23.44	-	-
	69.84	73.44	95.1
	52.55	48.44	108.5
	35.95	33.44	107.5

Linearity

Serum samples were serially diluted with Dilution Buffer and assayed.

Sample	Dilution	Observed (ng/ml)	Expected (ng/ml)	Recovery O/E (%)
1	-	59.50	-	-
	2x	29.90	29.75	100.5
	4x	14.41	14.87	96.8
	8x	7.57	7.44	101.7
2	-	52.62	-	-
	2x	26.37	26.31	100.2
	4x	13.01	13.15	98.9
	8x	6.77	6.58	102.9

Effect of sample matrix

Citrate, heparin and EDTA plasmas were compared to respective serum samples from the same 10 individuals. Results are shown below:




Figure 3: AFABP levels measured using Human AFABP ELISA from 10 individuals using serum, EDTA, citrate and heparin plasma, respectively.

Volunteer No.	Serum (ng/ml)	Plasma (ng/ml)		
		EDTA	Citrate	Heparin
1	6.87	6.05	4.86	6.41
2	13.51	12.18	10.93	12.16
3	14.94	13.28	11.71	9.74
4	12.67	12.04	11.30	12.55
5	20.33	14.53	13.82	17.72
6	10.46	9.43	8.94	9.68
7	16.09	14.21	12.30	12.36
8	13.94	12.30	11.90	13.94
9	14.31	12.90	10.3	13.80
10	7.55	6.92	6.64	7.42
Mean (ng/ml)	13.07	11.38	10.33	11.58
Mean Plasma/Serum (%)	-	87.1	79.1	88.6
Coefficient of determination R ²	-	0.91	0.84	0.91

QUANTITATIVE DETERMINATION OF HUMAN ADIPOCYTE FABP

Summary of protocol

- Reconstitute Master Standard and prepare set of Standards
- Reconstitute QCs and dilute them 10., dilute samples 10.
- Add 100 µl Standards, QCs and samples
- Incubate at RT for 1 hour/300 rpm
- Wash plate 5 times
- Add 100 µl Biotin Labelled antibody
- Incubate at RT for 1 hour/300 rpm
- Wash plate 5 times
- Add 100 µl Streptavidin-HRP Conjugate
- Incubate at RT for 30 min/300 rpm
- Wash plate 5 times
- Add 100 µl Substrate Solution
- Incubate at RT for 10 min
- Add 100 µl stop solution
- Read absorbance and calculate results

Related products

- RD291036200R Adipocyte FABP Mouse ELISA
- RD172036100 Adipocyte Fatty Acid Binding Protein Human *E. coli* Tag free
- RD165036050 Adipocyte Fatty Acid Binding Protein NATIVE, Human Adipose Tissue
- RD181037050 Adipocyte Fatty Acid Binding Protein Human, Rabbit Polyclonal Antibody
- RD184037100 Adipocyte Fatty Acid Binding Protein Human, Goat Polyclonal Antibody

References

1. Makowski L, Brittingham KC, Reynolds JM, Suttles J and Hotamisligil GS: The Fatty Acid-binding Protein, aP2, Coordinates Macrophage Cholesterol Trafficking and Inflammatory Activity. *J Biol Chem.* Apr 1; 280(13): 12888-12895 (2005).
2. Maeda K, Cao H, Kono K, Gorgun CZ, Furuhashi M, Uysal KT, Cao Q, Atsumi G, Malone H, Krishnan B, Minokoshi Y, Kahn BB, Parker RA and Hotamisligil GS: Adipocyte/macrophage fatty acid binding proteins control integrated metabolic responses in obesity and diabetes. *Cell Metabolism.* February; 1(2): 107-119 (2005).
3. Boord JB, Maeda K, Makowski L, Babaev VR, Fazio S, Linton MF, Hotamisligil GS: Combined adipocyte-macrophage fatty acid-binding protein deficiency improves metabolism, atherosclerosis, and survival in apolipoprotein E-deficient mice. *Circulation.* Sep 14; 110(11): 1492-1498 (2004).
4. Lehmann F, Haile S, Axen E, Medina C, Uppenberg J, Svensson S, Lundback T, Rondahl L, Barf T: Discovery of inhibitors of human adipocyte fatty acid-binding protein, a potential type 2 diabetes target. *Bioorg Med Chem Lett.* Sep 6; 14(17): 4445-4448 (2004).
5. Damcott CM, Moffett SP, Feingold E, Barmada MM, Marshall JA, Hamman RF, Ferrell RE: Genetic variation in fatty acid-binding protein-4 and peroxisome proliferator-activated receptor gamma interactively influence insulin sensitivity and body composition in males. *Metabolism.* Mar; 53(3): 303-9 (2004).
6. Jenkins-Kruchten AE, Bennaars-Eiden A, Ross JR, Shen WJ, Kraemer FB, Bernlohr DA: Fatty acid-binding protein-hormone-sensitive lipase interaction. Fatty acid dependence on binding. *J Biol Chem.* Nov 28; 278(48): 47636-47643 (2003).
7. Hertz AV, Bennaars-Eiden A, Bernlohr DA: Increased lipolysis in transgenic animals overexpressing the epithelial fatty acid binding protein in adipose cells. *J Lipid Res.* Dec; 43(12): 2105-2111 (2002).
8. Fu Y, Luo N, Lopes-Virella MF, Garvey WT: The adipocyte lipid binding protein (ALBP/aP2) gene facilitates foam cell formation in human THP-1 macrophages. *Atherosclerosis.* Dec; 165(2): 259-269 (2002).
9. Fisher RM, Hoffstedt J, Hotamisligil GS, Thorne A, Ryden M: Effects of obesity and weight loss on the expression of proteins involved in fatty acid metabolism in human adipose tissue. *Int J Obes Relat Metab Disord.* Oct; 26(10): 1379-1385 (2002).
10. Boord JB, Maeda K, Makowski L, Babaev VR, Fazio S, Linton MF, Hotamisligil GS: Adipocyte fatty acid-binding protein, aP2, alters late atherosclerotic lesion formation in severe hypercholesterolemia. *Arterioscler Thromb Vasc Biol.* Oct 1; 22(10): 1686-1691 (2002).
11. Scheja L, Makowski L, Uysal KT, Wiesbrock SM, Shimshek DR, Meyers DS, Morgan M, Parker RA, Hotamisligil GS: Altered insulin secretion associated with reduced lipolytic efficiency in aP2-/- mice. *Diabetes.* Oct; 48(10): 1987-1994 (1999).
12. Milner KL, van der Poorten D, Xu A, Bugianesi E, Kench JG, Lam KS, Chisholm DJ, George J: Adipocyte fatty acid binding protein levels relate to inflammation and fibrosis in nonalcoholic fatty liver disease. *Hepatology.* Jun; 49(6):1926-34 (2009).
13. Yeung DC, Xu A, Tso AW, Chow WS, Wat NM, Fong CH, Tam S, Sham PC, Lam KS: Circulating levels of adipocyte and epidermal fatty acid-binding proteins in relation to nephropathy staging and macrovascular complications in type 2 diabetic patients. *Diabetes Care.* Jan; 32(1):132-4 (2009).
14. Stejskal D, Karpisek M, Bronsky J: Serum adipocyte-fatty acid binding protein discriminates patients with permanent and temporary body weight loss. *J Clin Lab Anal.* 22(5):380-2 (2008).
15. Karpisek M, Stejskal D, Kotolova H, Kollar P, Janoutova G, Ochmanova R, Cizek L, Horakova D, Yahia RB, Lichnovska R, Janout V: Treatment with atorvastatin reduces serum adipocyte-fatty acid binding protein value in patients with hyperlipidaemia. *Eur J Clin Invest.* Aug; 37(8):637-42 (2007).
16. Tso AW, Xu A, Sham PC, Wat NM, Wang Y, Fong CH, Cheung BM, Janus ED, Lam KS: Serum adipocyte fatty acid binding protein as a new biomarker predicting the development of type 2 diabetes: a 10-year prospective study in a Chinese cohort. *Diabetes Care.* Oct; 30(10):2667-72 (2007).
17. Yeung DC, Xu A, Cheung CW, Wat NM, Yam MH, Fong CH, Chau MT, Lam KS: Serum. Adipocyte fatty acid-binding protein levels were independently associated with carotid atherosclerosis. *Arterioscler Thromb Vasc Biol.* Aug; 27(8):1796-802 (2007).
18. Xu A, Tso AW, Cheung BM, Wang Y, Wat NM, Fong CH, Yeung DC, Janus ED, Sham PC, Lam KS: Circulating adipocyte-fatty acid binding protein levels predict the development of the metabolic syndrome: a 5-year prospective study. *Circulation.* Mar 27; 115(12): 1537-1543 (2007).
19. Bronsky J, Karpisek M, Bronska E, Pechova M, Jancikova B, Kotolova H, Stejskal D, Prusa R, Nevoral J: Adiponectin, adipocyte fatty acid binding protein, and epidermal fatty acid binding protein: Proteins newly identified in human breast milk. *Clin Chem.* Sep; 52(9): 1763-1770 (2006).
20. von Eynatten M, Schneider JG, Allolio B: Adipocyte fatty acid binding protein (A-FABP) ist ein neuer Serum-Marker für Adipositas, Insulinresistenz und Metabolisches Syndrom. *Diabetologie.* 1: S23-S172 (2006).
21. Xu A, Wang Y, Xu JY, Stejskal D, Tam S, Zhang J, Wat NM, Wong WK, Lam KS: Adipocyte fatty acid-binding protein is a plasma biomarker closely associated with obesity and metabolic syndrome. *Clin Chem.* Mar; 52(3): 405-413 (2006).
22. Stejskal D and Karpisek M: Adipocyte fatty acid binding protein (A-FABP) in Caucasian population: a new marker of metabolic syndrome? *Eur J Clin Invest.* Sep; 36(9): 621-625 (2006).

Contact Information



BioVendor – Laboratorni medicina a.s.
 Karasek 1767/1, 621 00 Brno, Czech Republic
 Phone: +420 549 124 185, Fax: +420 549 211 460
 E-mail: info@biovendor.com

BioVendor GmbH
 Otto-Hahn-Straße 16, 34123 Kassel, Germany
 Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
 E-mail: infoEU@biovendor.com

BioVendor GesmbH
 Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
 Phone: +43 1 890 9025, Fax: +43 1 890 5163
 E-mail: infoAustria@biovendor.com

BioVendor, LLC
 128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
 Phone: +1-800-404-7807, Phone: +1-828-575-9250
 Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com

➤ www.biovendor.com

*Visit **www.biovendor.com**
 to find more information
 about BioVendor products.*