

# QUANTITATIVE DETERMINATION OF HUMAN IL-1 BETA


**NEW PRODUCT**

## *Human Interleukin-1 beta ELISA*

- › Sensitivity (0.4 pg/ml)
- › Very good analytical characteristics
- › Validated for human serum, plasma (EDTA, citrate, heparin) and saliva samples

**CYTOKINES AND CHEMOKINES  
AND RELATED MOLECULES  
IMMUNE RESPONSE, INFECTION  
AND INFLAMMATION  
ONCOLOGY  
SEPSIS**

# HUMAN INTERLEUKIN-1 BETA ELISA



## Introduction

IL-1 Beta, a polypeptide cytokine, represents one of the most important mediators of inflammation and host responses to infections [1].

Even low concentrations of IL-1 Beta cause fever, hypotension and production of additional proinflammatory chemokines/ cytokines, such as IL-6 [2]. IL-1 Beta exerts biological effects by binding the membrane-bound type I IL-1 receptor (IL-1R), which then associates with the IL-1-receptor accessory protein (IL-1RAcP) to form a complex capable of intracellular signaling [3]. This signalling controls expression of a number of inflammatory and catabolic genes [4].

Besides its favorable role in mediating host responses to microbial invasion, IL-1 Beta has also harmful effects [5]. IL-1 Beta can promote tumor invasiveness, tumor angiogenesis and metastasis [6]. IL-1 Beta also exacerbates damage during

chronic diseases and acute tissue injury [7]. Overexpression of IL-1 Beta was observed in the pathophysiological changes that occur during different diseases, such as rheumatoid arthritis, neuropathic pain, inflammatory bowel disease, osteoarthritis, multiple sclerosis and neurodegenerative diseases [8, 9, 10].

It was observed that IL-1 Beta impairs insulin-producing Beta-cell function [11]. Macrophage-derived IL-1 Beta production in insulin-sensitive organs leads to progression of inflammation and induction of insulin resistance in obesity [12].

Regarding other biofluids, it was found that IL-1 Beta is one of the most abundant cytokine in saliva. It was observed that salivary level of IL-1 Beta was higher in the patients with periodontitis compared to periodontally healthy subjects [13, 14].

# QUANTITATIVE DETERMINATION OF HUMAN IL-1 BETA

## BioVendor Human Interleukin-1 beta ELISA RD194559200R

### Intended use

The RD194559200R Human Interleukin-1 beta ELISA is a sandwich enzyme immunoassay for the quantitative measurement of native human Interleukin-1 Beta Protein.

➤ **It is intended for research use only**

- The total assay time is less than 4 hours
- The kit measures total human IL-1 Beta in serum, plasma (EDTA, citrate, heparin) and saliva
- Assay format is 96 wells
- Standard is purified E.coli protein
- Components of the kit are provided ready to use, concentrated or lyophilized

### Clinical application


- Cytokines and chemokines and related molecules
- Immune Response, Infection and Inflammation
- Oncology
- Sepsis

### Test principle

In the BioVendor Human Interleukin-1 Beta ELISA, standards and samples are incubated in a microtiter plate wells pre-coated with monoclonal anti-human IL-1 Beta antibody. After 60 minutes incubation and a washing, biotin-labelled monoclonal anti-human IL-1 Beta antibody is added and incubated with captured IL-1 Beta for 60 minutes. After another washing, the streptavidin-HRP conjugate is added. After 60 minutes incubation and the last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration IL-1 Beta. A standard curve is constructed by plotting absorbance values against concentrations of IL-1 Beta standards, and concentrations of unknown samples are determined using this standard curve.

#### HUMAN INTERLEUKIN-1 BETA ELISA CAT. NO.: RD194559200R

Assay format	Sandwich ELISA, Biotin-labelled antibody, 96 wells/kit
Samples	Serum, Plasma (EDTA, citrate, heparin) and Saliva
Standards	1.25 to 80 pg/ml
Limit of detection	0.4 pg/ml



# HUMAN INTERLEUKIN-1 BETA ELISA

## Precision

Intra-assay (Within-Run) (n=8)

Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
Saliva 1	168.77	4.5	2.7
Saliva 2	45.67	1.9	4.2

Inter-assay (Run-to-Run) (n=7)

Sample	Mean (pg/ml)	SD (pg/ml)	CV (%)
Serum 1	54.93	2.8	5.2
Serum 2	101.57	6.8	6.7

## Spiking recovery

Serum samples were spiked with different amounts of human IL-1 Beta and assayed.

Sample	Observed (pg/ml)	Expected (pg/ml)	Recovery O/E (%)
Serum 1	2.36	-	-
	8.09	9.86	82.0
	14.26	17.36	82.1
	26.77	32.36	82.7
Serum 2	3.24	-	-
	9.44	10.74	87.9
	14.87	18.24	81.5
	26.92	33.24	81.0
Saliva 1	139.59	-	-
	200.98	189.59	106.0
	248.48	239.59	103.7
	371.43	339.59	109.4
Saliva 2	74.18	-	-
	128.06	124.18	103.1
	164.94	174.18	94.7
	287.52	274.18	104.9

## Linearity

Serum, plasma and urine samples were serially diluted with Dilution Buffer and assayed.

Sample	Dilution	Observed (pg/ml)	Expected (pg/ml)	Recovery O/E (%)
Serum 1	-	79.66	-	-
	2x	40.05	39.83	100.6
	4x	20.34	19.91	102.1
	8x	9.43	9.96	94.7
Serum 2	-	56.22	-	-
	2x	28.56	28.11	101.6
	4x	14.99	14.04	106.7
	8x	7.64	7.03	108.7
Saliva 1	-	1080.81	-	-
	2x	506.27	540.41	93.7
	4x	264.56	270.20	97.9
	8x	140.78	135.10	104.2
Saliva 2	-	663.22	-	-
	2x	307.99	331.61	92.9
	4x	153.41	165.81	92.5
	8x	80.40	82.90	97.0

# QUANTITATIVE DETERMINATION OF HUMAN IL-1 BETA

## Summary of protocol

- Reconstitute Master Standard and prepare set of Standards
- Dilute samples (serum 3x, plasma 3x, saliva 20x)
- Add 100 µl Standards and samples
- Incubate at RT for 1 hour/300 rpm
- Wash plate 3 times
- Add 100 µl Biotin Labelled Antibody solution
- Incubate at RT for 1 hour/300 rpm
- Wash plate 3 times
- Add 100 µl Streptavidin - HRP Conjugate
- Incubate at RT for 1 hour/300 rpm
- Wash plate 3 times
- Add 100 µl Substrate Solution
- Incubate at RT for 10 min
- Add 100 µl Stop Solution
- Read absorbance and calculate results

## References

1. Dinarello, Charles A. „Immunological and inflammatory functions of the interleukin-1 family.“ *Annual review of immunology* 27 (2009): 519-550
2. Church, Leigh D., Graham P. Cook, and Michael F. McDermott. „Primer: inflammasomes and interleukin 1 $\beta$  in inflammatory disorders.“ *Nature Clinical Practice Rheumatology* 4.1 (2008): 34-42
3. Korherr, Christian, et al. „A critical role for interleukin-1 receptor accessory protein in interleukin-1 signaling.“ *European journal of immunology* 27.1 (1997): 262-267
4. Risbud, Makarand V., and Irving M. Shapiro. „Role of cytokines in intervertebral disc degeneration: pain and disc content.“ *Nature Reviews Rheumatology* 10.1 (2014): 44-56
5. Eder, Claudia. „Mechanisms of interleukin-1 $\beta$  release.“ *Immunobiology* 214.7 (2009): 543-553
6. Apte, Ron N., et al. „The involvement of IL-1 in tumorigenesis, tumor invasiveness, metastasis and tumor-host interactions.“ *Cancer and Metastasis Reviews* 25.3 (2006): 387-408
7. Lopez-Castejon, Gloria, and David Brough. „Understanding the mechanism of IL-1 $\beta$  secretion.“ *Cytokine & growth factor reviews* 22.4 (2011): 189-195
8. Dinarello, Charles A. „Biologic basis for interleukin-1 in disease.“ *Blood* 87.6 (1996): 2095-2147
9. Braddock, Martin, and Anthony Quinn. „Targeting IL-1 in inflammatory disease: new opportunities for therapeutic intervention.“ *Nature Reviews Drug Discovery* 3.4 (2004): 330-340
10. Dinarello, Charles A. „Therapeutic strategies to reduce IL-1 activity in treating local and systemic inflammation.“ *Current opinion in pharmacology* 4.4 (2004): 378-385
11. Mandrup-Poulsen, T., et al. „Cytokines cause functional and structural damage to isolated islets of Langerhans.“ *Allergy* 40.6 (1985): 424-429
12. Maedler, Kathrin, et al. „Interleukin-1 beta targeted therapy for type 2 diabetes.“ *Expert opinion on biological therapy* 9.9 (2009): 1177-1188
13. Gursoy, Ulvi Kahraman, et al. „Salivary interleukin-1 $\beta$  concentration and the presence of multiple pathogens in periodontitis.“ *Journal of clinical periodontology* 36.11 (2009): 922-927
14. Tobón-Arroyave, S. I., P.E. Jaramillo-González, and D. M. Isaza-Guzmán. „Correlation between salivary IL-1 $\beta$  levels and periodontal clinical status.“ *Archives of oral biology* 53.4 (2008): 346-352

## Contact Information



**BioVendor – Laboratorni medicina a.s.**  
Karasek 1767/1, 621 00 Brno, Czech Republic  
Phone: +420 549 124 185, Fax: +420 549 211 460  
E-mail: info@biovendor.com

**BioVendor GmbH**  
Otto-Hahn-Straße 16, 34123 Kassel, Germany  
Phone: +49 6221 4339 100, Fax: +49 6221 4339 111  
E-mail: infoEU@biovendor.com

➤ [www.biovendor.com](http://www.biovendor.com)

**BioVendor GesmbH**  
Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria  
Phone: +43 1 890 9025, Fax: +43 1 890 5163  
E-mail: infoAustria@biovendor.com

**BioVendor, LLC**  
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA  
Phone: +1-800-404-7807, Phone: +1-828-575-9250  
Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com

Visit [www.biovendor.com](http://www.biovendor.com)  
to find more information  
about BioVendor products.