

BioVendor

Research
and Diagnostic Products



HUMAN PARATHYROID HORMONE (PTH) ELISA

Product Data Sheet

Cat. No.: RIS003R

For Research Use Only

CONTENTS

1.	INTENDED USE	3
2.	STORAGE, EXPIRATION	3
3.	INTRODUCTION	3
4.	TEST PRINCIPLE	4
5.	PRECAUTIONS	5
6.	TECHNICAL HINTS	5
7.	REAGENT SUPPLIED	6
8.	MATERIAL REQUIRED BUT NOT SUPPLIED	7
9.	PREPARATION OF REAGENTS	7
10.	PREPARATION OF SAMPLES	8
11.	ASSAY PROCEDURE	8
12.	CALCULATIONS	9
13.	PERFORMANCE CHARACTERISTICS	10
14.	QUALITY CONTROL	12
15.	REFERENCE INTERVALS	12
16.	REFERENCES	13

**»» This kit is manufactured by:
BioVendor – Laboratorní medicína a.s.**

»» Use only the current version of Product Data Sheet enclosed with the kit!

1. INTENDED USE

The RIS003R Human Parathyroid hormone ELISA Immunoenzymetric assay for the measurement of human Intact Parathyroid Hormone (PTH) in serum and plasma. It is intended for research use only.

2. STORAGE, EXPIRATION

Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.

Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.

After reconstitution, calibrators and controls should be frozen immediately after use and kept at -20°C for 3 months. Avoid successive freeze thaw cycles.

The concentrated Wash Solution is stable at room temperature until expiration date.

Freshly prepared Working Wash solution should be used on the same day.

After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.

Alterations in physical appearance of kit reagents may indicate instability or deterioration.

3. INTRODUCTION

Biological activities

Human parathyroid hormone (hPTH) is a major physiological regulator of phosphocalcic metabolism. hPTH increases serum calcium concentrations by its actions on kidney (enhancing tubular Ca⁺⁺ reabsorption and phosphate excretion) and bone (stimulating osteoclastic activity and bone resorption). It indirectly affects intestinal absorption of Ca⁺⁺ by stimulating renal 1 α -hydroxylation of 25 hydroxyvitamin D. The release of PTH is controlled in a negative feedback loop by the serum concentration of Ca⁺⁺.

PTH is synthesized in the chief cells of the parathyroid glands and secreted as an 84 amino acid molecule called "intact PTH", which is the main bioactive product. This molecule is degraded by proteolytic cleavage between amino acids 33-37 at peripheral sites to form biologically active amino-terminal fragments and biologically inactive carboxyl-terminal fragments. The carboxyl-terminal fragments are cleared only by glomerular filtration, while the bioactive intact PTH and amino-terminal fragments are also metabolically degraded in the liver and other tissues. The half-life of the carboxyl-terminal fragments increases dramatically in patients with renal failure. Thus, the measurement of intact PTH correlates best with the hormone production and biological activity.

Clinical application

The measurement of intact hPTH is used to establish the diagnosis of primary hyperparathyroidism by demonstrating elevated serum levels of bioactive PTH. It allows documenting the occurrence of secondary hyperparathyroidism in patients with Vit.D deficiency, intestinal malabsorption, or renal failure. In this last situation, the absence of interference with the inactive carboxyl-terminal fragments is especially valuable. The specificity and high sensitivity of the assay also allows differentiating clearly low serum PTH levels in hypoparathyroidism or in tumor-induced hypercalcaemia.

4. TEST PRINCIPLE

The Biovendor hPTH-ELISA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplates. Calibrators and samples react with the capture polyclonal antibodies (PAb, goat anti 1-34 PTH fragment) coated on microtiter well. After incubation, the excess of antigen is removed by washing. Then monoclonal antibodies (MAb, mouse anti 44-68 PTH fragment) labelled with horseradish peroxidase (HRP) are added. After an incubation period allowing the formation of a sandwich: coated PABs – human PTH – MAb – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. The chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the PTH concentration.

A calibration curve is plotted and PTH concentration in samples is determined by interpolation from the calibration curve. The use of the ELISA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range

5. PRECAUTIONS

Safety

For research use only

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl, the chromogen contains TMB. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves

6. TECHNICAL HINTS

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling.

Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.

Use a clean plastic container to prepare the Wash Solution.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.


To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).

Prepare a calibration curve for each run, do not use data from previous runs.

Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

7. REAGENT SUPPLIED

Reagents	96 tests Kit	Color Code	Reconstitution			
 Microtiterplate with 96 anti PTH (polyclonal antibodies) coated wells	96 wells	blue	Ready for use			
<table border="1" data-bbox="229 593 459 667"> <tr> <td>Ab</td> <td>HRP</td> </tr> </table> Conjugate: HRP labelled anti-PTH (monoclonal antibodies) in TRIS-maleate buffer with bovine serum albumin, thymol and sheep serum	Ab	HRP	1 vial 11 ml	red	Ready for use	
Ab	HRP					
<table border="1" data-bbox="194 891 375 958"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero calibrator in human serum and thymol	CAL	0	1 vials lyophilized	yellow	Add 3.0 ml distilled water	
CAL	0					
<table border="1" data-bbox="194 1070 375 1137"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrator N = 1 to 5 (see exact values on vial labels) in human serum and thymol	CAL	N	5 vials lyophilized	yellow	Add 1.0 ml distilled water	
CAL	N					
<table border="1" data-bbox="225 1317 534 1391"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum with thymol	CONTROL	N	2 vials lyophilized	silver	Add 1.0 ml distilled water	
CONTROL	N					
<table border="1" data-bbox="194 1512 418 1594"> <tr> <td>INC</td> <td>BUF</td> </tr> </table> Incubation Buffer with EDTA, Benzamidin and azide (< 0.1%)	INC	BUF	1 vial 6 ml	blue	Ready for use	
INC	BUF					
<table border="1" data-bbox="225 1720 568 1794"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash Solution (NaCl-Tween20)	WASH	SOLN	CONC	1 vial 25 ml	green	Dilute 28 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="225 1877 512 1928"> <tr> <td>CHROM</td> <td>TMB</td> </tr> </table> Chromogenic Solution (TMB: Tetramethylbenzidine)	CHROM	TMB	1 vial 25 ml	white	Ready for use	
CHROM	TMB					

STOP SOLN	1 vial 25 ml	white	Ready for use
Stopping reagent: HCl 1.0N			

- Note:**
1. Use the zero calibrator for sample dilutions.
 2. 1 pg of the calibrator preparation is equivalent to 1 pg of a synthetic PTH (1-84) from the Japanese Peptide Institute.

8. MATERIAL REQUIRED BUT NOT SUPPLIED

1. High quality distilled water
2. Pipettes for delivery of: 50 µl, 100 µl, 200 µl, 1 ml, 2 ml and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Horizontal microtiterplate shaker capable of 700 rpm ± 100 rpm
6. Washer for microtiterplates
7. Microtiterplate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (monochromatic reading)
8. Optional equipment: The ELISA-AID™ necessary to read the plate according to polychromatic reading (see paragraph XI.A.) can be purchased from Robert Maciels Associates, Inc. Mass. 0.2174 USA.

9. PREPARATION OF REAGENTS

- Calibrators: Reconstitute the zero calibrator with 3.0 ml distilled water and other calibrators with 1 ml distilled water.
- Controls : Reconstitute the controls with 1 ml distilled water.
- Working Wash solution : Prepare an adequate volume of Working Wash solution by adding 27 volumes of distilled water to 1 volume of Wash Solution (28x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

10. PREPARATION OF SAMPLES

Blood samples should be promptly separated from the blood cells.

Serum and plasma must be kept at 2 - 8°C.

If the test is not run within 8 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.

Prior to use, all samples should be at room temperature. It is recommended to vortex the samples before use.

It is advisable to assay serum samples.

Do not use haemolysed samples.

11. ASSAY PROCEDURE

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
 2. Secure the strips into the holding frame.
 3. Pipette 50 µl of Incubation Buffer into all wells.
 4. Pipette 200 µl of each Calibrator, Control and Sample into the appropriate wells.
 5. Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
 6. Aspirate the liquid from each well.
 7. Wash the plate 4 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
 8. Pipette 100 µl of anti-PTH-HRP conjugate into all the wells.
 9. Incubate for 1 hour at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.
 10. Aspirate the liquid from each well.
 11. Wash the plate 4 times by:
 - Dispensing 0.4 ml of Wash Solution into each well
 - Aspirating the content of each well
 12. Pipette 100 µl of the Chromogenic Solution into each well within 15 minutes following the washing step.
 13. Incubate the microtiterplate for 30 minutes at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm, avoid direct sunlight.
 14. Pipette 200 µl of Stop Solution into each well.
- Read the absorbancies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section 12.

12. CALCULATIONS

Polychromatic Reading:

In this case, the ELISA-AID™ software will do the data processing.

The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).

A second reading is performed at 490 nm against the same reference filter.

The ELISA-AID™ Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.

The principle of polychromatic data processing is as follows:

$X_i = \text{OD at 450 nm}$

$Y_i = \text{OD at 490 nm}$

Using a standard unweighted linear regression, the parameters A & B are calculated: $Y = A \cdot X - B$

If $X_i < 3$ OD units, then X calculated = X_i

If $X_i > 3$ OD units, then X calculated = $(Y_i - B) / A$

A 4 parameter logistic curve fitting is used to build up the calibration curve.

The PTH concentration in samples is determined by interpolation on the calibration curve.

Bichromatic Reading

Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).

Calculate the mean of duplicate determinations.

On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each calibrator against the corresponding concentration of PTH (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.

Read the concentration for each control and sample by interpolation on the calibration curve.

Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

Typical data

The following data are for illustration only and should never be used instead of the real time calibration curve.

hPTH-ELISA		OD units Polychromatic model
Calibrator	0 pg/ml	0.050
	22 pg/ml	0.149
	70 pg/ml	0.344
	224 pg/ml	0.999
	666 pg/ml	2.721
	1400 pg/ml	4.483

13. PERFORMANCE CHARACTERISTICS

- Detection Limit**

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 2 pg/ml.

- Specificity**

Cross-reactive hormones or fragments were added to the zero calibrator, a high value calibrator (900 pg/ml) and a low value calibrator (100 pg/ml). The apparent PTH response was measured.

Cross-reactant	No significant interference up to
PTH 1-34 synthetic fragment, human	1000 pg/ml
PTH 44-68 synthetic fragment, human	20000 pg/ml
PTH 53-84 synthetic fragment, human	20000 pg/ml
PTH 73-84 synthetic fragment, human	100000 pg/ml
PTH-related protein 1-34 synthetic fragment, human	100000 pg/ml

- Precision**

Intra-assay (n=10)

Serum	<X> ± SD (pg/ml)	CV %
A	41.0 ± 0.5	1.1
B	594 ± 12	2.0

SD : Standard Deviation; CV: Coefficient of variation

Inter-assay (n=20)

Serum	<X> ± SD (pg/ml)	CV %
A	45.7 ± 3.3	7.1
B	381 ± 11.1	2.9

SD : Standard Deviation; CV: Coefficient of variation

- Accuracy**

Recovery test

Sample	Added PTH (pg/ml)	Recovered PTH (pg/ml)	Recovery (%)
Serum	371	333	90
Heparin plasma	371	347	93
EDTA plasma	371	350	94

Dilution test

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
Serum	1/1	-	955
	1/2	477	506
	1/4	239	229
	1/8	119	124
	1/16	60	63
	1/32	34	34

Samples were diluted with zero calibrator.

- Time delay between last calibrator and sample dispensing**

As shown hereafter, assay results remain accurate even when a sample is dispensed 60 minutes after the calibrators have been added to the coated wells.

TIME DELAY					
	T0	15 min	30 min	45 min	60 min
S1	158	146	166	144	137
S2	77	72	70	67	72
S3	328	320	323	342	357
S4	260	250	250	258	251

14. QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

15. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

	N	Median (pg/ml)	Range (pg/ml)
Normal patients	156	29	16 - 46
Hyperparathyroidism	64	291	106 - > 1000
Hypoparathyroidism	11	0	0 – 6.4

The range is based on 5% to 95% percentiles.

16. REFERENCES

»» References to PTH:


1. HABENER J.F., and POTTE J.T., Jr. (1978) **"Biosynthesis of parathyroid hormone"**. New Engl. J. Med., 299, 11:580 and 299, 12:635.
2. MARTIN K.J., HRUSKA K.A., FREITAG J.J., KLAH S. and SLOPOLOSKY E. (1979) **"The peripheral metabolism of parathyroid hormone"**. New Engl. J. Med., 301, 20:1092.
3. GOLTZMAN D., HENDERSON B. and LOVERIDGE N. "Cytochemical bioassay of PTH. (1980) **Characteristics of the assay and analysis of circulating hormone forms"**. J. Clin. Invest., 65:1309.
4. POTTS J.T. Jr., KRONENBERG H.M., ROSENBLATT M. (1982) **"Parathyroid hormone : Chemistry, biosynthesis and mode of action"**. Adv. Protein Chem., 323.
5. HACKENG W.H.L., LIPS P., NETELENBOS J.C. and LIPS C.J.M. (1986) **"Clinical implication of estimation of intact parathyroid hormone (PTH) versus total immunoreactive PTH in normal subjects and hyperparathyroid patients"**. J. Clin. Endocrinol. Metab., 63:447.
6. BOUILLON R., COOPMANS W., DE GROOTE D.E.H., RADOUX D., ELIARD P.H. (1990) **"Immunoradiometric assay of Parathyrin with polyclonal and monoclonal region specific antibodies"**. Clin. Chem., 36/2:271-276.


»» For more references on this product see our WebPages at www.biovendor.com

Assay procedure summary

	CALIBRATORS (μl)	SAMPLE(S) CONTROLS (μl)
Incubation buffer	50	50
Calibrators (0-5)	200	-
Samples, Controls	-	200
Incubate for 2 hours at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 4 times with 400 μ l of Wash Solution and aspirate.		
Anti-PTH-HRP	100	100
Incubate for 1 hour at room temperature with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 4 times with 400 μ l of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 30 min at room temperature with continuous shaking at 700 rpm.		
Stop Solution	200	200
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (versus 630 or 650 nm) and 490 nm (versus 630 or 650 nm)		

NOTES





BioVendor – Laboratorní medicína a.s.

Karasek 1767/1, 621 00 Brno, Czech Republic

Phone: +420-549-124-185, Fax: +420-549-211-460

E-mail: info@biovendor.com, sales@biovendor.com

Web: www.biovendor.com

There are BioVendor branches and distributors near you.

To find the office closest to you, visit www.biovendor.com/contact

