

Bicarbonate Assay Kit (CO₂)

Method: PEPC Enzymatic

| Cat .No. | Size | Instrument |
|----------|----------|--|
| GB460E | 5×20 ml | For Hitachi 717 & Shimadzu CL7200/8000 |
| GB460E/B | 6×80 ml | |
| GS461E | 5×20 ml | For Hitachi 917 & OlympusAU640/400/600 |
| GS461E/B | 8×60 ml | |
| GH461E | 5×20 ml | Hitachi 902 |
| GT461E | 5×20 ml | For TOSHIBA 40 |
| GX461E | 1×100 ml | For SYNCHRON CX4-5-7-9/ LX20/DXC600-800 |



INTENDED USE

For the *in vitro* quantitative determination of CO₂ in serum or plasma.

CLINICAL SIGNIFICANCE^[1]

Increased blood CO₂, (hypercapnia) causes respiratory acidosis. CO₂ rises with decreased alveolar ventilation due to diseases of the lungs or bronchial tree, or breathing CO₂ enriched air. Depression of the overall lung capacity by certain drugs may lead to retention of CO₂.

PRINCIPLE^[2, 3]



The reduction in absorbance at 340 nm caused by the oxidation of NADH is proportional to the bicarbonate concentration in the sample.

SPECIMEN COLLECTION AND PREPARATION^[4]

Serum or heparinized plasma may be used. EDTA, citrate and oxalate should not be used as anticoagulants, as they will affect results. Samples should be drawn on ice and analyzed within 1 hour. Samples should be kept tightly closed, as CO₂ will diffuse from the sample causing erroneous values (up to 6 mmol/hr).

REAGENT COMPOSITION

| Contents | Concentration of Solutions |
|-----------------------|----------------------------|
| Reagent 1 (R1) | |
| Tris Buffer | PH 7.5 |
| PEP | 12.5 mmol/L |
| NADH analog | 0.6 mmol/L |
| MDH | >4100 U/L |
| PEPC | >400 U/L |

STABILITY AND PREPARATION OF REAGENTS

All reagents are ready to use. Stable up to the expiry date when stored at 2-8°C. The reagent is stable for 28 days On-board the analyzer after opening and kept at 2-8°C.

ASSAY PROCEDURE


Test Procedure for Analyzers (HITACHI7170/917)

Assay Mode: 2 Point Rate, 4 - 10

Wave Length (main/sub): 405 nm/546 nm

Sample: 2 µl

R1: 200 µl



- Mix 2 µl sample with 200 µl R1 and incubate at 37 °C for 1 minute.
- Read initial absorbance and start timer simultaneously, read again after 1 and 2 minutes.

MATERIALS REQUIRED BUT NOT PROVIDED

Randox Assayed Multi-sera Level 2 (Cat .No.HN 1530) and Level 3 (Cat .No. HE 1532).

CALCULATION

$$\text{Concentration} = \frac{\Delta A_{\text{sample}}/\text{min}}{\Delta A_{\text{calibrator}}/\text{min}} \times \text{Calibrator value}$$

CALIBRATION

Recommend that this assay should be calibrated using Gcell Calibrator (Cat .No. GC-CO₂).

QUALITY CONTROL

Use Gcell multi-control serum or Randox control serum. Values obtained should fall within a specified range. If these values fall outside the range and, the following steps should be taken:

- Check instrument settings and light source.
- Check reaction temperature.
- Check expiration date of kit and contents.

NORMAL VALUE^[5]

Serum or plasma: Venous 22 - 29 mmol/L.

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

SPECIFIC PERFORMANCE CHARACTERISTICS

LINEARITY

The method is linear up to 50 mmol/L. Sample above this concentration should be diluted with 0.9% NaCl and reassay. Multiply the result by

dilution factor.

PRECISION

The CV of this test should be less than 10%.

| Intra assay precision | | |
|-----------------------|--------|--------|
| N=20 | Level1 | Level2 |
| Mean | 13.1 | 19.6 |
| SD | 0.24 | 0.19 |
| CV | 1.82% | 0.96% |
| Inter assay precision | | |
| N=5 | Level1 | Level2 |
| Mean | 16.2 | 19.0 |
| SD | 0.48 | 0.49 |
| CV | 2.98% | 2.62% |

SENSITIVITY

The minimum detectable level that can be distinguished from zero has been determined as 2 mmol/L.

INTERFERENCE










The main interference in this assay is CO₂ from the air or from the breath of the analyst. The assay is not affected by the following interfering substances at the indicated concentration:

| | |
|-------------------|------------|
| Hemoglobin: | 400 mg/dl |
| Direct bilirubin: | 40 mg/dl |
| Intralipid: | 1000 mg/dl |
| Ascorbic acid: | 50 mg/dl |

REFERENCES

1. Tietz, N. N., et al "Textbook of Clinical Chemistry" W. B. Saunders Co., 1986; 1172-1253.
2. Jacobs, N., et al "Laboratory Test Handbook" 2nd. ed., Williams and Wilkins 1990.
3. Forrester, R.L., Wataji, L.J., Silverman, D.A., Pierre K.J., Clin, Chem. 1976; 22/2: 243-245.
4. Young D.S., Effects of Drugs on Chemical Laboratory Tests, 3rd ed., AACCC Press 1990.
5. Norris, K.A., Atkinson, A.R., Smith, W.G., Clin. Chem. 1975; 21/8: 1093 - 1101.

INDEX OF SYMBOLS

| | |
|---|---|
|  | Manufacture |
|  | Catalogue Number |
|  | Lot number |
|  | Date of manufacture |
|  | Use by(Expiration date) |
|  | For In-Vitro Diagnostic use only |
|  | Stored at 2-8 °C |
|  | Attention:See instruction for use |
|  | Authorized Representative in the European Company |