

## Zinc Assay Kit (Zn)

Method: PAPS Reagent

Cat .No.	Size	Instrument
GB430E	R1: 2×40 ml R2: 1×20 ml	For Hitachi 717 & ShimadzuCL7200/8000
GS431E	R1: 2×40 ml R2: 1×20 ml	For Hitachi917 & OlympusAU640/400/600
GH431E	R1: 2×40 ml R2: 1×20 ml	For Hitachi902
GX431E	R1: 1×80 ml R2: 1×20 ml	For SYNCHRON CX4-5-7-9/LX20/DXC600-8
GD431E	R1: 24×4.3 ml R2: 6×4.3 ml	For DATE DEMENSION

### INTENDED USE

For the *in vitro* quantitative determination of Zinc in serum.

### CLINICAL SIGNIFICANCE<sup>[1]</sup>

Zinc is an essential trace metal, which is second only to Iron. It is present in Zinc metalloenzymes e.g. Carbonic anhydrase, alkaline phosphatase, R.N.A and D.N.A polymerases, thymidine kinase, carboxypeptidases and alcohol dehydrogenase.

### ASSAY PRINCIPLE<sup>[2]</sup>

Zinc present in the sample is chelated by Nitro-PAPS in the reagent. The formation of this complex is measured at a wavelength of 570 nm.

### SAMPLE COLLECTION AND PREPARATION

Serum.

Serum samples are stable for one week at 4°C.

### REAGENT COMPOSITON

Contents	Concentration of Solutions
<b>R1</b>	
borate buffer, pH 8.20	370 mmol/L
salicylaldoxime	12.5 mmol/L
dimethylglyoxime	1.25 mmol/L
surfactants	
preservatives	
<b>R2</b>	
Nitro-PAPS	0.40 mmol/L
<b>Calibrator</b>	200 µg/dl

### STABILITY AND PREPARATION OF REAGENTS

All reagents are ready to use.

Unopened kit: Up to the expiry date at 2-8°C.


The reagents are stable for 1 month after opening and kept at 2-8°C.

### ASSAY PROCEDURE

Test Procedure for Analyzers (HITACHI 917)

Assay Mode: End Point, 34

Wave Length (main): 570 nm



- Mix 12 µl sample with 200 µl R1 and incubate at 37°C for 5 minutes.
- Add 50 µl R2 into cuvette, mix and incubate for another 5 minutes at 37°C
- Measure the absorbance of the sample ( $A_{\text{sample}}$ ) and calibrator ( $A_{\text{calibrator}}$ ) against reagent blank.

### CALCULATION

$$\text{Concentration} = \frac{A_{\text{sample}} - A_{\text{blank}}}{A_{\text{calibrator}} - A_{\text{blank}}} \times \text{Calibrator value}$$

### CALIBRATION

Recommend that this assay should be calibrated using Gcell Calibrator (Cat .No: GC-Zn) .

### QUALITY CONTROL

Randox Assayed Multi-sera, Level 2 and Level 3 are recommended for daily quality control. Two levels of controls should be assayed at least once a day. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error, the following steps should be taken:

- Check instrument settings and light source.
- Check reaction temperature.
- Check expiration date of kit and contents.

### NORMAL VALUE

Serum :70-115 µg/dl (10.7-17.7 µmol/L)

### CONVERSION FACTORS

µmol/L × 6.51 = µg/dl

### SPECIFIC PERFORMANCE CHARACTERISTICS

#### LINEARITY

The method is linear up to 500 µg/dl. If the concentration exceeds the top value, further dilute the sample 1+1 with 0.9% NaCl solution and repeat assay. Multiply the result by 2.

#### PRECISION

The CV of the test should be less than 5%.

Intra assay precision		
N=20	Level 1	Level 2
Mean	24.76	20.55
SD	0.78	0.96
CV	3.13%	4.66%
Inter assay precision		
N=5	Level1	Level 2
Mean	24.51	21.07
SD	0.94	0.92
CV	3.84%	4.37%

**SENSITIVITY**

The minimum detectable concentration of Zn with an acceptable level of precision was determined as 6.92 ug/dl.

**INTERFERENCE**

Bilirubin  $\leq 50$  mg/dl, Ascorbic Acid  $\leq 200$  mg/dl,  
Heparin sodium  $\leq 100$  U/ml, without interference;  
Intralipid  $\leq 500$  mg/dl, Interference.

**CORRELATION**

This method (Y) was compared with another commercially available method (X) and the following linear regression equation obtained:

$Y=1.073X-1.609$ ,  $R^2=0.997$ ; 237 patient samples were analyzed.










**SAFETY PRECAUTIONS AND WARNINGS**

1. For in vitro diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
2. The reagents contain Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
3. Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
4. All specimens used in this test should be considered potentially infectious. Universal Precautions, as they apply at your facility, should be used for handling and disposing of materials during and after testing.

**REFERENCES**

1. M. Saito, T. Makino et al., Clinica Chimica Acta, 120(1982) 127-135.
2. R. Homster, B. Zak, Clin. Chem. 31/8, 1310-1313(1985).

**INDEX OF SYMBOLS**

	Manufacture
	Catalogue Number
	Lot number
	Date of manufacture
	Use by(Expiration date)
	For In-Vitro Diagnostic use only
	Stored at 2-8°C
	Attention:See instruction for use
	Authorized Representative in the European Company