

MYOGLOBIN ASSAY

Dual Vial Liquid Stable

Diazyme's Myoglobin Assay is an excellent cost effective cardiovascular test that is used to aid in the early detection of myocardial damage. The latex enhanced immunoturbidimetric methodology is highly precise with excellent correlations to existing commercial myoglobin tests and is designed to work on most modern high throughput clinical chemistry analyzers. The Myoglobin assay offers liquid stable reagent, calibrator, high and low end controls separately for added convenience.

DIAZYME MYOGLOBIN ASSAY ADVANTAGES

- Fast test results for a rapid turnaround time
- Liquid stable reagent, calibrator and controls are offered separately for added convenience
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

REGULATORY STATUS

510(k) Cleared



AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- Roche
- Hitachi



MYOGLOBIN ASSAY


Dual Vial Liquid Stable



ASSAY SPECIFICATIONS

Method	Latex Enhanced Immunoturbidimetric
Sample Type & Volume	<ul style="list-style-type: none">• Serum• Plasma<ul style="list-style-type: none">- Li-heparin- K₃EDTA Sample Volume 4 µL
Method Comparison	Deming Regression: N = 66 y-intercept = -5.141 Slope = 0.959 R ² = 0.9927 Regression Analysis: N = 66 y-intercept = -4.228 Slope: 0.9526 R ² : 0.9855 Samples Ranged From: 16.9 – 615.9 ng/mL
Assay Range	13.2 – 615.9 ng/mL
LOB	4.4 ng/mL
LOD	7.2 ng/mL
LOQ	13.2 ng/mL
Calibration Levels	5-Point Calibration
Reagent On-Board Stability	Opened: Eight weeks on board analyzer

Myoglobin Assay Procedure*



*Analyzer Dependent

Parameter questions for Myoglobin Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email support@diazyme.com

ASSAY PRECISION

The precision of the Diazyme Myoglobin Assay was evaluated according to CLSI EP5-A guideline. In the study, three levels of serum based controls containing approximately 66, 170, and 335 ng/mL of myoglobin, and three serum sample containing approximately 35, 150, and 414 ng/mL of myoglobin, respectively, were tested with 2 runs per day in duplicates over 20 working days. Results were calculated using the EP Evaluator software precision statistic template and summarized in the following tables:

	Control Level 1	Control Level 2	Control Level 3	Serum Level 1	Serum Level 2	Serum Level 3
N	80	80	80	80	80	80
Mean	65.97	175.8	337.0	37.78	148.6	414.3
SD	2.45	6.69	11.9	1.77	3.53	19.7
CV%	3.71%	3.87%	3.54%	4.69%	2.37%	4.80%

	Control Level 1	Control Level 2	Control Level 3	Serum Level 1	Serum Level 2	Serum Level 3
N	80	80	80	80	80	80
Mean	65.97	172.8	337.0	37.78	148.6	414.3
SD	3.37	7.37	14.9	1.97	5.32	21.8
CV%	5.10%	4.30%	4.40%	5.20%	3.58%	5.3%

ASSAY INTERFERENCE

The following substances do not interfere with this assay at the levels tested (less than 10% bias):

Ampicillin:	up to 152 µM	Hemoglobin:	up to 1000 mg/dL
Carbamazepine:	up to 0.13 mM	Bilirubin:	up to 40 mg/dL
Na2+-Cefoxitin:	up to 1549 µM	Conjugated Bilirubin:	up to 40 mg/dL
Ibuprofen:	up to 2425 µM	Triglycerides:	up to 1000 mg/dL
Cyclosporin:	up to 0.125 µM	Intralipid:	up to 125 mg/dL
Levodopa:	up to 30.4 mM	Ascorbic acid:	up to 176 mg/dL
Methyldopa:	up to 71 µM	Rheumatoid factor:	up to 100 IU/mL
Metronidazole:	up to 701 µM	Heparin:	up to 1.5 IU/mL
Rifampicin:	up to 78.1 µM	N-acetylcysteine:	up to 11.04 mM
Theophylline:	up to 222 µM	Acetylsalicylic acid:	up to 2.78 mM
Phenylbutazone:	up to 650 µM		
Valproic Acid, Sodium Salt:	up to 3.5 mM		

DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064

PO Box 85608, San Diego, CA 92186

Tel: 858-455-4768 888-DIAZYME

www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland

Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366

sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201, 1011 Halei Road, Zhangjiang Hi-tech Park

Shanghai, 201203, People's Republic of China

Tel: 086-21-51320668 Fax: 086-21-51320663

www.lanyuanbio.com service@lanyuanbio.com

