

HOMOCYSTEINE 3 REAGENT ENZYMATIC ASSAY

Three Reagent Liquid Stable

Diazyme's Homocysteine 3 Reagent Enzymatic Assay features convenient ready to use reagent, calibrators and controls for the quantitative determination of total L-homocysteine in serum or plasma. Diazyme's proprietary Enzyme Cycling methodology is an excellent choice for cost conscious laboratories of all sizes due to a wide variety of instrument specific packaging options. The assay requires minimal patient sample and provides fast, accurate and precise results. A wide variety of reliable instrument parameters means the assay is readily available for installation on most automated clinical chemistry analyzers.

DIAZYME HOMOCYSTEINE 3 REAGENT ASSAY ADVANTAGES

- Award winning Homocysteine recognized by the American Association of Clinical Chemistry (AACC) for outstanding contribution to scientific research
- Innovative enzyme cycling based technology for accurate and reliable results
- Excellent correlation to HPLC and immunochemical methods
- No "carry over" issues with iron or lipase reagents
- Test renal patients with confidence since there is no interference from cystathionine which affects some other less specific methods
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation saving time and reducing sample handling

REGULATORY STATUS


510(k) Cleared

Health Canada Registered



AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- **Roche**
- Modular P
- Integra
- Cobas
- Hitachi
- **Beckman**
- Synchron
- **Siemens**
- Dimension



HOMOCYSTEIN 3 REAGENT ENZYMATIC ASSAY


Three Reagent Liquid Stable



ASSAY SPECIFICATIONS

Method	Diazyme Patented Enzyme Cycling
Sample Type & Volume	<ul style="list-style-type: none">• Serum• Plasma- EDTA- Li-heparin Sample Volume 18 μ L
Method Correlation	N = 66 y-intercept = 0.87 Slope = 0.98 R ² = 0.976
Linear Range	up to 50 μ mol/L
LOD	<1.5 μ mol/L
Calibration Levels	5-Point Calibration
Reagent On-Board Stability	Opened: At least 100 days (Analyzer Dependent)

Homocysteine 3 Reagent Assay Procedure*



*Analyzer Dependent

Parameter questions for Enzymatic Homocysteine 3 Reagent Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email support@diazyme.com

1. Vilaseca et al. *Clin. Chem.* 43: 690-692 (1997)
2. Faure-Delanef et al. *Am. J. Hum. Genet.* 60: 999-1001 (1997)

ASSAY PRECISION

Precision studies were conducted according to the NCCLS EP-5 protocol. Within precisions (CV%) for three levels of Hcy controls are 2.2% for 7 μ M Hcy, 3.0% for 12 μ M Hcy and 1.8% for 29.5 μ M Hcy. Total imprecision for three levels of Hcy controls are 4.1% for 7 μ M Hcy, 5.9% for 12 μ M Hcy and 4.0% for 29.5 μ M Hcy.

HCY Concentration	7 μ M N = 40	12 μ M N = 80	29.5 μ M N = 80
Within-Run Imprecision CV%	2.2	3.0	1.8
Total Imprecision CV%	4.1	5.9	4.0

ASSAY INTERFERENCE

An interference study was performed by testing a serum sample spiked with varied concentrations of endogenous substances. The following substances normally present in the serum produced less than 10% deviation when tested at the stated concentrations:

NH ₄ Cl:	500 μ M
NaPi:	1 mM
NaF:	1 mM
Triglycerides:	2500 mg/dL
Bilirubin:	20 mg/dL
Hemoglobin:	1200 mg/dL
*Glutathione:	0.5 mM
Ascorbic Acid:	10 mM
L-Cysteine:	1 mM
S-Adenosylmethionine (SAM):	20 μ M
**Adenosine:	100 μ M
**Cystathionine:	100 μ M

* Glutathione was originally tested at 1 mM level, the interference was +13.5%. When retested at 0.5 mM level, the interference was less than 10%.

** The concentrations tested are about 5-10 times higher than the normal range of serum levels.

REFERENCE RANGE

In most of the U.S. clinical laboratories, 15 μ mol/L is used as the cut-off value for normal level of Hcy for adults.¹⁻² In Europe, 12 μ mol/L is used as the cut-off value. However, each laboratory is recommended to establish a range of normal values for the population in their region.

DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064
PO Box 85608, San Diego, CA 92186
Tel: 858-455-4768 888-DIAZYME
www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park
Shanghai, 201203, People's Republic of China
Tel: 086-21-51320668 Fax: 086-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

