

CYSTATIN C ASSAY

Dual Vial Liquid Stable

Cystatin C is an emerging renal Biomarker for eGFR and is recommended as part of the KDIGO 2012 guidelines for the early confirmation and diagnosis of Chronic Kidney Disease (CKD).¹ In addition to confirming CKD at earlier stages than is possible with serum Creatinine, studies suggest that Cystatin C may help facilitate kidney disease screening efforts in the elderly, and those with diabetes, hypertension, or cardiovascular disease.²⁻⁷ The Diazyme Cystatin C Assay is a cost effective dual vial liquid stable system which is directly traceable to (ERM/DA471/IFCC) the international standard reference material. The test utilizes Avian IGY antibodies to virtually eliminate some of the most common causes for interference in immunoassay's.

DIAZYME CYSTATIN C ASSAY ADVANTAGES

- Diazyme's Cystatin C Assay uses Avian IGY antibodies which are not interfered with by rheumatoid and HAMA factors for increased reliability
- The assay is traceable to (ERM/DA471/IFCC) the international standard reference material providing increased accuracy and performance
- Liquid stable reagent, calibrator and controls are offered separately for added convenience
- A wide range of instrument parameters are offered for facilitating and simplifying implementation

REGULATORY STATUS

510(k) Cleared **CE**

AVAILABLE INSTRUMENT SPECIFIC PACKAGING


- | | | |
|----------------|------------------|------------------|
| • Roche | • Beckman | • Siemens |
| - Hitachi | - Synchron | - Dimension |
| | - AU Series | |



ASSAY SPECIFICATIONS

Method	Latex Enhanced Immunoturbidometric (Avian IGY antibodies quantification at 540 nm)
Sample Type & Volume	<ul style="list-style-type: none"> • Serum • Plasma - Heparin - EDTA Sample Volume 3 µL
Method Correlation	N = 45 y-intercept = 0.0715 Slope = 0.999 R ² = 0.992
Linear Range	0.2 to 8.0 mg/L
LOB	0.04 mg/L Cystatin C
LOQ	0.19 mg/L Cystatin C
LOD	0.068 mg/L Cystatin C
Calibration Levels	5-Point Calibration
Traceability	Standard traceable to ERM-DA471/IFCC primary reference material
Reagent On-Board Stability	Unopened: 24 months when stored at 2-8°C Opened: 4 weeks when kept stored at 2-8°C on-board Hitachi 917

Cystatin C Assay Procedure*



*Analyzer Dependent

Parameter questions for Cystatin C Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email support@diazyme.com

- Kidney International Supplements (2013) 3, 136–150; doi:10.1038/kisup.2012.72; 5-8.
- Filler G, Bökenkamp A, Hofmann W, Le Bricon T, Martínez-Brú C, Grubb A. Cystatin C as a marker of GFR - history, indications, and future re-search. Clin Biochem 2005; 38: 1-8.
- Dharmidharka VR, Kwon C, Stevens G. Serum cystatin C is superior to serum creatinine as a marker of kidney function: a meta-analysis. Am J Kidney Dis 2002; 40: 221-226.
- Grubb A, Björk J, Bondesson P, Lindström V, Sterner G, Nyman U. Cystatin C estimates glomerular filtration rate better than creatinine clear-ance using the Cockcroft-Gault formula. Scand J Clin Lab Invest 2005; 65: 1-10.
- Thomassen SA, Johannessen IL, Erlandsen EJ, Abrahamsen J, Randers E. Serum cystatin C as a marker of the renal function in patients with spinal cord injury. Spinal Cord 2002; 40: 524-528.
- Christensson AG, Grubb AO, Nilsson JA, Norrgren K, Sterner G, Sundkvist G. Serum Cystatin C advantageous compared with serum creatinine in the detection of mild but not severe diabetic nephropathy. J Intern Med 2004; 256: 510-518.
- Alan H.B. Wu. Tietz Clinical Guide to Laboratory Tests. Fourth Ed. Saunders Elsevier, 11830 Westline Industrial Drive, St. Louis, Missouri 63146. 2006; 328-329.
- Grubb A, Nyman U, Björk J, et al. Simple Cystatin C-based prediction equations for glomerular filtration rate compared with the modification of diet in renal disease prediction equation for adults and the Schwartz and the Counahan-Barratt prediction equations for children. Clin Chem 2005; 51:1420-31.
- Inker L, Brad A., Chester, F., et al. "KDIGO US Commentary on the 2012 KDIGO Clinical Practice Guideline for the Evaluation and Management of CKD." Am J Kidney Dis. 2014; 63(5):713-735.

ASSAY PRECISION

The precision of the Diazyme Cystatin C Assay was evaluated according to Clinical Laboratory Standards Institute (formerly NCCLS) EP5-A guideline. In the study, three samples containing Cystatin C were tested on Hitachi 917 2 runs per day in duplicates over 20 working days.

Within Run Precision (S _r)	Level 1 0.9 mg/L Cystatin C	Level 2 2.5 mg/L Cystatin C	Level 3 5.4 mg/L Cystatin C
No. of Data Points	80	80	80
Mean (mg/L)	0.91	2.51	5.40
SD (mg/L)	0.03	0.06	0.11
CV%	3.5%	2.5%	2.0%

Within Laboratory Precision (S _T)	Level 1 0.9 mg/L Cystatin C	Level 2 2.5 mg/L Cystatin C	Level 3 5.4 mg/L Cystatin C
No. of Data Points	80	80	80
Mean (mg/L)	0.91	2.51	5.40
SD (mg/L)	0.04	0.08	0.25
CV%	4.6%	3.0%	4.6%


ASSAY INTERFERENCE

The following substances do not interfere with this assay at the levels tested (less than 10% bias):

- Hemoglobin: up to 1000 mg/dL
- Bilirubin: up to 40 mg/dL
- Bilirubin Conjugated: up to 40 mg/dL
- Triglycerides: up to 1000 mg/dL
- Ascorbic Acid: up to 176 mg/dL
- Rheumatoid Factor: up to 1000 IU/mL

ASSAY REFERENCE RANGE

The reference interval is 0.5 - 1.03 mg/L. However, each laboratory is recommended to establish a range of normal values for the population in their region.⁸⁻⁹



(n = 451) Population GFR = 83.93 × Cystatin C [mg/L]⁻¹ - (1.68)

DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064
 PO Box 85608, San Diego, CA 92186
 Tel: 858-455-4768 888-DIAZYME

www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland
 Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park
 Shanghai, 201203, People's Republic of China
 Tel: 086-21-51320668 Fax: 086-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

