

MYELOPEROXIDASE (MPO) ASSAY

A Sensitive Inflammatory Biomarker

Myeloperoxidase (MPO) is a hemoprotein present in leukocytes of blood circulation. It is well known that in published literature, elevated levels of plasma MPO is a sensitive indicator of inflammatory disorders.¹⁻⁶ MPO is involved in the oxidation of lipids contained within LDL particles, and its reaction products including hydrogen peroxide are involved in the initiation of systemic inflammation.⁷⁻¹⁰ Diazyme's Latex Enhanced Immunoturbidimetric MPO Assay is accurate, cost effective and designed to work on validated general chemistry analyzers.

DIAZYME MYELOPEROXIDASE (MPO) ASSAY ADVANTAGES

- The MPO assay has been designed to work on most modern high throughput general chemistry analyzers
- Faster reporting and improved workflow for research laboratories
- Automated parameters available for a wide range of clinical instrumentation
- Liquid stable format requires no reagent preparation


REGULATORY STATUS

USA: For Research Use Only



AVAILABLE INSTRUMENT SPECIFIC PACKAGING

- Roche
- Hitachi



MYELOPEROXIDASE (MPO) ASSAY


Dual Vial Liquid Stable



ASSAY SPECIFICATIONS

Method	Latex Enhanced Immunoturbidimetric Assay
Sample Type & Volume	<ul style="list-style-type: none">• Plasma- Lithium Heparin- EDTA Sample Volume 15 μ L
Method Comparison	N = 54 y-intercept = 35.4 pmol/L Slope = 1.01 R ² = 0.98 Samples Ranged From: 31.4 ng/mL to 715.6 ng/mL
Linear Range	83 to 5000 pmol/L
LOQ	12 ng/mL
Calibration Levels	5-Point Calibration

Myeloperoxidase (MPO) Assay Procedure*



*Analyzer Dependent

For a list of validated parameters please contact Diazyme technical support at 858.455.4768 or email support@diazyme.com

1. Nilsson L et al. (1988) Activation of inflammatory system during cardiopulmonary bypass. Scand J Thorac Cardiovasc Surg. 22: 51-3
2. Heinecke JW et al. (1999) Mechanisms of oxidative damage by myeloperoxidase in atherosclerosis and other inflammatory disorders. J Lab Clin Med 133: 321-5
3. Podil'chak MD and Terletskaia LM (1988) Clinical value of determining myeloperoxidase and alkaline phosphatase activity of the leukocytes in patients with suppurative inflammatory processes. Klin Khir 59-60
4. G. Azzimondi G Re. et al (1997) Plasma lipoperoxidative markers in ischaemic stroke suggest brain embolism. European Journal of Emergency Medicine 4, 5-9
5. Luigi M. B. et al. (1996) Intracellular Neutrophil Myeloperoxidase is reduced in unstable angina and acute myocardial infarction, but its reduction is not related ischemia. JACC Vol. 27, No.3: 611-6.
6. Jessie Shih et al. (2008) Effect of collection tube type and preanalytical handling on myeloperoxidase concentrations. Clin. Chem. 54:6 1076-1079.
7. Podrez EA, Schmitt D, Hoff HF et al.: Myeloperoxidase-generated reactive nitrogen species convert LDL into an atherogenic form in vitro. J. Clin. Invest. 103, 1547-1560 (1999).
8. Naruko T, Ueda M, Haze K et al.: Neutrophil infiltration of culprit lesions in acute coronary syndromes. Circulation 106, 2894-2900 (2002).
9. Buffon A, Biasucci LM, Liuzzo G et al.: Widespread coronary inflammation in unstable angina. N. Engl. J. Med. 347, 5-12 (2002).
10. Sugiyama S, Okada Y, Sukhova GK et al.: Macrophage myeloperoxidase regulation by granulocyte macrophage colony-stimulating factor in human atherosclerosis and implications in acute coronary syndromes. Am. J. Pathol. 158, 879-891 (2001).

ASSAY PRECISION

The simple precision of the Diazyme MPO Immunoassay was evaluated. In the study, two levels of MPO controls containing 534 pmol/L (77ng/mL) and 3824 pmol/L (551 ng/mL) MPO respectively were tested with 15 duplicates in one run.

	Level 1: 534 pmol/L	Level 2: 3824 pmol/L
Number of Data Points	15	15
Mean (U/L)	534	3824
SD (U/L)	15	158
CV (%)	2.7%	4.1%

ASSAY INTERFERENCE

The substances normally present in the plasma were tested. Less than 10% deviation was produced when tested up to the concentrations shown below:

Ascorbic Acid:	10 mM
Bilirubin, free:	40 mg/dL
Bilirubin, conjugated:	40 mg/dL
Hemoglobin:	200 mg/dL
Triglyceride:	270 mg/dL
Rheumatoid:	75 IU/mL

DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064
PO Box 85608, San Diego, CA 92186
Tel: 858-455-4768 888-DIAZYME

www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park
Shanghai, 201203, People's Republic of China
Tel: 086-21-51320668 Fax: 086-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

