

# ***N-ACETYL- $\beta$ -D-GLUCOSAMINIDASE (NAG) ASSAY***

***Three Vial Liquid Stable***

Diazyme's NAG Assay is a cost effective test that uses a "highly sophisticated and powerful colorimetric substrate" that is not affected by urine color.<sup>1</sup> Each kit is supplied with a calibrator set for added convenience and controls are available separately. The test offers a wide range of instrument parameters for facilitating and simplifying implementation in the laboratory. Increased NAG levels in urine are usually an early indication of renal disease and can serve as a valuable renal monitoring test in disorders such as nephritic syndrome, glomerulonephritis, drug abuse associated nephrotoxicity, diabetes-associated nephropathy, hypertension and urinary tract infections.

## ***DIAZYME NAG ASSAY ADVANTAGES***

- Fast test results (under 5.5 minutes) for a rapid turnaround time
- Liquid stable reagent, calibrator and controls are offered separately for added convenience
- Wide range of instrument parameters available for facilitating and simplifying implementation
- Liquid stable format requires no reagent preparation, saving time and reducing sample handling

## ***REGULATORY STATUS***

510(k) Exempt



## ***AVAILABLE INSTRUMENT SPECIFIC PACKAGING***


- Roche
- Hitachi



## ASSAY SPECIFICATIONS

<b>Method</b>	Colorimetric (Enzymatic cleavage of a colorimetric substrate)
<b>Sample Type &amp; Volume</b>	• Urine Sample Volume 10 µL
<b>Linear Range</b>	Up to 200 U/L
<b>LOQ</b>	1.64 U/L
<b>Calibration Levels</b>	1-Point Calibration
<b>Reagent On-Board Stability</b>	Unopened: 24 months Opened: 1 month when stored at 2-8°C

### NAG Assay Procedure\*



#### \*Analyzer Dependent

Parameter questions for NAG Assay should be addressed to Diazyme technical support. Please call 858.455.4768 or email [support@diazyme.com](mailto:support@diazyme.com)

## ASSAY PRECISION

In the study, two levels of NAG controls and one NAG urine sample containing 40.9 U/L, 124.0 U/L and 9.64 U/L NAG respectively were tested on a Hitachi 917 in one run with 20 in replicates.

### Within-Run Precision:

	Sample 1	Sample 2	Sample 3
<b>N</b>	20	20	20
<b>Mean</b>	38.99	119.71	9.68
<b>Std. Dev.</b>	0.39	1.16	0.41
<b>CV %</b>	0.99%	0.97%	4.23%

## ASSAY INTERFERENCE

The common urine interfering substances triglyceride, ascorbic acid, free bilirubin, and conjugated bilirubin showed no significant interference ( $\geq 10\%$ ) up to the concentrations summarized below.

Triglyceride:	1000 mg/dL
Ascorbic Acid:	0.500 mg/dL
Bilirubin:	5 mg/dL
Bilirubin Conjugated:	5 mg/dL

## DIAZYME LABORATORIES

12889 Gregg Court, Poway, CA 92064  
PO Box 85608, San Diego, CA 92186  
Tel: 858-455-4768 888-DIAZYME

[www.diazyme.com](http://www.diazyme.com) [sales@diazyme.com](mailto:sales@diazyme.com)

## DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Deutschland  
Tel. +49 (0) 351 886 3300 Fax +49 (0) 351 886 3366  
[sales@diazyme.de](mailto:sales@diazyme.de)

## SHANGHAI DIAZYME CO., LTD.

Room 201,1011 Halei Road, Zhangjiang Hi-tech Park  
Shanghai, 201203, People's Republic of China  
Tel: 086-21-51320668 Fax: 086-21-51320663  
[www.lanyuanbio.com](http://www.lanyuanbio.com) [service@lanyuanbio.com](mailto:service@lanyuanbio.com)

