

TruQuick™ cTnI 10T

A rapid test for the diagnosis of myocardial infarction (MI) to detect cardiac Troponin I (cTnI) qualitatively in whole blood, serum or plasma.

REF TQ3410

IVD

Rx Only

INTENDED USE

TruQuick cTnI is a rapid chromatographic immunoassay for the qualitative detection of human cardiac Troponin I in whole blood, serum or plasma as an aid in the diagnosis of myocardial infarction (MI).

SUMMARY AND EXPLANATION OF THE TEST

Cardiac Troponin I (cTnI) is a protein found in cardiac muscle with a molecular weight of 22.5 kDa.¹ Troponin I is part of a three subunit complex comprising of Troponin T and Troponin C. Along with troponomyosin, this structural complex forms the main component that regulates the calcium sensitive ATPase activity of actomyosin in striated skeletal and cardiac muscle.² After cardiac injury occurs, Troponin I is released into the blood 4-6 hours after the onset of pain. The release pattern of cTnI is similar to CK-MB, but while CK-MB levels return to normal after 72 hours, Troponin I remains elevated for 6-10 days, thus providing for a longer window of detection for cardiac injury. The high specificity of cTnI measurements for the identification of myocardial damage has been demonstrated in conditions such as the perioperative period, after marathon runs, and blunt chest trauma.³ cTnI release has also been documented in cardiac conditions other than acute myocardial infarction (AMI) such as unstable angina, congestive heart failure, and ischemic damage due to coronary artery bypass surgery.⁴ Because of its high specificity and sensitivity in the myocardial tissue, Troponin I has recently become the most preferred biomarker for myocardial infarction.⁵

TruQuick cTnI is a simple test that utilizes a combination of anti-cTnI antibody-coated particles and capture reagent to detect cTnI in whole blood, serum or plasma. The minimum detection level is 0.5 ng/mL.

BIOLOGICAL PRINCIPLES

TruQuick cTnI is a qualitative, membrane based immunoassay for the detection of cardiac Troponin I (cTnI) in whole blood, serum or plasma. In this test procedure, capture reagent is immobilized in the test line region of the test. After specimen is added to the specimen area of the cassette, it reacts with anti-cTnI antibody-coated colloid gold particles in the test. This mixture migrates chromatographically along the length of the test and interacts with the immobilized capture reagent. The test format can detect cardiac Troponin I (cTnI) in specimens. If the specimen contains cardiac Troponin I (cTnI), a colored line will appear in the test line region, indicating a positive result. If the specimen does not contain cardiac Troponin I (cTnI), a colored line will not appear in this region, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Test Cassettes: The test contains anti-cTnI antibody-coated colloid gold particles and capture reagent coated on the membrane. Each cassette is packaged in a foil pouch.
- Buffer: A buffered solution containing Proclin 300 as a preservative. The Buffer is supplied in a dropper vial ready for use.
- Droppers
- Package insert

MATERIALS NOT PROVIDED

- Specimen Collection Containers
- Centrifuge
- Timer

For fingerstick whole blood

- Lancets
- Heparinized capillary tubes and dispensing bulb

PRECAUTIONS

1. For in vitro diagnostic use only. Do not use after expiration date.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Do not use test if pouch is damaged.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
5. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
6. The used test should be discarded according to local regulations.
7. Humidity and temperature can adversely affect results.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazard and Precautionary Statements.

SHELF LIFE AND STORAGE

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30 C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

SPECIMEN COLLECTION AND PREPARATION


1. TruQuick cTnI can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
2. To collect **Fingerstick Whole Blood specimens:**
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.

- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
 - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube:**
- Touch the end of the capillary tube to the blood until filled to approximately 75 μ L. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen area of the Test Cassette.
- Add the Fingerstick Whole Blood specimen to the test by using **hanging drops:**
- Position the patient's finger so that the drop of blood is just above the specimen area of the Test Cassette.
 - Allow three hanging drops of fingerstick whole blood to fall into the center of the specimen area on the Test Cassette, or move the patient's finger so that the hanging drop touches the center of the specimen area. Avoid touching the finger directly to the specimen area.
3. Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear non-hemolyzed specimens.
 4. Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8 C for up to three days. For long term storage, specimens should be kept below -20 C. Whole blood collected by venipuncture should be stored at 2-8 C if the test is to be run within one day of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
 5. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
 6. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

TEST PROCEDURE

Allow the test, specimen, Buffer and/or controls to reach room temperature (15-30 C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the Test Cassette from the sealed pouch and use it as soon as possible. Place the cassette on a clean and level surface.
2. For **Serum or Plasma** specimen:
 - Hold the dropper vertically and transfer **2 drops of serum or plasma (approximately 50 μ L)** to the specimen area, then **add 1 drop of Buffer (approximately 40 μ L)**, and start the timer. See illustration below.
- For **Venipuncture Whole Blood** specimen:
 - Hold the dropper vertically and transfer **3 drops of whole blood (approximately 75 μ L)** to the specimen area, then **add 1 drop of Buffer (approximately 40 μ L)**, and start the timer. See illustration below.
- For **Fingerstick Whole Blood** specimen:
 - To use a capillary tube: Fill the capillary tube and transfer **approximately 75 μ L of fingerstick whole blood specimen** to the specimen area of Test Cassette, then **add 1 drop of Buffer (approximately 40 μ L)** and start the timer. See illustration below.
 - To use hanging drops: Allow **3 hanging drops of fingerstick whole blood specimen (approximately 75 μ L)** to fall into the specimen area of Test Cassette, then **add 1 drop of Buffer (approximately 40 μ L)** and start the timer. See illustration below.
3. Wait for the colored line(s) to appear. Read results at 10 minutes. Do not interpret the result after 20 minutes.



INTERPRETATION OF RESULTS

(Please refer to the illustration above.)

POSITIVE: * **Two lines appear.** One colored line should be in the control line region (C) and another colored line should be in the test line region (T).

***NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of cardiac Troponin I (cTnI) present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

NEGATIVE: **One colored line appears in the control line region (C).** No line appears in the test line region (T).

INVALID: **Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact Meridian's Technical Services Department at 800-343-3858 or your local distributor.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

A procedural control is included in the test. A colored line appearing in the control line region(C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

External controls are not supplied with this kit; however, it is recommended that positive and negative external controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

TruQuick cTnI has been compared with a leading commercial cTnI EIA test, demonstrating an overall correlation of 99.1%.

LIMITATIONS OF THE PROCEDURE

1. TruQuick cTnI is for in vitro diagnostic use only. This test should be used for the detection of Troponin I in whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in cTnI can be determined by this qualitative test.
2. TruQuick cTnI will only indicate the qualitative level of cTnI in the specimen and should not be used as the sole criteria for the diagnosis of myocardial infarction.
3. TruQuick cTnI cannot detect less than 0.5 ng/mL of cTnI in specimens. A negative result at any time does not preclude the possibility of myocardial infarction.
4. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
5. Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor (RF) may affect expected results. Even if the test results are positive, further clinical evaluation should be considered with other clinical information available to the physician.
6. There is a slight possibility that some whole blood specimens with very high viscosity or which have been stored for more than one day may not run properly on the Test Cassette. Repeat the test with a serum or plasma specimen from the same patient using a new test cassette.

SPECIFIC PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

TruQuick cTnI was evaluated with a leading commercial cTnI EIA test using clinical specimens. The results show that the sensitivity of TruQuick cTnI is 99.4% and the specificity is 99.0% relative to the leading EIA test.

Method	EIA		Total Result		
	Results	Positive		Negative	
TruQuick cTnI	Positive	172	5	177	
	Negative	1	472	473	
Total Result			173	477	650

Sensitivity: 172/173 = 99.4% (95% CI*: 96.8% ~ 99.9%);

Specificity: 472/477 = 99.0% (95% CI*: 97.6% ~ 99.7%);

Correlation: (172+472)/(172+1+5+472) = 99.1%(95% CI*: 98.0% ~ 99.7%). *Confidence Intervals

REPRODUCIBILITY

Intra-Assay Precision

Within-run precision was determined by using 15 replicates of five specimens: a negative, cTnI 1.0 ng/mL positive, cTnI 5.0 ng/mL positive, cTnI 10 ng/mL positive and cTnI 40 ng/mL positive. The negative, cTnI 1.0 ng/mL positive, cTnI 5.0 ng/mL positive, cTnI 10 ng/mL positive and cTnI 40 ng/mL positive values were correctly identified > 99% of the time.

Inter-Assay Precision

Between-run precision was determined by 15 independent assays on the same five specimens: a negative, cTnI 1.0 ng/mL positive, cTnI 5.0 ng/mL positive, cTnI 10 ng/mL positive and cTnI 40 ng/mL positive specimens. Three different lots of TruQuick cTnI were tested over a three day period using negative, cTnI 1.0 ng/mL positive, cTnI 5.0 ng/mL positive, cTnI 10 ng/mL positive and cTnI 40 ng/mL positive specimens. The specimens were correctly identified > 99% of the time.

CROSSREACTIVITY

TruQuick cTnI was tested using samples from patients diagnosed with 10,000 ng/mL Skeletal Troponin I, 2,000 ng/mL Troponin T, 20,000 ng/mL Cardiac Myosin, HBsAg, HBeAg, HBeAb, HBeAb, syphilis, anti-HIV, anti-*H. pylori*, Mononucleosis, anti-CMV, anti-Rubella and anti-Toxoplasmosis. The results showed no crossreactivity.

ANALYTICAL SENSITIVITY

The limit of detection for TruQuick cTnI is 0.5 mg/mL of cardiac troponin I antigen in plasma.

TESTS FOR INTERFERING SUBSTANCES

The following potentially interfering substances were added to cTnI negative and positive specimens.

Acetaminophen 20 mg/dL	Caffeine 20 mg/dL
Acetylsalicylic Acid 20 mg/dL	Genistic Acid 20 mg/dL
Ascorbic Acid 20 mg/dL	Albumin 10,500 mg/dL
Creatine 200 mg/dL	Hemoglobin 1,000 mg/dL
Bilirubin 1,000 mg/dL	Oxalic Acid 600 mg/dL
Cholesterol 800 mg/dL	Triglycerides 1,600 mg/dL


None of the substances interfered in the assay at the concentration tested.

REFERENCES

1. Adams JE, et al. Biochemical markers of myocardial injury. Is MB creatine kinase the choice for the 1990s? Immunoassay Circulation. AMA. 1993;88:750-763.
2. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. J Biol Chem. 1991;266:966.
3. Adams JE, Sicard GA, Allen BT, Bridwell KH, Lenke LG, Davila-Roman VG, et al. Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. N Engl J Med. 1994;330:670-674.
4. Hossein-Nia M, et al. Cardiac troponin I release in heart transplantation. Ann Thorac Surg. 1996;61:227.
5. Alpert JS, et al. Myocardial infarction redefined, Joint European Society of Cardiology American College of Cardiology: J Am Coll Cardio. 2000;36(3):959.

SNTQ3410












REV. 05/17

 Manufactured By	Meridian Bioscience, Inc. Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 USA Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
--	--

SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols

	Use By	CONTROL +	Positive control
LOT	Batch Code	CONTROL -	Negative control
IVD	In vitro diagnostic medical device	EC REP	Authorized Representative in the European Community
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent
REF	Catalogue number		Do not freeze
	Consult Instructions for Use	BUF RXN	Reaction Buffer
	Manufacturer		For IVD Performance Evaluation Only
	Contains sufficient for $\lt;-\gt;$ tests	SOLN STOP	Stopping Solution
	Temperature limitation	CONJ ENZ	Enzyme Conjugate
SN	Serial number	CONTROL	Assay Control
TEST	Test Device	REAG	Reagent
	Date of manufacture	BUF WASH	Wash Buffer
BUF	Buffer		Warning
CONJ	Conjugate	DIL SPE	Specimen Diluent (or Sample Diluent)
SUBS	Substrate	BUF WASH 20X	Wash Buffer Concentration: 20X
RUO	Research Use Only	DET REAG	Detection Reagent
IUO	Investigational Use Only	R_x Only	Prescription Use Only
	Do not use if package is damaged		

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.