

## TruQuick™ MYO 10T

A rapid test to detect Myoglobin qualitatively in whole blood, serum or plasma.

REF TQ3810

IVD

Rx Only

### INTENDED USE

TruQuick MYO is a rapid chromatographic immunoassay for the qualitative detection of human Myoglobin in whole blood, serum or plasma as an aid in the diagnosis of myocardial infarction (MI).

### SUMMARY AND EXPLANATION OF THE TEST

Myoglobin (MYO) is a heme-protein normally found in skeletal and cardiac muscle with a molecular weight of 17.8 kDa. It constitutes about two percent of total muscle protein and is responsible with transporting oxygen within the muscle cells.<sup>1</sup> When the muscle cells are damaged, Myoglobin is released to the blood rapidly due to its relatively small size. Following the death of tissue associated with MI, Myoglobin is one of the first markers to rise above normal levels. The level of Myoglobin increases measurably above baseline within two to four hours post-infarct, peaking at nine to 12 hours and returning to baseline within 24-36 hours.<sup>2, 3</sup> A number of reports suggest the measurement of Myoglobin as a diagnostic aid in confirming the absence of myocardial infarction with negative predictive values of up to 100% reported at certain time periods after onset of symptoms.<sup>4</sup>

TruQuick MYO is a simple test that utilizes a combination of anti-Myoglobin antibody-coated particles and capture reagents to qualitatively detect Myoglobin in whole blood, serum or plasma. The minimum detection level is 50 ng/mL.

### BIOLOGICAL PRINCIPLES

TruQuick MYO is a qualitative, membrane based immunoassay for the detection of Myoglobin in whole blood, serum or plasma. The membrane is pre-coated with specific capture antibodies in the test line regions of the test. During testing, the whole blood, serum or plasma specimen reacts with the particle coated with specific antibodies. The mixture migrates upward on the membrane chromatographically by capillary action to react with specific capture reagents on the membrane to generate a colored line. The presence of this colored line in the specific test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Test Cassettes: The test contains anti-Myoglobin antibody-conjugated colloid gold particles and capture reagents coated on the membrane.
- Buffer: A buffered solution containing Proclin 300 as a preservative. The Buffer is supplied in a dropper vial ready to use.
- Droppers
- Package insert

### MATERIALS NOT PROVIDED

- Specimen Collection Containers
- Centrifuge
- Timer

### For fingerstick whole blood

- Lancets
- Heparinized capillary tubes and dispensing bulb

### PRECAUTIONS

1. For in vitro diagnostic use only. Do not use after expiration date.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Do not use test if pouch is damaged.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
5. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
6. The used test should be discarded according to local regulations.
7. Humidity and temperature can adversely affect results.

### HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazard and Precautionary Statements.

### SHELF LIFE AND STORAGE

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30 C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

### SPECIMEN COLLECTION AND PREPARATION

1. TruQuick MYO can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
2. To collect **Fingerstick Whole Blood specimens**:
  - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
  - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
  - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
  - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.

Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube**:

- Touch the end of the capillary tube to the blood until filled to approximately 75 µL. Avoid air bubbles.
- Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen area of the test cassette.


Add the Fingerstick Whole Blood specimen to the test by using **hanging drops**:

- Position the patient's finger so that the drop of blood is just above the specimen area of the Test Cassette.
  - Allow three hanging drops of fingerstick whole blood to fall into the center of the specimen area on the Test Cassette, or move the patient's finger so that the hanging drop touches the center of the specimen area. Avoid touching the finger directly to the specimen area.
3. Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear nonhemolyzed specimens.
  4. Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8 C for up to three days. For long term storage, specimens should be kept below -20 C. Whole blood collected by venipuncture should be stored at 2-8 C if the test is to be run within two days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
  5. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
  6. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

### TEST PROCEDURE

Allow the test, specimen, Buffer and/or controls to reach room temperature (15-30 C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the Test Cassette from the sealed pouch and use it within one hour. Place the cassette on a clean and level surface.
  - A. For **Serum or Plasma** specimen:
    - Hold the dropper vertically and transfer **2 drops of serum or plasma (approximately 50 µL)** to the specimen area, then **add 1 drop of Buffer (approximately 40 µL)**, and start the timer. See illustration below.
  - A. For **Venipuncture Whole Blood** specimen:
    - Hold the dropper vertically and transfer **3 drops of whole blood (approximately 75 µL)** to the specimen area, then **add 1 drop of Buffer (approximately 40 µL)**, and start the timer. See illustration below.
  - A. For **Fingerstick Whole Blood** specimen:
    - To use a capillary tube: Fill the capillary tube **and transfer approximately 75 µL of fingerstick whole blood specimen** to the specimen area of Test Cassette, then **add 1 drop of Buffer (approximately 40 µL)** and start the timer. See illustration below.
    - To use hanging drops: Allow **3 hanging drops of fingerstick whole blood specimen (approximately 75 µL)** to fall into the specimen area of Test Cassette, then **add 1 drop of Buffer (approximately 40 µL)** and start the timer. See illustration below.
3. Wait for the colored line(s) to appear. Read results at 10 minutes. Do not interpret the result after 20 minutes.



### INTERPRETATION OF RESULTS

(Please refer to the illustration above.)

**POSITIVE:** A colored line in the control line region (C) and the presence of one colored line in the test line regions indicates a positive result. This indicates that the concentration of Myoglobin is above the minimum detection level.

**\*NOTE:** The intensity of the color in the test line region will vary depending on the concentration of Myoglobin, present in the specimen. Therefore, any shade of color in the test line region should be considered positive.

**NEGATIVE:** One colored line appears in the control line region (C). No line appears in the test line region (T). This indicates that the concentration of Myoglobin is below the minimum detection level.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact Meridian's Technical Services Department at 800-343-3858 or your local distributor.

### QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

Internal procedural controls are included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

External controls are not supplied with this kit; however, it is recommended that positive and negative external controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

### EXPECTED VALUES

TruQuick MYO has been compared with a leading commercial Myoglobin EIA test, demonstrating an overall correlation of 97.5% with Myoglobin EIA.

### LIMITATIONS OF THE PROCEDURE

1. TruQuick MYO is for in vitro diagnostic use only. This test should be used for the detection of Myoglobin in whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in Myoglobin can be determined by this qualitative test.
2. TruQuick MYO will only indicate the qualitative level of Myoglobin in the specimen and should not be used as the sole criteria for the diagnosis of myocardial infarction.
3. TruQuick MYO cannot detect less than 50 ng/mL Myoglobin in specimens. A negative result at any time does not preclude the possibility of myocardial infarction.
4. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
5. Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor (RF) may affect expected results. Even if the test results are positive, further clinical evaluation should be considered with other clinical information available to the physician.
6. There is a slight possibility that some whole blood specimens with very high viscosity or which have been stored for more than two days may not run properly on the Test Cassette. Repeat the test with a serum or plasma specimen from the same patient using a new Test Cassette.

### SPECIFIC PERFORMANCE CHARACTERISTICS

#### Sensitivity and Specificity

TruQuick MYO has been evaluated with a leading commercial Myoglobin EIA test using clinical specimens. The results show that relative to leading EIA tests, TruQuick MYO shows > 99.9% sensitivity and 97.2% specificity for Myoglobin.

#### TruQuick MYO vs. EIA

Method	EIA		Total Result
	Positive	Negative	
TruQuick MYO	54	11	65
	0	379	379
Total Result	54	390	444

Sensitivity: 54/54 = > 99.9% (95% CI\*: 94.6%~100.0%);

Specificity: 379/390 = 97.2% (95% CI\*: 95.0%~98.6%);

Correlation: (54+379)/(54+11+379) = 97.5% (95% CI\*: 95.6%~98.8%). \*Confidence Intervals

### REPRODUCIBILITY

#### Intra-Assay Precision

Within-run precision was determined by using 15 replicates of serum and plasma specimens of the following: Myoglobin specimen levels at 0 ng/mL, 50 ng/mL, 100 ng/mL, 200 ng/mL and 400 ng/mL. The specimens were correctly identified > 99% of the time.

#### Inter-Assay Precision

Between-run precision was determined by three independent assays on serum and plasma for the same five specimens: 0 ng/mL, 50 ng/mL, 100 ng/mL, 200 ng/mL and 400 ng/mL of Myoglobin. Three different lots of TruQuick MYO were tested using these specimens. The specimens were correctly identified > 99% of the time.

### CROSSREACTIVITY

TruQuick MYO was tested with samples containing HBsAg, HBsAb, HBeAg, HBeAb, HBeAb, HBeAb, syphilis, anti-HIV, anti-*H. pylori*, Mononucleosis, anti-CMV, anti-Rubella and anti-Toxoplasmosis. The samples showed no crossreactivity.

### TESTS FOR INTERFERING SUBSTANCES

The following potentially interfering substances were added to Myoglobin negative and positive specimens.

Acetaminophen 20 mg/dL	Caffeine 20 mg/dL
Acetylsalicylic acid 20 mg/dL	Genistic acid 20 mg/dL
Ascorbic acid 20 mg/dL	Albumin 10,500 mg/dL
Creatine 200 mg/dL	Hemoglobin 1,000 mg/dL
Bilirubin 1,000 mg/dL	Oxalic acid 600 mg/dL
Cholesterol 800 mg/dL	Triglycerides 1,600 mg/dL

None of the substances interfered in the assay at the concentration tested.

### REFERENCES

1. Wong SS. Strategic utilization of cardiac markers for diagnosis of acute myocardial infarction. Ann Clin Lab Sci. 1996;26:301-12.
2. Kagen LJ. Myoglobin methods and diagnostic uses. CRC Crit. Rev. Clin. Lab. Sci. 1978;2:273.
3. Chapelle JP, et al. Serum myoglobin determinations in the assessment of acute myocardial infarction. Eur. Heart J. 1982;3:122.

4. Hamfelt A, et al. Use of biochemical tests for myocardial infarction in the county of Vasternorrland, a clinical chemistry routine for the diagnosis of myocardial infarction. Scand J Clin Lab Invest Suppl. 1990;200:20.



SNTQ3810












REV. 06/17

 Manufactured By	Meridian Bioscience, Inc. Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 USA Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. l. Via dell' Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

**SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols**

	Use By	<b>CONTROL +</b>	Positive control
<b>LOT</b>	Batch Code	<b>CONTROL -</b>	Negative control
<b>IVD</b>	In vitro diagnostic medical device	<b>EC REP</b>	Authorized Representative in the European Community
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices	<b>SMP   PREP   DIL   SPE</b>	Sample Preparation Apparatus containing Sample Diluent
<b>REF</b>	Catalogue number		Do not freeze
	Consult Instructions for Use	<b>BUF RXN</b>	Reaction Buffer
	Manufacturer		For IVD Performance Evaluation Only
	Contains sufficient for <math>\leq 10</math> tests	<b>SOLN STOP</b>	Stopping Solution
	Temperature limitation	<b>CONJ   ENZ</b>	Enzyme Conjugate
<b>SN</b>	Serial number	<b>CONTROL</b>	Assay Control
<b>TEST</b>	Test Device	<b>REAG</b>	Reagent
	Date of manufacture	<b>BUF   WASH</b>	Wash Buffer
<b>BUF</b>	Buffer		Warning
<b>CONJ</b>	Conjugate	<b>DIL   SPE</b>	Specimen Diluent (or Sample Diluent)
<b>SUBS</b>	Substrate	<b>BUF   WASH   20X</b>	Wash Buffer Concentration: 20X
<b>RUO</b>	Research Use Only	<b>DET   REAG</b>	Detection Reagent
<b>IUO</b>	Investigational Use Only	<b>R. Only</b>	Prescription Use Only
	Do not use if package is damaged		

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.