# V E T C H R O M A T 4

Note: Please refer to the table below to identify various symbols

	Sufficient <n> tests
	Read instruction for use
	Use by date
	Batch code
	Catalog number
	Caution
	Manufacturer
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse

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




- Turn the power on and Please follow procedure on the Screening of Vet chroma™ reader. This picture is for reference only.
  - STEP 1 : Insert the cartridge.
  - STEP 2 : Enter the required information.
  - STEP 3 & STEP 4 : Follow the picture below.

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# Vet chroma™ T4

## INTENDED USE

Vet chroma™ T4 is a Fluorescence Immunoassay (FIA) for the quantitative determination of total thyroxine (T4) in canine, feline serum/plasma. It is useful as an aid in management and monitoring of thyroid disorder.

## INTRODUCTION

Thyroxine (T4) is one of two major hormones produced by the thyroid gland (the other is called triiodothyronine, or T3). T4 is a useful marker for the diagnosis of hypothyroidism and hyperthyroidism in dog and cat.

## WARNINGS AND PRECAUTIONS

- For in Vitro use only.
- Carefully follow the instructions and procedures described in this ‘‘Instruction for use’’.
  - Sample
    - The sample type: Serum or plasma (Heparin) of Dog and Cat.
    - Recommend the fresh samples and avoid direct sunlight.
    - If the test could not be performed within 24 hours, serum or plasma should be immediately frozen below −20℃. The freezing storage of sample up to 3 months does not affect the quality of results.
    - Once the sample was frozen, it should be thawed one time and only for test, because repeated freezing and thawing can result in the changed test values.
  - Sample with severe hemolytic and hyperlipidemia cannot be used and should be recollected.
  - Do not cross use collection tool of sample between specimen and specimen. Because it avoid contaminating each other. It must be use the clean and new tools at the time of use.

- Reagents and cartridges
  - Just before use, allow the cartridge, detection buffer and sample to be at room temperature for approximately 30 minutes.
  - The cartridge and detection buffer can be stored at room temperature for up to 12 hours.
  - Test components that were left at room temperature over 12 hours may affect test results.
  - Test components that were left at room temperature frequently will expire sooner.
  - Used detection buffer tubes, pipette tips and cartridges should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
  - Do not reuse. A detection buffer and capillary cap should be used for processing one sample only. So should a cartridge.
  - Lot numbers of all the test components (cartridge, Lot chip and detection buffer) must match each other.
  - Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield misleading of test result(s).
  - The cartridge should remain sealed in its original pouch before use. Do not use the cartridge, if is damaged or already opened.
  - The detection buffer contains sodium azide for preservatives. If you exposure to skin or eye with reagent, please rinse immediately with clean water and ask for medical advice.

- Reader
  - Vet chroma™ T4 as well as the instrument for Vet chroma™ tests should be used away from vibration and/or magnetic field. During normal usage, it can be noted that instrument for Vet chroma™ tests may produce minor vibration.

## LIMITATION OF THE TEST SYSTEM

- The test may yield false positive/negative result(s) due to the flow things.
  - The cross–reactions and/or non–specific adhesion of certain sample components to the capture/de–tector antibodies.
  - The non–responsiveness of the antigen to the antibodies is most common where the epitope is masked by some unknown components, so as not to be detected or captured by the antibodies.
  - The instability or degradation of the antigen with time and/or temperature may cause the false negative as it makes antigen unrecognizable by the antibodies.
  - Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.

- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

## DEMAND COMPONENT

- Instrument for Vet chroma™ T4 tests
  - Vet chroma™ Reader FPRR025 (Please refer to the ‘‘Instrument for Vet chroma™ tests operation manual’’ for complete information and operating instructions)

## MATERIALS SUPPLIED

- REF) CFPA–5–1
- Cartridge Box:
  - Cartridges
  - Detection Buffer tubes
  - Lot Chip
  - Capillary caps
  - Instruction for Use

### STORAGE AND STABILITY

- The cartridge and detection buffer is stable for 24 months (while sealed in an aluminum foil pouch) if stored at 2–8 ℃.
- After the cartridge pouch is opened, the test should be performed immediately.

### QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and can be performed at regular intervals.
- Quality control tests could be used for re–calibration of product.
- Control materials are not provided with Vet chroma™ T4. For more information regarding obtaining the control materials, contact ANIVET Inc.’ s Sales Division for assistance. (Please refer to the instruction for use of control material.)

## TEST SETUP

- Check the contents of Vet chroma™ T4: Sealed Cartridge, Detection Buffer Tubes and Lot Chip and capillary caps.
- Ensure that the lot number of the cartridge matches that of the Lot chip as well as the detection buffer.
- Keep the sealed cartridge (if stored in refrigerator) and the detection buffer tube at room temperature for at least 30 minutes just prior to the test.
- Place the cartridge on a clean, dust–free and flat surface.

## TEST PROCEDURE

- Please follow procedure on the Screening of Vet chroma™ reader.
  - Make a puncture on the top of the detection buffer tube by inserting an empty capillary cap.
  - Transfer 50µl of sample (Canine, Feline serum/plasma/control) using a micropipette to a tube containing the Detection buffer tube (punched tube).
  - Assemble the capillary cap and the tube into one. Shake the 5 times or more until the sample out of the capillary cap by inversion. Incubate the sample mixture at room temperature for 5 minutes.
  - After 5minutes, remove the cap off the top of assembled tube. Discard 2 drops of reagent onto the paper towel before applying to the cartridge. Load only 2 drops of the mixture onto the sample well of the cartridge. And Leave the Cartridge on the Vet chroma™ holder.
  - Instrument for Vet chroma™ tests will start automatically scanning the sample–loaded cartridge after 10minutes.
  - Read the test result on the display screen of the instrument for Vet chroma™ tests.

## INTERPRETATION OF TEST RESULT

- Instrument for Vet chroma™ tests calculates the test result automatically and displays T4 concentration of the test sample in terms of nmol/L and µg/dL
- Working range & Reference Values:
  - Working range: 8–90 nmol/L (0.62 – 6.99 µg/dL)
  - Unit conversion factor: 1 µg/dL = 12.87 nmol/L
  - Reference values:

		nmol/L	µg/dt
Canine(Dog)	Low	<15	<1.17
	Borderline	15-23	1.17-1.79
	Normal	23-50	1.79-3.89
	High	>50	>3.89
Feline(Cat)	Low	<15	<1.17
	Normal	15-40	1.17-3.11
	Borderline	40-50	3.11-3.89
	High	>50	>3.89

- These recommendations are for guidance only.

- The aim of this test is to measure circulating total T4 concentrations in dogs and cats.
- Interpretation of the result by the veterinarian should always. Take into account the history, clinical examination, and any further diagnostic test results, as no diagnostic method is 100% accurate.

### Manufactured by

### ANIVET DIAGNOSTICS Incorporated

43, Geodudanji 1–gil, Dongnae–myeon, Chuncheon–si, Gang–won–do, 24398 Republic of Korea
Tel: +(82) –33–243–1411 / Fax: +(82) –33–243–9373

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