



Revised 16 Sept. 2010 rm (Vers. 5.1)

USA:RUO

INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of Hepatitis B Virus "e" Antigen and Antibody in human plasma and sera.

For Research Use only.

PRINCIPLE OF THE TEST

HBeAg:

HBeAg, if present in the sample, is captured by a specific monoclonal antibody, in the 1st incubation.

In the 2nd incubation, after washing, a tracer, composed of a mix of two specific anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP), is added to the microplate and binds to the captured HBeAg.

The concentration of the bound enzyme on the solid phase is proportional to the amount of HBeAg in the sample and its activity is detected by adding the chromogen/substrate in the 3^{rd} incubation.

The presence of HBeAg in the sample is determined by means of a cut-off value that allows for the semiquantitative detection of the antigen.

HBeAb:

Anti HBeAg antibodies, if present in the sample, compete with a recombinant HBeAg preparation for a fixed amount of an anti HBeAg antibody, coated on the microplate wells.

The competitive assay is carried out in two incubations, the first with the sample and recHBeAg, and the second with a tracer, composed of two anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP).

The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti HBeAg antibodies in the sample and its activity is detected by adding the chromogen/substrate in the third incubation.

The concentration of HBeAg specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti HBeAg antibodies.





Revised 16 Sept. 2010 rm (Vers. 5.1)

USA:RUO

COMPONENTS

The kit contains reagents for total 96 tests.

1. Microplate

n° 1 coated microplate

12 strips of 8 breakable wells coated with anti HBeAg specific monoclonal antibody, postcoated with bovine serum proteins and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at $2^{\circ}C - 8^{\circ}C$

2. Negative Control

1 x 2.0 ml/vial. Ready to use control. It contains bovine serum, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The negative control is colorless.

3. Antigen Positive Control

1 x 1.0 ml/vial. Ready to use control. It contains 2% bovine serum albumin, non infectious recombinant HBeAg, 100 mM tris buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The positive control is green color coded.

4. Antibody Positive Control

1 x 1.0 ml/vial. Ready to use control. It contains 2% bovine serum albumin, human anti HbeAg positive plasma at about 10 PEI U/ml, 100 mM tris buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The label is red colored. The positive control is yellow color coded.

5. Antigen Calibrator

n° 1 vial. Lyophilized calibrator for HBeAg.

To be dissolved with EIA grade water as reported in the label. It contains fetal bovine serum, non infectious recombinant HBeAg at 1 PEI U/ml +/-10%, 0.02% gentamicine sulphate and 0.1% Kathon GC as preservatives.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

6. Antibody Calibrator

n° 1 vial. Lyophilized calibrator for anti HBeAg a ntibody.

To be dissolved with EIA grade water as reported in the label.

It contains fetal bovine serum, positive plasma at 0.25 PEI U/ml +/-10%, 0.02% gentamicine sulphate and 0.1% Kathon GC as preservatives. The label is red colored.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

7. Wash buffer concentrate

1 x 60 ml/bottle. 20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.





Revised 16 Sept. 2010 rm (Vers. 5.1)

USA:RUO

8. Enzyme conjugate

1 x 16 ml/vial. Ready to use conjugate. It contains Horseradish peroxidase conjugated with a mix of monoclonal antibodies to HBeAg, 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives. The reagent is red color coded.

9. HBe Antigen

 1×10 ml/vial. Ready to use reagent. It contains recombinant HBeAg, fetal bovine serum, buffered solution pH 8.0+/-0.1, 0.1% Kathon GC and 0.09% sodium azide as preservatives. The reagent is blue color coded.

10. Chromogen/Substrate

1 x 16 ml/vial. Ready-to-use component. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H_2O_2 . *Note: To be stored protected from light as sensitive to strong illumination.*

11. Sulphuric Acid

1 x 15 ml/vial. It contains 0.3 M H₂SO₄ solution. Attention: Irritant (Xi R36/38; S2/26/30)

- **12. Plate sealing foils** n°2
- **13. Package insert** n°1

MATERIALS REQUIRED BUT NOT PROVIDED

- 1. Calibrated Micropipettes (150ul, 100ul and 50ul) and disposable plastic tips.
- 2. EIA grade water (double distilled or deionized, charcoal treated to remove oxidizing chemicals used as disinfectants).
- 3. Timer with 60 minute range or higher.
- 4. Absorbent paper tissues.
- 5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
- 6. Calibrated ELISA microwell reader with 450nm (reading) and if possible with 620-630nm (blanking) filters.
- 7. Calibrated ELISA microplate washer.
- 8. Vortex or similar mixing tools.

SPECIMEN: PREPARATION AND RECOMMENDATIONS

- 1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- 2. Avoid any addition of preservatives; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
- 3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results.
- 4. Hemolyzed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.





Revised 16 Sept. 2010 rm (Vers. 5.1)

- 5. Sera and plasma can be stored at $+2^{\circ}$ C 8° C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20° C for several months. Any frozen samples should not be freeze/thawed more than once as this may generate particles that could affect the test result.
- 6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8 μ filters to clean up the sample for testing.

PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

1. Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call the customer service. Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at $+2^{\circ}$ C - 8° C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Antigen Positive Control:

Ready to use. Mix well on vortex before use.

4. Antibody Positive Control:

Ready to use. Mix well on vortex before use.

5. Antigen Calibrator:

Add the volume of ELISA grade water, reported on the label, to the lyophilized powder; let fully dissolve and then gently mix on vortex.

Note: The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.

6. Antibody Calibrator:

Add the volume of ELISA grade water, reported on the label, to the lyophilized powder; let fully dissolve and then gently mix on vortex.

Note: The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.

7. Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200 ml and mixed gently end-over-end before use.

During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles. *Note: Once diluted, the wash solution is stable for 1 week at* $+2^{\circ}C - 8^{\circ}C$.

8. Enzyme conjugate:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

9. HBe Antigen:

Ready to use. Mix well on vortex before use.

USA:RUO





USA:RUO

DRG[®] HBe Ag/Ab ELISA (EIA-4162)

Revised 16 Sept. 2010 rm (Vers. 5.1)

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

10. Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

11. Sulphuric Acid:

Ready to use. Mix well on vortex before use. Attention: Irritant (Xi R36/38; S2/26/30) Legenda: R 36/38 = Irritating to eyes and skin.

S 2/26/30 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

- 1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of +2%.
- 2. The ELISA incubator has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- 3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests. 4-5 washing cycles (aspiration + dispensation of 350 µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
- 4. Incubation times have a tolerance of +5%.
- 5. The ELISA reader has to be equipped with a reading filter of 450nm and ideally with a second filter (620-630 nm) for blanking purposes Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
- 6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.





Revised 16 Sept. 2010 rm (Vers. 5.1)

USA:RUO

PRE ASSAY CONTROLS AND OPERATIONS

- 1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- 2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB+H₂O₂) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
- 3. Dilute all the content of the 20x concentrated Wash Solution as described above.
- 4. Dissolve the Calibrator as described above and gently mix.
- 5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- 6. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
- 7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- 8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- 9. Check that the micropipettes are set to the required volume.
- 10. Check that all the other equipment is available and ready to use.

In case of problems, do not proceed further with the test and advise the supervisor.

ASSAY PROCEDURE

The assay has to be performed according to the procedure given below, taking care to maintain the same incubation time for all the samples being tested.

HBe Antigen

- 1. Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- 2. Leave the A1 well empty for blanking purposes.
- 3. Pipette 100 µl of the Negative Control in triplicate, 100 µl of the Antigen Calibrator in duplicate and then 100 µl of the Antigen Positive Control in single.
- 4. Then dispense $100 \ \mu l$ of samples in the proper wells.
- 5. Check for the presence of samples in wells by naked eye (there is a marked colour difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0.100).
- 6. Incubate the microplate for $60 \text{ min at } +37^{\circ}\text{C}$.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

7. When the first incubation is finished, wash the microwells as previously described (section 9.3)





USA: RUO

DRG[®] HBe Ag/Ab ELISA (EIA-4162)

Revised 16 Sept. 2010 rm (Vers. 5.1)

8. Dispense 100 µl Enzyme Conjugate in all wells, except for A1, used for blanking operations.

Important note: Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.

- 9. Check that the reagent has been dispensed properly and then incubate the microplate for 60 min at +37°C.
- 10. When the second incubation is finished, wash the microwells as previously described (section 9.3)
- 11. Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

Important note: Do not expose to strong direct light as a high background might be generated.

- 12. Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive control and positive samples will turn from clear to blue.
- 13. Pipette 100 μl Sulphuric Acid into all the wells using the same pipetting sequence as in step 11. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
- 14. Measure the color intensity of the solution in each well, as described in section 9.5 using a 450 nm filter (reading) and if possible a 620-630nm filter (background subtraction), blanking the instrument on A1.

HBe Antibody

- 1. Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- 2. Leave the A1 well empty for blanking purposes.
- 3. Pipette 50 µl of the Negative Control in triplicate, 50 µl of the Antibody Calibrator in duplicate and then 50 µl of the Antibody Positive Control in single.
- 4. Then dispense 50 μ l of samples in the proper wells.
- 5. Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0.100).
- 6. Dispense then 50 μ l of HBe Antigen in all the wells, except for A1.
- 7. Incubate the microplate for 60 min at +37°C.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- 8. When the first incubation is finished, wash the microwells as previously described (section 9.3)
- 9. Finally proceed as described for the HBeAg assay from point 8 to the last one.

Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.

2. Reading should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.



Revised 16 Sept. 2010 rm (Vers. 5.1)

USA:RUO

ASSAY SCHEME

HBe antigen test	
Controls and calibrator	100 µl
Samples	100 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme Conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂ mix	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

HBe antibody test

Controls and calibrator	50 µl
Samples	50 µl
Neutralising antigen	50 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzymatic conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂ mixture	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm



Revised 16 Sept. 2010 rm (Vers. 5.1)

An example of dispensation scheme is reported below:

Microplate 1 2 3 5 4 6 7 8 9 10 11 12 S2 BLK Α NC **S**3 B С NC S4 NC S5 D CAL E S6 CAL **S**7 F G PC **S**8 **S**1 **S**9 Η Legenda: BLK = Blank NC = Negative Control CAL = Calibrator PC = Positive Control S = Sample

CALCULATION OF THE CUT-OFF

The results are calculated by means of a cut-off value determined with the following formula:

HBeAg:

NC + 0.100 = Cut-Off (Co)

The value found for the test is used for the interpretation of results as described in the next paragraph.

HBeAb:

(NC + PC) / 3 = Cut-Off (Co)

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Note: S = OD450nm of the sample , Co = cut-off value



USA:RUO

9



Revised 16 Sept. 2010 rm (Vers. 5.1)



USA: RUO

REFERENCES

- 1. Engvall E. and Perlmann P. J. Immunochemistry, 8, 871-874, 1971
- 2. Engvall E. and Perlmann P.. J.Immunol. 109, 129-135, 1971
- 3. Remington J.S. and Klein J.O.. In "Infectious diseases of the fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
- 4. Volk W.A. In "Essential of Medical Microbiology". 2nd ed. pp 729, G.B.Lippincott Company, Philadelphia, New York, S.Josè, Toronto.
- 5. Snydman D.R., Bryan J.A. and Dixon R.E. Ann.Int.Med., 83, pp 838, 1975.
- 6. Barker L.F., Gerety R.J., Lorenz D.E., Viral Hepatitis. 581-587, 1978.
- 7. Cossart Y.. Brit.Med.Bull. 28, pp 156, 1972
- 8. Lander J.J., Alter H. and Purcell R.. J.Immunol. 106, pp 1066, 1971
- 9. Mushawar I.K., Dienstag J.L., Polesky H.F. et al.. Ann.J.Clin.Pathol.. 76, pp 773, 1981.
- 10. Ling C.M., Mushawar I.K. et al.. Infection and Immunity, 24: 235, 1979.
- 11. Mushawar I.K., Overby L.R. et al.. J.Med.Virol..2: 77, 1978
- 12. Aldershville J., Frosner G.G. et al.. J.Med.Dis., 141: 293, 1980
- 13. Magnius L.O., Lindhom A. et al., J.Am.Med.Assoc., 231: 356, 1975
- 14. Krugman S., Overby L.R. et al.. N.Engl.J.Med.. 300: 101, 1979