SALIVA ELISAS

Speichel Analytik

DRG

Testosterone
Progesterone
Estradiol
17alpha-OH Progesterone
DHEA
Estriol
Cortisol
Lactate

DRG

DRG ELISAS

SALIVA Immunoassays Speichel Analytik

Sensitive ELISAS
for detection of free steroids
and other parameters in saliva
samples

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Benefits:

- Non invasive sampling, can be done also by patients
- Measurement of the free (unbound) fraction of steroids
- Much higher diagnostic relevance with reproducible and reliable results, compared to serum analytics
- Simple and patient-friendly measurement of hormone profiles (Diurnal-, monthly **profiles**)
- Easy sampling in babies and kids
- The measurement of steroids in **veterinary** saliva samples is possible without any further pre-treatment

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FDA cleared

Cortisol

- Sensitivity: 0.537 ng/mL
- Standard range: 2 - 80 ng/mL

17 α -OH Progesterone

- Sensitivity: 2.5 pg/mL
- Standard range: 10 - 1000 pg/mL

Estriol

- Sensitivity: 3 pg/mL
- Standard range: 5 - 4000 pg/mL

FDA cleared

Progesterone

- Sensitivity: 3.8 pg/mL
- Standard range: 10 - 5000 pg/mL

FDA cleared

Testosterone

- Sensitivity: 1.9 pg/mL
- Standard range: 10 - 5000 pg/mL

Cortisol HS (high sensitive)

- Sensitivity: 0.012 ng/mL
- Standard range: 0.1 - 30 ng/mL

FDA cleared

Estradiol

- Sensitivity: 0.39 pg/mL
- Standard range: 1 - 100 pg/mL

Estrone*

- High Sensitivity
- Short incubation time

DHEA

- Sensitivity: 2.19 pg/mL
- Standard range: 10 - 1440 pg/mL

* Coming soon

Sports Medicine/Sportmedizin

Lactate*

Lactate is generated by lactate dehydrogenase (LDH) under hypoxic or anaerobic conditions. Monitoring lactate level is, therefore, useful in the study of cellular and animal physiology, providing an important index of the balance between tissue oxygen demand and utilization.

DRG Lactate Assay Kit is based on lactate dehydrogenase catalyzed oxidation of lactate, in which the formed NADH is coupled to the formazan (WST) / phenazine methosulfate (PMS)

* Coming soon

Reagent. The intensity of the product color, measured at 450 nm, is proportional to the lactate concentration in the sample.

Cat. No.:	SLV-4764
Incubation time:	30/10 min.
Testprinciple:	Quantitative/ Competitive binding
Standard Range:	1 - 12 mmol/L
Specimen Sample Volumen:	20 μ L Saliva
Sensitivity:	0.08 mmol/L

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Sampling device

As saliva contains only tiny amounts of protein non-polar analytes such as steroids tend to adsorb to various plastic surfaces. Glass would tend to be the best material for saliva sampling and for the most part Polyethylene, although it is still in use in many laboratories worldwide, should be avoided. However there is an inexpensive polypropylene device available from DRG (Sali-Tubes). It does not show any interference in saliva even in Progesterone testing, which is the most critical steroid.

Blood in saliva

Traces of blood might interfere with the results of saliva testing. It is therefore important that the salivary samples should not show any red colour. In cases of visible blood traces in saliva the sample should be discarded. After rinsing the device several times with tap water a new saliva sample can be taken but only after waiting for approximately 15 minutes.

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
DRG Diagnostics

DRG Instruments GmbH, founded in 1973 by Dr. Geacintov, subsidiary of DRG Intl. Inc., USA, is a manufacturer of ELISAS. The DRG Group operates through a network of DRG subsidiaries in Germany, Poland, Russia and the Czech Republic and through distributors worldwide.

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ELISAS that perform


DRG develops and manufactures diagnostic ELISA test kits for use in clinical and research laboratories. The experience of our production and management team guarantees to provide high quality products, competitive prices and excellent customer service.

DRG works to DIN EN ISO 9001:2000, ISO 13485:2007 and ISO 13485:2003 under CMDCAS standard, certified by TÜV Rheinland Product Safety GmbH, an indication of our commitment to customer service, quality control and improved health care.

DRG Instruments GmbH, Germany
Frauenbergstraße 18
D-35039 Marburg
Tel. +49 (0) 64 21/170 00,
Fax +49 (0) 64 21/17 00 50
Internet: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by

DRG International Inc. USA
1167 U.S. Highway 22 East
Mountainside, N.J. 07092 USA
Phone: +1 (908) 233-2079
Fax +1 (908) 233-0758
Internet: www.drg-international.com
E-mail: corp@drg-international.com