# AUTOSTAT™II Anti-dsDNA

Anti-dsDNA

**REF** FGA01

## PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-dsDNA enzyme linked immunosorbent assay (ELISA) specific for antibodies to double stranded DNA

## INTENDED USE

Enzyme-linked immunosorbent assay method for the quantitative determination of specific IgG autoantibodies to dsDNA in human serum. The results of the anti-dsDNA assay can be used as an aid in the diagnosis of auto-immune diseases including Systemic Lupus Erythematosus (SLE). Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the HYCOR HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

## INTRODUCTION

Anti-double-stranded (ds) DNA antibodies are found along with a wide variety of antibodies recognising other nuclear components in Systemic Lupus Erythematosus (SLE). SLE is a disease of unknown aetiology, being characterised by widespread organ involvement including renal, CNS and skin as well as joint disease. Anti-dsDNA antibodies have been suggested to be directly involved in the disease process, being included in the ARA criteria for the diagnosis of the disease. This involvement is thought to be due to either dsDNA-anti-dsDNA complexes which have been detected in some patients or to the cross-reaction of anti-dsDNA antibodies with other molecules such as the phospholipids or glycosaminoglycans of connective tissues and cell surfaces. In SLE patients, levels of anti-dsDNA antibodies can be used in the monitoring of disease activity, since rapid increases have been reported to precede exacerbation of the disease and falls in antibody levels often accompany drug therapy. Anti-dsDNA antibodies are not detected in the sera of patients with drug induced lupus.

## Alternative Test Comparison

A wide range of immunochemical techniques have been used for the detection of anti-dsDNA antibodies. These include immunofluorescence, radioimmunoassay, haemagglutination and ELISA.

## Immunofluorescence

Immunofluorescence detection of anti-dsDNA antibodies has, in recent years, centred on the use of *Crithidia luciliae*, a protozoan with circular dsDNA. This assay has the disadvantage of being semi-quantitative.

## Radioimmunoassay

Radioimmunoassays which use [<sup>14</sup>C]- or [<sup>125</sup>I]-labelled DNA depend on the precipitation of the DNA when incubated in serum containing anti-DNA antibodies. These assays can be subject to false positive results if some other agent which is able to bind to DNA is present in the serum. It is also possible that anti-DNA antibodies of low affinity will not be detected by these techniques. Problems inherent with the handling of radiochemicals are also encountered.

## Agglutination Assays

Agglutination assays which depend on the agglutination of DNA coated particles, such as animal erythrocytes or inert chemical particles, suffer from the same disadvantages as above. In addition, the use of animal erythrocytes makes possible the detection of heterophilic antibodies leading to false positives.

## Enzyme Immunoassays

Solid phase immunoassays or ELISA's where the DNA is absorbed to a plastic surface, such as a microtitre plate, will have the advantage that antibody binding is measured directly. In addition, these assays can be modified to allow selective measurement of individual classes or subclasses of immunoglobulin. The major drawback to these assays has been the difficulty in attaching DNA to the solid phase, a difficulty which has been successfully overcome in the Autostat™II assay system.


## PRINCIPLE OF THE PROCEDURE

The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgG monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow. The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

COLORLESS

↓  
BLUE

ADD STOP SOLUTION

↓  
YELLOW

- |                     |
|---------------------|
| 1. Antigen          |
| 2. Analyte          |
| 3. Enzyme Conjugate |
| 4. Substrate        |

## MATERIALS PROVIDED

- Mouse anti-human IgG conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are colored pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakapart wells. The wells are coated with purified dsDNA from calf thymus.

- Standards: Four vials containing 1ml of standard. The standards are calibrated to International Reference Wo/80, are in IU/ml and contain human antisera. The standards contain 0.09% sodium azide as a preservative. The concentrations allocated to the standards are

Standard 1	50 IU/ml
Standard 2	200 IU/ml
Standard 3	500 IU/ml
Standard 4	1000 IU/ml

- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.

- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H<sub>2</sub>SO<sub>4</sub> stop solution.

## WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue. Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

## MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

## STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

## INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

A 'hook effect' may only be seen with very high samples which are above the assay range. No hook effect is seen up to 2686 IU/ml. This is a sample specific effect.

## SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

## LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out SLE or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat™II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that anti-dsDNA antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

## REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/50 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat™II kits are supplied in a ready-to-use format.

## MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/50 in sample diluent by adding 10µl to 490µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

## HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HYTEC automated instrument.

Dilutions of controls and samples are performed automatically by the HYTEC instrument.

Assay results are calculated by the HYTEC system automatically.

## CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or quantitatively:

### Qualitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2}}{\text{OD of standard 2}} \times \text{OD of sample or control}$$

The concentration of standard 2 is 200 IU/ml.

This method provides a qualitative result only

### Quantitative Calculation:

Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

### AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.062	0.036	0.049	
Standard 1	0.187	0.202	0.194	
Standard 2	0.618	0.621	0.620	
Standard 3	1.291	1.151	1.221	
Standard 4	2.291	2.115	2.203	
Sample 1	0.884	0.784	0.834	275 IU/ml
Sample 2	0.090	0.083	0.086	13 IU/ml

### INTERPRETATION OF RESULTS

A total of 161 samples were assayed. 34 of the samples were from normal subjects and 127 from people with an autoimmune condition. The cut-offs shown below were based on 81 non-dsDNA positive samples tested.

IU/ml dsDNA	Negative	Equivocal	Positive
	<40	40 - 60	>60

It is recommended that equivocal samples be retested using a subsequent sample.

### AUTOMATED PERFORMANCE CHARACTERISTICS

#### Analytical Specificity and Sensitivity

Of a panel of 'normal' asymptomatic individuals all samples gave a negative result. No cross-reactivity with anti-thyroglobulin, anti-TPO, anti-GPC, various ENAs, ANCA, anti-GBM or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 13 IU/ml.

#### Relative Specificity and Sensitivity

The automated assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Automated</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	28	4	4
	Equivocal	1	1	2
	Negative	4	9	108

Relative sensitivity = 87.5 %  
 Relative specificity = 96.4 %  
 Overall agreement = 94.4 %

Equivocal results were omitted from the above calculations.

#### Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	335	365	457
%cv	5.2	6.9	3.9

Inter-assay variation:	A	B	C
x	431	585	1156
%cv	9.4	8.4	9.8

### QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

### MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.058	0.059	0.058	
Standard 1	0.238	0.258	0.248	
Standard 2	0.806	0.874	0.840	
Standard 3	1.748	1.709	1.729	
Standard 4	2.536	2.524	2.530	
Sample 1	0.688	0.696	0.692	162.2
Sample 2	0.106	0.115	0.111	13.6

### INTERPRETATION OF RESULTS

A total of 200 samples were assayed. 50 of the samples were from normal subjects and 150 from people with an autoimmune condition. The cut-offs shown below were based on 107 non-dsDNA positive samples tested.

IU/ml dsDNA	Negative	Equivocal	Positive
	<40	40 - 60	>60

It is recommended that equivocal samples be retested using a subsequent sample.

### MANUAL PERFORMANCE CHARACTERISTICS

#### Analytical Specificity and Sensitivity

Of a panel of 'normal' asymptomatic individuals all samples gave a negative result. No cross-reactivity with anti-thyroglobulin, anti-TPO, anti-GPC, various ENAs, ANCA, anti-GBM or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 4.0 IU/ml.

#### Relative Specificity and Sensitivity

The manual assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Manual</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	63	3	1
	Equivocal	5	5	11
	Negative	0	2	110

Relative sensitivity = 98.4 %  
 Relative specificity = 100 %  
 Overall agreement = 99.4 %

Equivocal results were omitted from the above calculations.

#### Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	93.9	215.3	337.2
%cv	5.5	5.4	3.6

Inter-assay variation:	A	B	C
x	93.9	215.3	337.2
%cv	12.3	14.9	8.1

**Hycor Biomedical**  
 7272 Chapman Ave  
 Garden Grove  
 California 92841  
 U.S.A.  
 +1 800 382 2527 (Tel)  
 +1 714 901 1264 (Fax)  
[www.hycorbiomedical.com](http://www.hycorbiomedical.com)

**Advena Ltd.**  
 Pure Offices  
 Plato Close  
 Warwick CV34 6WE  
 U.K.  
 +44(0) 1926 800153 (Tel)

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### BIBLIOGRAPHY

- Dubois EL. (1974) Lupus Erythematosus. University of Southern California Press.
- Hughes GRV. et al. (1971) Anti-DNA activity in systemic lupus erythematosus: a diagnostic and therapeutic guide. Ann. Rheum. Dis., **30**, 259-264.
- Whinfield JB. et al (1975) Role of DNA-anti-DNA complexes in systemic lupus erythematosus. Scan. J. Rheumatol., **11**, 59-62.
- Eilat D. (1985) Cross-reactions of anti-DNA antibodies and the central dogma of lupus nephritis. Immunol. Today., **6**, Issue 4.
- Swaak AJK (1982) Prognostic value of anti-dsDNA in SLE. Ann. Rheum. Dis., **41**, 388-395.
- Brinkman K, Termaat R, Van den Brink H, Berden J (1991) The specificity of the anti-dsDNA ELISA. J. Immunological Methods., **139**, 91-100.
- Smeenk, Brinkman K, Termaat R, Van den Brink H (1991) Anti-dsDNA: choice of assay in relation to clinical value. Rheumatology International, Spring 1991.

**Hycor Biomedical**  
7272 Chapman Ave  
Garden Grove  
California 92841  
U.S.A.  
+1 800 382 2527 (Tel)  
+1 714 901 1264 (Fax)  
[www.hycorbiomedical.com](http://www.hycorbiomedical.com)

**Advena Ltd.**  
Pure Offices  
Plato Close  
Warwick CV34 6WE  
U.K.  
+44(0) 1926 800153 (Tel)

### Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual	Hytec
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	4 IU/ml	13 IU/ml
Relative specificity/Spécificité relative/Spesificitat/Specificite/ Especificidad	100%	96.4%
Relative sensitivity/Sensibilit� relative/ Empfindlichkeit / Sensibilit�/ Sensibilidad relativa	98.4%	87.5%
Intra-assay variation/ la reproductibilit� intra-lot/ Intra Testvariationen/Riproducibilit� intralotto/Variaci�n intraan�lisis	A 5.5%	A 5.2%
	B 5.4%	B 6.9%
	C 3.6%	C 3.9%
Inter-assay variation/ la reproductibilit� inter-lot/ Intra Testvariationen /Riproducibilit� interlotto/ Variaci�n interan�lisis	A 12.3%	A 9.4%
	B 14.9%	B 8.4%
	C 8.1%	C 9.8%

	English	Français	Deutsch	Italiano	Esp�ol
<b>LOT</b>	Batch / Lot code	Code du lot	Chargenbezeichnung	Codice del lotto	Codigo de lote
Expiry / Use by	Expiry / Use by	Date du expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
Manufactured by	Manufactured by	Fabriqu� par	Hergestellt von	Fabbricato da	Fabricado por
Store at	Store at	Conserver �	Lagerung bei	Conservare a	Conservar a
See instructions for use	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
<b>IVD</b>	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	S�lo para el uso diagn�stico in vitro
<b>CONTROL -</b>	Negative control	Contr�le n�gatif	Negativ Kontrolle	Controllo negativo	Control negativo
<b>CONTROL +</b>	Positive control	Contr�le positif	Positiv Kontrolle	Controllo positivo	Control positivo
<b>REF</b>	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	N�mero de cat�logo
<b>ECIREP</b>	Authorized representative	Repr�sentant Autoris�	Autorisierter Repr�sentant	Rappresentante Autorizzato	Representante Autorizado
$\Sigma$	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enth�lt ausreichend Reagenzien f�r <n> Tests	Contenuto sufficiente per <n> test	Contiene lo necesario para <n> test
<b>CAL</b>	Calibrator	Calibrateur	Kalibrator	Calibratore	Calibrador
<b>CONJ</b>	Conjugate	Conjug�	Konjugat	Conjugate	Conjugar
<b>MT PLATE</b>	Microtiter plate	Plaque de microtitration	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
<b>SAMPLEDIL</b>	Sample diluent	Diluant	Proben verd�nner	Diluire prova	Diluir muestra
<b>BUF WASH 20X</b>	Wash buffer Concentrate	Solution de lavage Concentr�	Waschpuffer Konzentrat	lavari di tampone concentrare	lavado de defensa concentrar
<b>SUBS</b>	Substrate	Substrat	Substrat	Substrato	Substrato
<b>SOLN STOP</b>	Stop solution	Solution d'arr�t	Stop Pl�sung	terminare la soluzione	Parar la soluc ion