AUTOSTAT™II Anti-Tg Anti-Thyroglobulin

REF FGA03

PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-Thyroglobulin enzyme linked immunosorbent assay (ELISA) specific for antibodies to thyroglobulin (Tg).

INTENDED USE

Enzyme-linked immunosorbent assay method for the quantitative determination of specific IgG autoantibodies to thyroglobulin in human serum. The results of the anti-Tg assay can be used as an aid in the diagnosis of auto-immune diseases including Hashimoto's Thyroiditis and Graves' Disease. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the Hycor HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

INTRODUCTION

Anti-thyroglobulin antibodies are characteristically present in the serum of patients with goitrous and atrophic thyroiditis or primary thyrotoxicosis. They are also found at lower levels in sera of patients with other thyroid diseases.

The highest concentrations occur in patients with the 'fibrous' and oxyphil variants of Hashimoto's goitre, whereas the juvenile forms of this disease are usually associated with lower levels. High levels are also common in primary myxoedema.

Almost 100% of patients with these diseases or with Hashitoxicosis have anti-Thyroglobulin, anti-Thyroid Peroxidase (TPO) antibodies or both. A small number of Graves' disease patients produce these antibodies but at lower levels. As these two autoantibodies vary independently of each other, and many patients have only anti-TPO antibodies (fewer have only anti-Thyroglobulin antibodies), it is recommended that both types are measured. An anti-Thyroid Peroxidase (TPO) test is also available from Hycor Biomedical.

Alternative Test Comparison

Shortly after the discovery of anti-thyroglobulin antibodies in the sera of patients with Hashimoto's disease, a second group of antibodies was found, using a complement fixation test which reacted with the microsomal fraction of thyroid tissue. Immunofluorescence and haemagglutination tests have subsequently been developed for anti-microsomal antibodies.

These assays are, however, limited by their semi-quantitative and subjective nature. In addition, haemagglutination titres may not be an accurate guide to the concentration of antibody in the serum because of differences in the agglutination properties of different immunoglobulin subclasses. Some sera contain a non-specific inhibitor of agglutination which can produce false negative results.

In recent years, a number of ELISA's have been developed for the determination of anti-thyroglobulin antibodies. This technique has the advantage of being sensitive, specific, quantitative and objective. The Autostat™II anti-Thyroglobulin is also suitable for use on automated immunoassay pipetting systems.

PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgG monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

COLORLESS

↓
BLUE

ADD STOP SOLUTION

↓
YELLOW

- | |
|---------------------|
| 1. Antigen |
| 2. Analyte |
| 3. Enzyme Conjugate |
| 4. Substrate |

MATERIALS PROVIDED

- Mouse anti-human IgG conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakapart wells. The wells are coated with purified human thyroglobulin.

- Standards: Four vials containing 1ml of standard. The standards are calibrated to IRP 65/93, are in IU/ml and contain human antisera. The standards contain 0.09 sodium azide as a preservative. The concentrations allocated to the standards are

Standard 1	150 IU/ml
Standard 2	500 IU/ml
Standard 3	2000 IU/ml
Standard 4	5000 IU/ml

- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.

- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H₂SO₄ stop solution.

WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing sodium azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue.

Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

A 'hook effect' may only be seen with very high samples which are above the assay range. No hook effect is seen up to 14,937 IU/ml. This is a sample specific effect.

SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out thyroid disease or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat™II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that anti-Thyroglobulin antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat™II kits are supplied in a ready-to-use format.

MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HYTEC automated instrument.

Dilutions of controls and samples are performed automatically by the HYTEC instrument.

Assay results are calculated by the HYTEC system automatically.

CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or quantitatively:

Qualitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2}}{\text{OD of standard 2}} \times \text{OD of sample or control}$$

The concentration of standard 2 is 500 IU/ml.

This method provides a qualitative result only

Quantitative Calculation:

Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.062	0.060	0.061	
Standard 1	0.274	0.304	0.289	
Standard 2	0.878	0.821	0.850	
Standard 3	1.977	1.967	1.972	
Standard 4	2.570	2.512	2.541	
Sample 1	0.869	0.845	0.857	504 IU/ml
Sample 2	1.363	1.447	1.405	826 IU/ml

INTERPRETATION OF RESULTS

A total of 120 samples were assayed. 25 of the samples were from normal subjects and 95 from people with an autoimmune condition. The cut-offs shown below were partly based on the statistical evaluation of 76 of these samples, none of which were anti-Tg positive.

IU/ml	Negative	Equivocal	Positive
Tg	<225	225 - 325	>325

It is recommended that equivocal samples be retested using a subsequent sample.

AUTOMATED PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of a panel of 25 'normal' asymptomatic individuals all samples gave a negative result. No cross-reactivity with Anti-Sm, anti-RNP, anti-Jo-1, anti-dsDNA, Rheumatoid Factor, anti-Gastric Parietal Cell, anti-Cardiolipin or ANCA was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 13.8 IU/ml.

Relative Specificity and Sensitivity

The automated assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Automated</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	28	1	1
	Equivocal	1	0	1
	Negative	1	1	86
Relative sensitivity		= 96.6 %		
Relative specificity		= 98.9 %		
Overall agreement		= 98.3 %		

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	1769	270	2216
%cv	9.3	3.4	6.3
Inter-assay variation:	A	B	C
x	1405	9407	1641
%cv	8.2	5.1	10.0

QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.077	0.074	0.076	
Standard 1	0.232	0.231	0.232	
Standard 2	0.488	0.581	0.535	
Standard 3	1.388	1.370	1.379	
Standard 4	2.284	2.233	2.259	
Sample 1	1.051	1.067	1.059	1327.8
Sample 2	0.079	0.081	0.080	4.0

INTERPRETATION OF RESULTS

A total of 120 samples were assayed. 25 of the samples were from normal subjects and 95 from people with an autoimmune condition. The cut-offs shown below were partly based on the statistical evaluation of 76 of these samples, none of which were anti-Tg positive.

IU/ml	Negative	Equivocal	Positive
Tg	<225	225 - 325	>325

It is recommended that equivocal samples be retested using a subsequent sample.

MANUAL PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of a panel of 25 'normal' asymptomatic individuals all samples gave a negative result. No cross-reactivity with Anti-Sm, anti-RNP, anti-Jo-1, anti-dsDNA, Rheumatoid Factor, anti-Gastric Parietal Cell, anti-Cardiolipin or ANCA was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 26.0 IU/ml.

Relative Specificity and Sensitivity

The manual assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Manual</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	69	1	0
	Equivocal	3	9	3
	Negative	0	1	94
Relative sensitivity		= 100 %		
Relative specificity		= 100 %		
Overall agreement		= 100 %		

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	370.6	569.7	2466.6
%cv	3.8	2.8	4.6
Inter-assay variation:	A	B	C
x	370.6	569.7	2466.6
%cv	9.9	12.3	9.8

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)

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BIBLIOGRAPHY

- Doniach D. et al. (1979) Goitrous autoimmune thyroiditis (Hashimoto's disease). Clin. Endocrinol. Metab., 8, 63-80.
- Goodburn R et al. (1982) The preparation on thyroid microsomal antigen for use in the indirect micro ELISA for the detection of an anti-thyroid microsomal antibody. Clin. Chem. Acta., 119, 291-297.
- Roman SH et al. (1984) Enzyme linked immunosorbent microassay and haemagglutination compared for detection of thyroglobulin and thyroid microsomal autoantibodies. Clin. Chem., 30, 246-251.
- Bigazzi PE et al. (1986) Antibodies to tissue specific endocrine, gastrointestinal and neurological antigens. In: Rose NR, Friedman H, Fahey JL. eds. Manual of Clinical Laboratory Immunology, 3rd Edition pp. 762-770.
- Tunbridge WM et al. (1977) The spectrum of thyroid disease in a community: The Wickham Survey. Clin. Endocrinology, 7, 481-493.

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual	Hytec
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	26.0 IU/ml	13.8 IU/ml
Relative specificity/Spécificité relative/Spécificite/ Especificidad	100%	98.9%
Relative sensitivity/Sensibilité relative/ Empfindlichkeit / Sensibilità/ Sensibilidad relativa	100%	96.6%
Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intralotto/Variación intraanálisis	A 3.8%	A 9.3%
	B 2.8%	B 3.4%
	C 4.6%	C 6.3%
Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis	A 9.9%	A 8.2%
	B 12.3%	B 3.4%
	C 9.8%	C 10.0%

	English	Français	Deutsch	Italiano	Español
LOT	Batch / Lot code	Code du lot	Chargenbezeichnung	Codice del lotto	Código de lote
Expiry / Use by	Expiry / Use by	Date du expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
Manufactured by	Manufactured by	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por
Store at	Store at	Conserver à	Lagerung bei	Conservare a	Conservar a
See instructions for use	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
IVD	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro
CONTROL -	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo
CONTROL +	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo
REF	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	Número de catálogo
ECIREP	Authorized representative	Représentant Autorisé	Autorisierter Repräsentant	Rappresentante Autorizzato	Representante Autorizado
Σ	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enthält ausreichend Reagenzien für <n> Tests	Contenuto sufficiente per <n> test	Contiene lo necesario para <n> test
CAL	Calibrator	Calibrateur	Kalibrator	Calibratore	Calibrador
CONJ	Conjugate	Conjugué	Konjugat	Conjugate	Conjugar
MT PLATE	Microtiter plate	Plaque de microtitration	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
SAMPLEDIL	Sample diluent	Diluant	Proben verdüner	Diluente prova	Diluir muestra
BUF WASH 20X	Wash buffer concentrate	Solution de lavage Concentré	Waschpuffer Konzentrat	lavari di tampone concentrare	lavado de defensa concentrar
SUBS	Substrate	Substrat	Substrat	Substrato	Substrato
SOLN STOP	Stop solution	Solution d'arrêt	Stop Lösung	terminare la soluzione	Parar la solución