AUTOSTAT™II Anti-GPC

Anti-Gastric Parietal Cell

REF FGA10

PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-Gastric Parietal Cell enzyme linked immunosorbent assay (ELISA) specific for antibodies to Gastric Parietal Cell (GPC).

INTENDED USE

Enzyme-linked immunosorbent assay method for the semi-quantitative determination of specific IgG autoantibodies to GPC in human serum. The results of the anti-GPC assay can be used as an aid in the diagnosis of auto-immune diseases including Pernicious Aamaemia. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the Hycor HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

INTRODUCTION

Circulating antibodies to intracellular antigens of gastric parietal cells, antibodies to intrinsic factor (type I) and antibodies to intrinsic factor B12 complex (type II) are all a common feature of pernicious anaemia (PA). Employing an indirect immunofluorescence assay, antibodies to gastric parietal cells are detected in more than 95% of patients with pernicious anaemia. Type I and II antibodies to intrinsic factor are technically difficult to detect, making anti-GPC antibodies the test of choice.

Anti-gastric parietal cell antibodies are also found in 60% of atrophic gastric and 22% of gastric ulceration cases without detectable anaemia because these conditions occur well before the onset of pernicious anaemia.

The occurrence of anti-GPC antibodies often arises in the presence of other autoimmune diseases, notably chronic thyroiditis, where 30% of patients with PA also possess at least one of the anti-thyroid antibodies. Other cases of serology overlap have been found in adrenocortical complications, pituitary and parathyroid abnormalities, and in patients with insulin dependent diabetes mellitus. Relatives of patients with PA were seen to be anti-GPC positive (47%) and anti-thyroid antibody positive (67%). The incidence of anti-GPC antibodies rises with age and is more common in women than men.

The Autostat™II Anti-Gastric Parietal Cell ELISA uses purified H⁺/K⁺ ATPase antigen coated onto microplate wells which, when used together with a standard curve, offers a sensitive, specific, semi-quantitative assay without associated background 'noise' due to tissue variability inherent in the classical immunofluorescence tests.

PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgG monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

COLORLESS

↓
BLUE

ADD STOP SOLUTION

↓
YELLOW

- | |
|---------------------|
| 1. Antigen |
| 2. Analyte |
| 3. Enzyme Conjugate |
| 4. Substrate |

MATERIALS PROVIDED

- Mouse anti-human IgG conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakpart wells. The wells are coated with GPC antigen purified from porcine stomach.

- Standards: Four vials containing 1ml of standard. The standards are calibrated arbitrary units and contain human antisera. The standards contain 0.09% sodium azide as a preservative. The concentrations allocated to the standards are

Standard 1	100 u/ml
Standard 2	200 u/ml
Standard 3	500 u/ml
Standard 4	1000 u/ml

- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.

- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H₂SO₄ stop solution.

WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue. Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

A 'hook effect' may only be seen with very high samples which are above the assay range. No hook effect is seen up to 9975 u/ml. This is a sample specific effect.

SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out pernicious anaemia or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the AutostatTMII kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that anti-GPC antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in AutostatTMII kits are supplied in a ready-to-use format.

MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HYTEC automated instrument. Dilutions of controls and samples are performed automatically by the HYTEC instrument.

Assay results are calculated by the HYTEC system automatically.

CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or semi-quantitatively:

Qualitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2}}{\text{OD of standard 2}} \times \text{OD of sample or control}$$

The concentration of standard 2 is 200 u/ml.

This method provides a qualitative result only

Semi-Quantitative Calculation:

Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.050	0.058	0.054	
Standard 1	0.570	0.567	0.568	
Standard 2	0.956	0.947	0.952	
Standard 3	1.641	1.562	1.602	
Standard 4	2.125	2.097	2.111	
Sample 1	0.060	0.095	0.078	5.3 u/ml
Sample 2	1.000	1.096	1.048	221.4 u/ml

INTERPRETATION OF RESULTS

A total of 139 samples were assayed. 40 of the samples were from normal subjects and 99 from people with an autoimmune condition. The cut-offs shown below were partly based on the statistical evaluation of 105 of these samples, all of which were anti-GPC negative.

u/ml	Negative	Equivocal	Positive
GPC	<100	100 - 125	>125

It is recommended that equivocal samples be retested using a subsequent sample.

AUTOMATED PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of a panel of 41 'normal' asymptomatic individuals 2 samples gave a positive result. No cross-reactivity with anti-dsDNA, anti-TPO, various ENAs, anti-pANCA, anti-cANCA or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 6.7 u/ml.

Relative Specificity and Sensitivity

The automated assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Automated</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	20	3	2
	Equivocal	2	0	9
	Negative	1	1	101
Relative sensitivity	=	90.9 %		
Relative specificity	=	99.0 %		
Overall agreement	=	97.6 %		

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	218.2	422.2	684.8
%cv	2.3	13.1	4.1
Inter-assay variation:	A	B	C
x	219.2	414.7	661.6
%cv	2.6	5.3	3.7

QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

Manual and HYTEC Correlation

A linear regression analysis of the correlation between the manual and the HYTEC GPC assays was performed. 141 samples were run on both kits. The r^2 value is 0.916.

MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.047	0.041	0.044	
Standard 1	0.479	0.475	0.477	
Standard 2	1.009	0.976	0.992	
Standard 3	1.448	1.323	1.385	
Standard 4	1.998	1.901	1.950	
Sample 1	0.877	0.842	0.859	172 u/ml
Sample 2	0.404	0.399	0.401	75 u/ml

INTERPRETATION OF RESULTS

A total of 139 samples were assayed. 40 of the samples were from normal subjects and 99 from people with an autoimmune condition. The cut-offs shown below were partly based on the statistical evaluation of 105 of these samples, all of which were anti-GPC negative.

u/ml	Negative	Equivocal	Positive
GPC	<100	100 - 125	>125

It is recommended that equivocal samples be retested using a subsequent sample.

MANUAL PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of a panel of 41 'normal' asymptomatic individuals 2 samples gave a positive result. No cross-reactivity with anti-dsDNA, anti-TPO, various ENAs, anti-pANCA, anti-cANCA or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 10.33 u/ml.

Relative Specificity and Sensitivity

The manual assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

		<u>Manual</u>		
		Positive	Equivocal	Negative
<u>Predicate</u>	Positive	23	0	2
	Equivocal	2	1	8
	Negative	1	1	101
Relative sensitivity	=	92.0 %		
Relative specificity	=	99.0 %		
Overall agreement	=	97.6 %		

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	131.8	341.5	638.2
%cv	3.0	4.5	6.1
Inter-assay variation:	A	B	C
x	131.8	341.5	638.2
%cv	10.0	9.4	14.6

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

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BIBLIOGRAPHY

- Bigazzi PE, Burek CL, Rone NR. (1986) Antibodies to tissue specific endocrine gastrointestinal and neurological antigens. In Rone NR., Freldman H., Fahey JL., Eds. Manual of Clinical Laboratory Immunology, 3rd edition, 762-770.
- Karlsson FA, Burman P, Loof L, March S. (1988) Major parietal cell antigen in autoimmune gastritis with pernicious anaemia is the acid producing H+/K+ATPase of stomach. J. Clin. Invest., **81**, 475-479.
- Whittingham S, Mackay IR (1985) Pernicious anaemia and gastric atrophy. In Rose NR., Mackay IR Eds, The Autoimmune Diseases, Orlando Academic Press.
- Callaghan JM, et al. (1993) Alpha and beta subunits of the gastric H+/K+ATPase are concordantly targeted by parietal cell antibodies associated with autoimmune gastritis. Autoimmunity, 1993, **16**, 289-295.

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual	Hytec
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	10.33 u/ml	6.7 u/ml
Relative specificity/Spécificité relative/Spesifitat/Specificite/ Especificidad	99.0%	99.0%
Relative sensitivity/Sensibilité relative/ Empfindlichkeit / Sensibilità/ Sensibilidad relativa	92.0%	90.9%
Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intra-lot/Variación intraanálisis	A 3.0% B 4.5% C 6.1%	A 2.3% B 13.1% C 4.1%
Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis	A 10.0% B 9.4% C 14.6%	A 2.6% B 5.3% C 3.7%

	English	Français	Deutsch	Italiano	Español
LOT	Batch / Lot code	Code du lot	Chargenbezeich- nung	Codice del lotto	Código de lote
Expiry / Use by	Date du expiration	Date de expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
Manufactu- red by	Fabriqué par	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por
Store at	Conservé à	Conservé à	Lagerung bei	Conservare a	Conservar a
See instructions for use	Consulter les instructions d'utilisation	Consulter les instructions d'utilisation	Gebrauchsanweis- ung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
IVD	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro
CONTROL -	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo
CONTROL +	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo
REF	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	Número de catálogo
ECIREP	Authorized representa- tive	Représentant Autorisé	Autorisierter Représentant	Rappresentante Autorizzato	Representante Autorizado
Σ	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enthält ausreichend Reagenzien für <n> Tests	Contenuto sufficiente per <n> test	Contiene lo necesario para <n> test
CAL	Calibrator	Calibrateur	Kalibrator	Calibratore	Calibrador
CONJ	Conjugate	Conjugué	Konjugat	Conjugate	Conjugar
MT PLATE	Microtiter plate	Plaque de microtitratio- n	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
SAMPLEDIL	Sample diluent	Diluant	Proben verdüner	Diluire prova	Diluir muestra
BUF WASH 20X	Wash buffer Concentrat- e	Solution de lavage Concentré	Waschpuffer Konzentrat	lavari di tampone concentrare	lavado de defensa concentrar
SUBS	Substrate	Substrat	Substrat	Substrato	Substrato
SOLN STOP	Stop solution	Solution d'arrêt	Stop Lösung	terminare la soluzione	Parar la soluc- ion