AUTOSTAT™II β-2 GPI IgG

β-2 Glycoprotein-I IgG

REF FGA22

PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-β-2 Glycoprotein-I IgG enzyme linked immunosorbent assay (ELISA) specific for IgG antibodies to β-2 Glycoprotein-I (β-2 GPI).

INTENDED USE

Enzyme-linked immunosorbent assay method for the semi-quantitative determination of specific IgG autoantibodies to β-2 GPI in human serum. The results of the anti-β-2 GPI assay can be used as an aid in the diagnosis of auto-immune diseases associated with elevated levels of anti-β-2 GPI antibodies including anti-phospholipid syndrome. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the Hycor HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

INTRODUCTION

Anti-phospholipid antibodies (aPL) are a group of immunoglobulins associated with binding of negatively charged phospholipids including cardiolipin and phosphatidyl serine. High serum levels of anti-phospholipid antibodies are frequently found in patients with recurrent arterial and venous thromboembolism, recurrent foetal loss and thrombocytopenia, the main features of anti-phospholipid syndrome (APS). These symptoms are also often present in cases of Systemic Lupus Erythematosus (SLE), when elevated levels of aPL are detected.

It has been reported that β2-Glycoprotein I (β2-GPI), a 50 kD plasma protein, is required for binding of anti-cardiolipin and anti-phospholipid antibodies and was more recently found to be the target for a specific group of aPL. The conformational epitope for anti-β2GPI antibodies is thought to develop when β2-GPI interacts with a lipid membrane composed of negatively charged phospholipids or, more recently reported, when it is adsorbed to the surface of a polystyrene microtitre plate.

Antibodies targeting β-2-GPI specifically have been found to be the dominant autoantibody type in anti-phospholipid positive sera and is widely accepted as a clear diagnostic indicator of APS and of SLE associated with APS.

All 3 isotypes (IgG, IgM and IgA) have been shown to be associated with APS. IgA anti-B2GPI in particular has been correlated significantly with venous thrombosis in these diseases. Although the levels of each isotype of B2GPI antibodies can correlate significantly with the other Isotypes, it is advised that all 3 isotypes are measured to ensure maximum sensitivity.

To increase the efficiency of clinical laboratory testing for APS, anti-cardiolipin IgM and IgG, β-2-GPI IgM and IgA ELISAs are also available from Hycor Biomedical Ltd.

PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgG monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

COLORLESS

↓
BLUE

ADD STOP SOLUTION

↓
YELLOW

- | |
|---------------------|
| 1. Antigen |
| 2. Analyte |
| 3. Enzyme Conjugate |
| 4. Substrate |

MATERIALS PROVIDED

- Sheep anti-human IgG conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakpart wells. The wells are coated with purified β-2 GPI from human plasma.

- Standards: Four vials containing 1ml of standard. The standards are in u/ml and contain human antisera. The standards contain 0.09% sodium azide as a preservative. The concentrations allocated to the standards are

Standard 1	10 u/ml
Standard 2	25 u/ml
Standard 3	50 u/ml
Standard 4	100 u/ml

- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.

- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H₂SO₄ stop solution.

WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue. Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C.

The opened kit should be used within one month.

INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

Within the limitations of the assay dynamic range, the equipment used and a dilution factor of 1/100, no hook effect was observable or of significance in this assay.

SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out anti-phospholipid syndrome or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat™II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that anti-β-2 GPI antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat™II kits are supplied in a ready-to-use format.

MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HYTEC automated instrument. Dilutions of controls and samples are performed automatically by the HYTEC instrument.

Assay results are calculated by the HYTEC system automatically.

CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or quantitatively:

Qualitative Calculation:

Only the assay blank and standard 2 should be run. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2} \times \text{OD of sample or control}}{\text{OD of standard 2}}$$

The concentration of standard 2 is 25 u/ml.

This method provides a qualitative result only

Semi-Quantitative Calculation:

An assay blank and standards 1, 2, 3 and 4 should be run. Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.095	0.092	0.094	
Standard 1	0.296	0.300	0.298	
Standard 2	0.651	0.635	0.643	
Standard 3	1.212	1.202	1.207	
Standard 4	2.247	2.336	2.292	
Sample 1	0.616	0.606	0.611	23.5 u/ml
Sample 2	1.487	1.550	1.518	64.8 u/ml

INTERPRETATION OF RESULTS

A total of 160 samples were assayed. 70 of the samples were from normal subjects, 50 samples were from patients with confirmed Anti-Phospholipid Syndrome and 40 samples from patients with other autoimmune conditions. The cut-off shown below was based on testing of 70 normal samples.

u/ml	Negative	Equivocal	Positive
β-2 GPI IgG	< 12	12 - 16	> 16

MANUAL PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of a panel of 70 'normal' samples, 1 sample gave a positive result.

A panel of 50 samples from patients clinically diagnosed with Anti-Phospholipid Syndrome were tested. 25 out of 50 samples resulted positive (50%).

Ten samples from patients with syphilis were tested, one gave a positive result.

40 samples from patients with other autoimmune conditions were tested with the following results;

DISEASE STATE	No. POSITIVE
Systemic Lupus Erythematosus (n = 20)	0
Rheumatoid Arthritis (n = 5)	0
ANCA/GBM (n = 11)	1
Thyroid Disease (n = 4)	0

The sensitivity of the assay was established by calculation of the mean plus three standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 2.76 u/ml.

Relative Specificity and Sensitivity

The assay was compared with another commercially available test using the above panel of APS samples, other disease state samples and 'normal' samples. The results are shown below:

Predicate Device		Autostat II		
		Positive	Negative	Equivocal
Predicate Device	Positive	17	1	0
	Negative	10	128	4

Relative sensitivity	94.4%
Relative specificity	92.8%
Overall agreement	92.9%

Reproducibility

Three samples across the assay range were tested in duplicate in a total of twenty assays.

Inter-assay variation	A	B	C
Mean concentration	12.3	52.6	24.4
%CV	14.3	10.7	9.2

Intra-assay variation	A	B	C
Mean concentration	12.3	52.6	24.4
Mean %CV	6.0	1.8	2.1

AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.050	0.056	0.053	
Standard 1	0.171	0.165	0.168	
Standard 2	0.373	0.386	0.380	
Standard 3	0.786	0.797	0.792	
Standard 4	1.723	1.743	1.733	
Sample 1	0.321	0.317	0.319	20.3 u/ml
Sample 2	0.212	0.216	0.214	12.3 u/ml

Reproducibility

Three samples across the assay range were tested in duplicate in a total of twenty assays.

Inter-assay variation	A	B	C
Mean concentration	63.9	25.3	21.3
%CV	9.2	12.0	6.6

Intra-assay variation	A	B	C
Mean concentration	63.9	25.3	21.3
Mean %CV	5.0	2.9	2.9


QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit, which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

Manual and Automated Correlation.

A linear regression analysis of the correlation between the manual and automated β-2 GPI IgG assays was performed. 58 samples that were run on both systems were used. The correlation coefficient is 0.93. 95% confidence intervals for the slope and intercept are 0.8162 to 1.007 and 0.8865 to 4.667 respectively.

BIBLIOGRAPHY

- Galli M et al (1990) Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein cofactor. *Lancet*. 335: 1544-1547.
- Hughes GRV (1993) The Antiphospholipid Syndrome: Ten Years On. *Lancet*. 342: 341-344
- Matsuura E & Koike T (1996) Beta-2-Glycoprotein I Autoantibodies. In *Autoantibodies*. Ed Peter AB and Schoenfeld Y.
- Tsutsumi A et al (1996) Antibodies to beta-2-glycoprotein I and clinical manifestations in patients with systemic lupus erythematosus. *Arthritis Rheum*. 39(9): 1466-74.
- Guerin J et al (1997) Antibodies to B2-Glycoprotein I; a specific marker for the antiphospholipid syndrome. *Clin Exp Immunol* 109:304.
- Sanmarco M et al (1997) Prevalence and clinical significance of IgG isotype anti-β2-Glycoprotein I antibodies in antiphospholipid syndrome: a comparative study with anticardiolipin antibodies. *J Lab Clin Med*. 129 (5): 499-506.
- Gomez-Pacheco L et al (1999) Serum anti-beta-2-glycoprotein I and anticardiolipin antibodies during thrombosis in systemic lupus erythematosus patients. *Am J Med*. 106 (4): 417-23.
- Rand JH (2002) Molecular Pathogenesis of the Antiphospholipid Syndrome. *Circulation Research*. 90: 29-37.
- Isotypes of anti-Beta2-glycoprotein I antibodies: association with thrombosis in patients with systemic lupus erythematosus. *J. Rheumatol* 2001 Mar;28(3) 520-4

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

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BIBLIOGRAPHY

- Galli M et al (1990) Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein cofactor. *Lancet*. 335: 1544-1547.
- Hughes GRV (1993) The Antiphospholipid Syndrome: Ten Years On. *Lancet*. 342: 341-344
- Matsuura E & Koike T (1996) Beta-2-Glycoprotein I Autoantibodies. In *Autoantibodies*. Ed Peter AB and Schoenfeld Y.
- Tsutsumi A et al (1996) Antibodies to beta-2-glycoprotein I and clinical manifestations in patients with systemic lupus erythematosus. *Arthritis Rheum*. 39(9): 1466-74.
- Guerin J et al (1997) Antibodies to B2-Glycoprotein I; a specific marker for the antiphospholipid syndrome. *Clin Exp Immunol* 109:304.
- Sanmarco M et al (1997) Prevalence and clinical significance of IgG isotype anti-β-2-Glycoprotein I antibodies in antiphospholipid syndrome: a comparative study with anticardiolipin antibodies. *J Lab Clin Med*. 129 (5): 499-506.
- Gomez-Pacheco L et al (1999) Serum anti-beta-2-glycoprotein I and anticardiolipin antibodies during thrombosis in systemic lupus erythematosus patients. *Am J Med*. 106 (4): 417-23.
- Rand JH (2002) Molecular Pathogenesis of the Antiphospholipid Syndrome. *Circulation Research*. 90: 29-37.
- Isotypes of anti-Beta2-glycoprotein I antibodies: association with thrombosis in patients with systemic lupus erythematosus. *J. Rheumatol* 2001 Mar;28(3) 520-4

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual	Hytec
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	2.76 u/ml	1.937 u/ml
Relative specificity/Spécificité relative/Spesifitat/Specificite/ Especificidad	92.8%	94.9%
Relative sensitivity/Sensibilité relative/ Empfindlichkeit / Sensibilità/ Sensibilidad relativa	94.4%	100%
Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intra-lot/Variación intraanálisis	A 6.0%	A 5.0%
	B 1.8%	B 2.9%
	C 2.1%	C 2.9%
Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis	A 14.3%	A 9.2%
	B 10.7%	B 12.0%
	C 9.2%	C 6.6%

	English	Français	Deutsch	Italiano	Español
LOT	Batch / Lot code	Code du lot	Chargenbezeichn ung	Codice del lotto	Código de lote
	Expiry / Use by	Date du expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
	Manufactu red by	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por
	Store at	Conserver à	Lagerung bei	Conservare a	Conservar a
	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweis ung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
IVD	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro
CONTROL -	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo
CONTROL +	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo
REF	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	Número de catálogo
ECREP	Authorized representa tive	Représentant Autorisé	Autorisierter Représentant	Rappresentante Autorizzato	Representante Autorizado
	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enthält ausreichend Reagenzien für <n> Tests	Contenuto sufficiente per <n> test	Contiene lo necesario para <n> test
CAL	Calibrator	Calibreur	Kalibrator	Calibratore	Calibrador
CONJ	Conjugate	Conjugué	Konjugat	Conjugate	Conjugar
MT PI ATF	Microtiter plate	Plaque de microtitratio n	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
SAMPLEDIL	Sample diluent	Diluant	Proben verdünnner	Diluire prova	Diluir muestra
BUF WASH 20X	Wash buffer Concentrat e	Solution de lavage Concentré	Waschpuffer Konzentrat	lavari di tampone concentrare	lavado de defensa concentrar
SUBS	Substrate	Substrat	Substrat	Substrato	Substrato
SOLN STOP	Stop solution	Solution d'arrêt	Stop Lösung	terminare la soluzione	Parar la soluc ion