# AUTOSTAT™II RF IgA

## Rheumatoid Factor IgA

**REF** FGA06

### PROPRIETARY AND COMMON NAMES

Hycor Biomedial Autostat™II Rheumatoid Factor IgA enzyme linked immunosorbent assay (ELISA) specific for Rheumatoid Factor IgA (RF IgA).

### INTENDED USE

Enzyme-linked immunosorbent assay method for the quantitative determination of specific IgA Rheumatoid Factor in human serum. The results of the RF IgA assay can be used as an aid in the diagnosis of Rheumatoid Arthritis when supported by other laboratory and clinical investigations. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the HYCOR HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

### INTRODUCTION

Immunoglobulins which bind to the Fc region of IgG (rheumatoid factors) are characteristically present in the serum of rheumatoid arthritis (RA) patients, although they may also be found at lower frequency in certain other diseases. Initial tests for rheumatoid factor (RF) were based on the agglutination of IgG-coated red blood cells or latex particles, but because these methods detect predominantly IgM RF, the presence of rheumatoid factors of other classes (which are much poorer agglutinins) was not recognised for some time. IgG, IgA, IgD and IgE rheumatoid factors have now been identified.

The presence of IgA RF in the serum of rheumatoid arthritis (RA) patients was first shown by immunoelectrophoresis of serum proteins eluted from insoluble IgG. Although usually found in association with IgM RF, IgA RF can occur in seronegative patients. Its frequency in RA has been estimated at 56-88%, with several reports of high incidence in patients with Sjogren's syndrome. Systemic Lupus Erythematosus (SLE) and bacterial endocarditis patients can also possess IgA RF.

In a prospective study of early RA patients, raised IgA RF levels were found to be significantly associated with the later development of erosive disease. Consequently, those patients with high levels received more treatment with specific drugs. A similar relationship between IgA RF and erosions in hands was found in a retrospective study which has led to the suggestion that IgA RF may be related to the development of erosions. In contrast there is no such association with IgM RF levels. Thus a test capable of determining IgA RF levels can be useful in patient prognosis and management.

### PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgA monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

### MATERIALS PROVIDED

- Mouse anti-human IgA conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.
- Microplate: One microplate is supplied which contains 12 strips of 8 breakapart wells. The wells are coated with human IgG.
- Standards: Four vials containing 1ml of standard. The standards are calibrated to IRP 64/2, are in IU/ml and contain human antisera. The standards contain 0.09% sodium azide as a preservative. The concentrations allocated to the standards are
 

|            |            |
|------------|------------|
| Standard 1 | 12.5 IU/ml |
| Standard 2 | 50 IU/ml   |
| Standard 3 | 100 IU/ml  |
| Standard 4 | 200 IU/ml  |
- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.
- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.
- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.
- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.
- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.
- Stop solution: One bottle containing 20ml of 0.25M H<sub>2</sub>SO<sub>4</sub> stop solution.

### WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue. Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

## MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

## STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

## INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

A 'hook effect' may only be seen with very high samples which are above the assay range. No hook effect is seen up to 1155.2 IU/ml. This is a sample specific effect.

## SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

## LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out rheumatoid arthritis or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat<sup>TM</sup>II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that RF IgA occurs at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

The use of this device in juveniles has not been established.

## REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat<sup>TM</sup>II kits are supplied in a ready-to-use format.

## MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

## HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HYTEC automated instrument. Dilutions of controls and samples are performed automatically by the HYTEC instrument. Assay results are calculated by the HYTEC system automatically.

## CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or quantitatively:

### Qualitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2}}{\text{OD of standard 2}} \times \text{OD of sample or control}$$

The concentration of standard 2 is 50 IU/ml.

This method provides a qualitative result only

### Quantitative Calculation:

Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

### AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

|               | OD    | OD    | mean OD | Result     |
|---------------|-------|-------|---------|------------|
| Reagent Blank | 0.071 | 0.071 | 0.071   |            |
| Standard 1    | 0.197 | 0.273 | 0.235   |            |
| Standard 2    | 0.948 | 0.926 | 0.937   |            |
| Standard 3    | 1.393 | 1.419 | 1.406   |            |
| Standard 4    | 2.021 | 1.936 | 1.978   |            |
| Sample 1      | 0.879 | 0.922 | 0.900   | 47.9 IU/ml |
| Sample 2      | 1.493 | 1.489 | 1.491   | 82 IU/ml   |

### INTERPRETATION OF RESULTS

A total of 119 samples were assayed. 33 of the samples were from normal subjects and 86 from people with an autoimmune condition. The cut-offs shown below were based on 81 RF IgA negative samples tested.

| IU/ml  | Negative | Equivocal | Positive |
|--------|----------|-----------|----------|
| RF IgA | <20      | 20 - 30   | >30      |

It is recommended that equivocal samples be retested using a subsequent sample.

### AUTOMATED PERFORMANCE CHARACTERISTICS

#### Analytical Specificity and Sensitivity

No cross-reactivity with anti-TPO, andti-Tg, anti-GPC, various ENAs, anti-cardiolipin, anti-GBM or ANCA was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 1.4 IU/ml.

#### Relative Specificity and Sensitivity

The automated assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

| <u>Predicate</u>     | <u>Automated</u> |           |          |
|----------------------|------------------|-----------|----------|
|                      | Positive         | Equivocal | Negative |
| Positive             | 36               | 5         | 2        |
| Negative             | 1                | 0         | 75       |
| Relative sensitivity | = 94.7 %         |           |          |
| Relative specificity | = 98.7 %         |           |          |
| Overall agreement    | = 97.4 %         |           |          |

Equivocal results were omitted from the above calculations.

#### Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

| Intra-assay variation: | A    | B     | C     |
|------------------------|------|-------|-------|
| x                      | 49.3 | 132.9 | 99.7  |
| %cv                    | 1.5  | 6.8   | 2.6   |
| Inter-assay variation: | A    | B     | C     |
| x                      | 57.2 | 89.7  | 117.8 |
| %cv                    | 3.3  | 3.7   | 4.7   |

### QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

### Manual and Automated Correlation.

A linear regression analysis of the correlation between the manual and the automated RF IgA assays was performed. 144 samples that were run on both kits were used. The r<sup>2</sup> value is 0.960.

### Clinical Studies.

The Autostat<sup>TM</sup> II RF IgA kit was evaluated in one internal and two external clinical studies. The results are summarised below.

| Patient Group           | No of Samples | No RF IgA Positive | % Positive |
|-------------------------|---------------|--------------------|------------|
| Normals                 | 113           | 2                  | 1.8        |
| Rheumatoid Arthritis    | 107           | 79                 | 73.8       |
| Rheumatological Samples | 478           | 97                 | 20.3       |
| Other Autoimmune        | 36            | 2                  | 5.6        |
| SLE                     | 6             | 3                  | 50.0       |

### MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

|               | OD    | OD    | mean OD | Result     |
|---------------|-------|-------|---------|------------|
| Reagent Blank | 0.105 | 0.103 | 0.104   |            |
| Standard 1    | 0.504 | 0.504 | 0.504   |            |
| Standard 2    | 1.242 | 1.263 | 1.253   |            |
| Standard 3    | 1.648 | 1.627 | 1.638   |            |
| Standard 4    | 2.012 | 2.062 | 2.037   |            |
| Sample 1      | 1.494 | 1.528 | 1.511   | 61.2 IU/ml |
| Sample 2      | 0.345 | 0.342 | 0.344   | 10.4 IU/ml |

### INTERPRETATION OF RESULTS

A total of 119 samples were assayed. 33 of the samples were from normal subjects and 86 from people with an autoimmune condition. The cut-offs shown below were based on 81 RF IgA negative samples tested.

| IU/ml  | Negative | Equivocal | Positive |
|--------|----------|-----------|----------|
| RF IgA | <20      | 20 - 30   | >30      |

It is recommended that equivocal samples be retested using a subsequent sample.

### MANUAL PERFORMANCE CHARACTERISTICS

#### Analytical Specificity and Sensitivity

No cross-reactivity with anti-TPO, andti-Tg, anti-GPC, various ENAs, anti-cardiolipin, anti-GBM or ANCA was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 0.6 IU/ml.

#### Relative Specificity and Sensitivity

The manual assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

| <u>Predicate</u>     | <u>Manual</u> |           |          |
|----------------------|---------------|-----------|----------|
|                      | Positive      | Equivocal | Negative |
| Positive             | 33            | 8         | 2        |
| Negative             | 1             | 0         | 75       |
| Relative sensitivity | = 94.3 %      |           |          |
| Relative specificity | = 98.7 %      |           |          |
| Overall agreement    | = 97.3 %      |           |          |

Equivocal results were omitted from the above calculations.

#### Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

| Intra-assay variation: | A    | B    | C    |
|------------------------|------|------|------|
| x                      | 17.5 | 47.6 | 93.9 |
| %cv                    | 7.7  | 3.0  | 4.0  |
| Inter-assay variation: | A    | B    | C    |
| x                      | 17.5 | 47.6 | 93.9 |
| %cv                    | 6.6  | 7.1  | 7.4  |

**Hycor Biomedical**  
 7272 Chapman Ave  
 Garden Grove  
 California 92841  
 U.S.A.  
 +1 800 382 2527 (Tel)  
 +1 714 901 1264 (Fax)  
 www.hycorbiomedical.com

**Advena Ltd.**  
 Pure Offices  
 Plato Close  
 Warwick CV34 6WE  
 U.K.  
 +44(0) 1926 800153 (Tel)

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#### BIBLIOGRAPHY

- Heimer R, Levin FM. (1966) On the distribution of rheumatoid factor among the immunoglobulins. *Immunochemistry*, **3**, 1-10.
- Dunne JV et al. (1979) IgA rheumatoid factor in the serum and saliva of patients with rheumatoid arthritis and Sjogren's syndrome. *Ann. Rheum. Dis.*, **38**, 161-166.
- Bampton JLM et al. (1985) Measurement of rheumatoid factors by an enzyme-linked immunosorbent assay (ELISA) and comparison with other methods. *Ann. Rheum. Dis.*, **44**, 14-19.
- Gioud-Paquet M. et al. (1987) IgM rheumatoid factor (RF), IgA RF, IgE RF and IgG RF detected by ELISA in rheumatoid arthritis. *Ann. Rheum. Dis.*, **46**, 65-71.
- Teitsson I et al. (1984) Prospective study of early rheumatoid arthritis (I). Association of the rheumatoid factor isotypes with disease activity. *Ann. Rheum. Dis.*, **43**, 673-678.
- Prospective study of early rheumatoid arthritis (II). Association of the rheumatoid factor isotypes with disease activity. *Ann. Rheum. Dis.*, **46**, 679-685.
- Amason JA et al. (1987) Relation between bone erosions and rheumatoid factor isotypes. *Ann. Rheum. Dis.*, **46**, 380-384.

**Hycor Biomedical**  
7272 Chapman Ave  
Garden Grove  
California 92841  
U.S.A.  
+1 800 382 2527 (Tel)  
+1 714 901 1264 (Fax)  
[www.hycorbiomedical.com](http://www.hycorbiomedical.com)

**Advena Ltd.**  
Pure Offices  
Plato Close  
Warwick CV34 6WE  
U.K.  
+44(0) 1926 800153 (Tel)

### Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

|  | Manual    | Hytec     |
|--|-----------|-----------|
| Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica              | 0.6 IU/ml | 1.4 IU/ml |
| Relative specificity/Spécificité relative/Spesifitat/Specificité/Especificidad   | 98.7%     | 98.7%     |
| Relative sensitivity/Sensibilité relative/ Empfindlichkeit / Sensibilità/ Sensibilidad relativa                                  | 94.3%     | 94.7%     |
| Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intralotto/Variación intraanálisis    | A 7.7%    | A 1.5%    |
|  | B 3.0%    | B 6.8%    |
|  | C 4.0%    | C 2.6%    |
| Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis | A 6.6%    | A 3.3%    |
|  | B 7.1%    | B 3.7%    |
|  | C 7.4%    | C 4.7%    |

|                          | English   | Français   | Deutsch  | Italiano   | Español  |
|--------------------------|---|--|--|--|--|
| <b>LOT</b>               | Batch / Lot code                                  | Code du lot                                      | Chargenbezeichnung   | Codice del lotto                                   | Codigo de lote                                       |
| Expiry / Use by          | Date du expiration                                | Verwendbar bis                                   | Utilizzare entro   | Fecha de caducidad                                 |  |
| Manufactured by          | Fabrique par                                      | Hergestellt von                                  | Fabbricato da  | Fabricado por                                      |  |
| Store at                 | Conservar à                                       | Lagerung bei                                     | Conservare a   | Conservar a  |  |
| See instructions for use | Consulter les instructions d'utilisation          | Gebrauchsanweisung beachten                      | Consultare le istruzioni per l'uso                           | Consulte las instrucciones de uso                  |  |
| <b>IVD</b>               | For in vitro diagnostic use only                  | Pour utilisation in vitro uniquement             | Nur zur in vitro Verwendung                                  | Per uso diagnostico in vitro                       | Sólo para el uso diagnóstico in vitro                |
| <b>CONTROL -</b>         | Negative control                                  | Contrôle négatif                                 | Negativ Kontrolle  | Controllo negativo                                 | Control negativo                                     |
| <b>CONTROL +</b>         | Positive control                                  | Contrôle positif                                 | Positiv Kontrolle  | Controllo positivo                                 | Control positivo                                     |
| <b>REF</b>               | Catalogue number                                  | Code produit                                     | Bestellnummer  | Numero di catalogo                                 | Número de catálogo                                   |
| <b>ECIREP</b>            | Authorized representative                         | Représentant Autorisé                            | Autorisierter Repräsentant                                   | Rappresentante Autorizzato                         | Representante Autorizado                             |
| $\Sigma$                 | Contains sufficient for <math>\Sigma</math> tests | Contenu suffisant pour <math>\Sigma</math> tests | Enthält ausreichend Reagenzien für <math>\Sigma</math> Tests | Contenuto sufficiente per <math>\Sigma</math> test | Contiene lo necessario para <math>\Sigma</math> test |
| <b>CAL</b>               | Calibrator  | Calibrateur                                      | Kalibrator   | Calibratore  | Calibrador   |
| <b>CONJ</b>              | Conjugate   | Conjugué   | Konjugat   | Conjugate  | Conjugar   |
| <b>MT PLATE</b>          | Microtiter plate                                  | Plaque de microtitration                         | Mikrotiterplatte   | Piatto Microtiter                                  | Plato Microtiter                                     |
| <b>SAMPLEDIL</b>         | Sample diluent                                    | Diluant  | Proben verdüner  | Diluire prova                                      | Diluir muestra                                       |
| <b>BUF WASH 20X</b>      | Wash buffer Concentrate                           | Solution de lavage Concentré                     | Waschpuffer Konzentrat                                       | lavari di tampone concentrare                      | lavado de defensa concentrar                         |
| <b>SUBS</b>              | Substrate   | Substrat   | Substrat   | Substrato  | Substrato  |
| <b>SOLN STOP</b>         | Stop solution                                     | Solution d'arrêt                                 | Stop Pflösung  | terminare la soluzione                             | Parar la solución                                    |