AUTOSTAT™II Anti-Gliadin IgA

IVD

Anti-Gliadin IgA

REF FGA20

PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-Gliadin IgA enzyme linked immunosorbent assay (ELISA) specific for IgA antibodies to gliadin.

INTENDED USE

Enzyme-linked immunosorbent assay method for the semi-quantitative determination of specific IgA autoantibodies to gliadin in human serum. The results of the anti-Gliadin IgA assay can be used as an aid in the diagnosis of diseases associated with elevated levels of anti-Gliadin antibodies including Coeliac Disease. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

This device has been validated for use with the HYCOR HYTEC automated EIA instrument. In addition, a manual assay protocol is described below.

For *in vitro* diagnostic use only.

INTRODUCTION

Anti-Gliadin antibodies are found in patients with Coeliac disease. This is a relatively common disease and can affect between one in three hundred to one in two thousand adults depending on the country. Gliadin is a mixture of around 50 components and has a molecular weight of 16-40kD. Gliadin can be divided into four major fractions according to electrophoretic mobility. Studies have shown that these fractions as well as prolamins contained in barley, oats and rye have toxic effects on patients with Coeliac disease.

In Coeliac disease the destruction of the gut villous structure increases the permeability of the absorptive area to macromolecules. This leads to increased antibody production to food proteins such as Gliadin which is the ethanol soluble fraction of wheat gluten.

The diagnosis of Coeliac disease currently involves taking jejunal biopsies as outlined in the European Society for Paediatric Gastroenterology and Nutrition (ESPGAN). Generally three tests are performed on children and two on adults. The stresses to the patient of undergoing several biopsies and changing their diet can be severe. There is therefore a great requirement for a less invasive method of testing for Coeliac disease especially in children.

Methods for detecting anti-gliadin antibodies include indirect immunofluorescence (IIF), mixed reverse (solid phase) passive antiglobulin hemadsorption, radioimmunoassay and enzyme linked immunosorbent assay (ELISA).

PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgA monoclonal antibody (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

MATERIALS PROVIDED

- Sheep anti-human IgA conjugate: One vial containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakpart wells. The wells are coated with purified gliadin from wheat gluten.

- Standards: Four vials containing 1ml of standard. The standards are calibrated to arbitrary units and contain human antisera. The standards contain 0.09% sodium azide as a preservative. The concentrations allocated to the standards are

Standard 1	5 u/ml
Standard 2	20 u/ml
Standard 3	50 u/ml
Standard 4	100 u/ml

- Positive control: One vial containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is color coded blue.

- Wash buffer: One bottle containing 50ml of wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H₂SO₄ stop solution.

WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to HIV 1 and 2, as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide (NaN₃) and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing NaN₃, flush away using copious amounts of water. Dispose with caution as NaN₃ can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue. Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

MATERIALS REQUIRED BUT NOT SUPPLIED

The items listed below may not all be required if using the HYTEC assay system. However, if using the HYTEC assay system additional reagents may be required. Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

INTERFERENCES AND HOOK EFFECT

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

A 'hook effect' may only be seen with very high samples which are above the assay range. No hook effect is seen up to 1813.8 u/ml. This is a sample specific effect.

SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out Coeliac Disease or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat™II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

A common occurrence noted in Coeliac Disease is IgA deficiency. This potential for false negatives should be borne in mind when interpreting the IgA anti-gliadin result. In these cases the anti-gliadin IgG result should be considered.

It should be noted that anti-Gliadin antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat™II kits are supplied in a ready-to-use format.

MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

HYTEC INSTRUMENT ASSAY PROTOCOL

Please refer to the HYTEC automated EIA system procedure manual for specific instructions on using this assay with the HY•TEC automated instrument. Dilutions of controls and samples are performed automatically by the HYTEC instrument. Assay results are calculated by the HYTEC system automatically.

CALCULATION OF RESULTS

The results of the assay may be calculated either qualitatively or semi-quantitatively:

Qualitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard 2}}{\text{OD of standard 2}} \times \text{OD of sample or control}$$

The concentration of standard 2 is 20 u/ml.

This method provides a qualitative result only

Semi-Quantitative Calculation:

Plot the blank corrected optical densities (OD's) of the standards against the concentration values, using a linear y-axis (OD) and a logarithmic x-axis (concentration). Use the reagent blank as the zero standard. The concentration value of the patient samples can then be determined from this calibration curve.

Alternatively, use a 4-parameter logistic curve fitting for the standard curve and for calculating results, using a log scale for the x-axis and a linear scale for the y-axis.

AUTOMATED TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.274	0.282	0.278	
Standard 1	0.522	0.571	0.546	
Standard 2	1.105	1.086	1.095	
Standard 3	1.728	1.668	1.698	
Standard 4	2.148	2.092	2.120	
Sample 1	1.147	1.201	1.174	21.9 u/ml
Sample 2	0.359	0.369	0.364	2.1 u/ml

INTERPRETATION OF RESULTS

A total of 161 samples were assayed. 37 of the samples were from normal subjects, 75 from untreated confirmed or suspected coeliac patients (both adult and paediatric) and 49 from people with another autoimmune condition. The age range of the normal samples tested was 18 to 60 years. The cut-offs shown below were based on 37 adult normal samples tested.

u/ml	Negative	Equivocal	Positive
Gliadin IgA	<5	5 - 6	>6

This cut-off was verified in a separate study with 50 paediatric samples (age range 3 months to 16 years).

It is recommended that equivocal samples be retested using a subsequent sample.

AUTOMATED PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

No cross-reactivity with anti-dsDNA, anti-TPO, anti-GPC, various ENAs, anti-cardiolipin, anti-GBM or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 1.0 u/ml.

Relative Specificity and Sensitivity

The automated assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

	Positive	Positive	<u>Automated</u> Equivocal	Negative	<u>Predicate</u>
	Positive	20	0	1	
	Negative	0	1	139	
Relative sensitivity	= 95.2 %				
Relative specificity	= 100 %				
Overall agreement	= 99.4 %				

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	27.4	8.3	44.6
%cv	4.1	8.1	4.2
Inter-assay variation:	A	B	C
x	27.4	16.3	59.2
%cv	9.1	6.3	6.5

QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

Manual and Automated Correlation.

A linear regression analysis of the correlation between the manual and the automated Gliadin IgA assays was performed. 136 samples that were run on both kits were used. The r² value is 0.991.

MANUAL TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

	OD	OD	mean OD	Result
Reagent Blank	0.264	0.257	0.261	
Standard 1	0.433	0.451	0.442	
Standard 2	0.883	0.894	0.888	
Standard 3	1.439	1.459	1.449	
Standard 4	1.969	1.927	1.948	
Sample 1	0.928	0.921	0.924	21.1 u/ml
Sample 2	1.656	1.659	1.658	44.6 u/ml

INTERPRETATION OF RESULTS

A total of 162 samples were assayed. 37 of the samples were from normal subjects, 76 from untreated confirmed or suspected coeliac patients (both adult and paediatric) and 49 from people with another autoimmune condition. The age range of the normal samples tested was 18 to 60 years. The cut-offs shown below were based on 37 adult normal samples tested.

u/ml	Negative	Equivocal	Positive
Gliadin IgA	<5	5 - 6	>6

This cut-off was verified in a separate study with 50 paediatric samples (age range 3 months to 16 years).

It is recommended that equivocal samples be retested using a subsequent sample.

MANUAL PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

No cross-reactivity with anti-dsDNA, anti-TPO, anti-GPC, various ENAs, anti-cardiolipin, anti-GBM or Rheumatoid Factor was seen.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 20 replicates of the zero standard which gave a value of 0.3 u/ml.

Relative Specificity and Sensitivity

The manual assay was compared to another commercially available test and was found to be substantially equivalent. The results are shown below:

	Positive	<u>Manual</u> Equivocal	Negative	<u>Predicate</u>
	Positive	19	0	2
	Negative	0	0	141
Relative sensitivity	= 90.5 %			
Relative specificity	= 100 %			
Overall agreement	= 98.8 %			

Equivocal results were omitted from the above calculations.

Reproducibility

Intra- and inter- assay variation were checked using a number of samples.

Intra-assay variation:	A	B	C
x	4.9	16.8	74.2
%cv	3.7	2.0	4.1
Inter-assay variation:	A	B	C
x	4.9	16.8	74.2
%cv	6.8	3.1	5.0

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

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BIBLIOGRAPHY

- Catassi C. (1996) Gliadin Antibodies. Autoantibodies, Edited by JB Peter and Y Shoenfeld, Elsevier Science Publications.
- Merridew SR, Wilson DV, Williams EJ. (1995) Anti-gliadin antibody measurement by chemiluminescence ELISA in the diagnosis of coeliac disease. J. Clin. Pathology., **48**, 509-512.
- Berger R, Schmidt G. (1996) Evaluation of six anti-gliadin antibody assays. J.Immunol. Methods., **191**, 77-86.
- Vogelsang H, Genser D, Wyatt J, Lochs H, Ferenci P, Granditsch G, Penner G. (1995) Screening for Coeliac Disease: A prospective study on the value of non-invasive tests. Am. J. Gastroenterology, Vol. **90**, No. 3, 394-398.
- Lindh E, Ljunghall S, Larsson K, Lavo B. (1992) Screening for antibodies against gliadin in patients with osteoporosis. J. Internal Medicine, **231**, 403-406.
- Burgin-Wolff A, Gaze H, Hadziselimovic F, Huber H, Lentze MJ, Nussle D, Reymond-Berthet C (1991) Anti-gliadin and anti-endomysium antibody determination for coeliac disease. Arc. Dis. Childhood, **66**, 941-947.

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual	Hytec
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	0.3 u/ml	1.0 u/ml
Relative specificity/Spécificité relative/Spesifitat/Specificite/ Especificidad	100%	100%
Relative sensitivity/Sensibilité relative/ Empfindlichkeit / Sensibilità/ Sensibilidad relativa	90.5%	95.2%
Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intralotto/Variación intraanálisis	A 3.7%	A 4.1%
	B 2.0%	B 8.1%
	C 4.1%	C 4.2%
Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis	A 6.8%	A 9.1%
	B 3.1%	B 6.3%
	C 5.0%	C 6.5%

	English	Français	Deutsch	Italiano	Español
LOT	Batch / Lot code	Code du lot	Chargenbezeichnung	Codice del lotto	Código de lote
Expiry / Use by	Expiry / Use by	Date de expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
Manufactured by	Manufactured by	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por
Store at	Store at	Conservé à	Lagerung bei	Conservare a	Conservar a
See instructions for use	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
IVD	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro
CONTROL -	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo
CONTROL +	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo
REF	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	Número de catálogo
ECIREP	Authorized representative	Représentant Autorisé	Autorisierter Repräsentant	Rappresentante Autorizzato	Representante Autorizado
Σ	Contains sufficient for <math>\lt;math> <td>Contenu suffisant pour <math>\lt;math> <td>Enthält ausreichend Reagenzien für <math>\lt;math> <td>Contenute sufficiente per <math>\lt;math> <td>Contiene lo necesario para <math>\lt;math> </td></td></td></td>	Contenu suffisant pour <math>\lt;math> <td>Enthält ausreichend Reagenzien für <math>\lt;math> <td>Contenute sufficiente per <math>\lt;math> <td>Contiene lo necesario para <math>\lt;math> </td></td></td>	Enthält ausreichend Reagenzien für <math>\lt;math> <td>Contenute sufficiente per <math>\lt;math> <td>Contiene lo necesario para <math>\lt;math> </td></td>	Contenute sufficiente per <math>\lt;math> <td>Contiene lo necesario para <math>\lt;math> </td>	Contiene lo necesario para <math>\lt;math>
CAL	Calibrator	Calibrateur	Kalibrator	Calibratore	Calibrador
CONJ	Conjugate	Conjugué	Konjugat	Conjugate	Conjugar
MT PLATE	Microtiter plate	Plaque de microtitration	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
SAMPLEDIL	Sample diluent	Diluant	Proben verdünner	Diluire prova	Diluir muestra
BUF	Wash buffer	Solution de lavage	Waschpuffer	lavari di tampone	lavado de defensa
WASH	Concentrate	Concentré	Konzentrat	concentrare	concentrar
ZOX					
SUBS	Substrate	Substrat	Substrat	Substrato	Substrato
SOLN	Stop solution	Solution d'arrêt	Stop Lösung	terminare la soluzione	Parar la solución
STOP					