AUTOSTAT™II ACA Isotype Anti-Cardiolipin Isotype

REF FGA37

PROPRIETARY AND COMMON NAMES

Hycor Biomedical Autostat™II Anti-Cardiolipin (ACA) Isotype enzyme linked immunosorbent assay (ELISA) specific for IgG, IgM and IgA antibodies to cardiolipin.

INTENDED USE

Enzyme-linked immunosorbent assay method for the semi-quantitative determination of IgG, IgM and IgA autoantibodies to cardiolipin in human serum. The results of the anti-cardiolipin assay can be used as an aid in the diagnosis of auto-immune diseases including anti-phospholipid syndrome. Levels of these autoantibodies are one indicator in a multi-factorial diagnostic regime.

For *in vitro* diagnostic use only.

INTRODUCTION

Anti-cardiolipin and anti-phosphatidyl serine antibodies along with others, such as the Lupus Anticoagulant (LAC), belong to the family of anti-phospholipid antibodies. Anti-cardiolipin antibodies are circulating serum antibodies often associated with recurrent arterial and venous thromboembolism, recurrent foetal loss and thrombocytopenia.

These symptoms are often present in cases of Systemic Lupus Erythematosus (SLE) and in many other conditions, both of an autoimmune and of a non-autoimmune nature. Some studies show that over 50% of SLE patients have one or more classes of anti-cardiolipin antibodies. The presence of this antibody serves as a marker for the risk of a thromboembolic event. Those SLE patients exhibiting high levels of these autoantibodies are 4 times more likely to have such an event than those not expressing the autoantibodies.

Anti-cardiolipin autoantibodies can be of any combination of the IgM, IgA and IgG classes. IgG antibodies are the most prevalent class of autoantibody and the class with the greatest clinical correlation. Samples found to have IgG levels in the 'high' anti-cardiolipin band are the most likely to display overt clinical symptoms. However, IgA and IgM autoantibodies are often found either alone or in association with the IgG class. Measurement of all three classes is therefore recommended.

Anti-cardiolipin are the most commonly measured anti-phospholipid antibodies. However, it is known that some patients with infectious diseases and other autoimmune disorders show some antibody activity against cardiolipin.

Patients with anti-phospholipid syndrome may also have autoantibodies which react with other phospholipids. These include antibodies to phosphatidyl serine. While anti-cardiolipin are the most commonly measured, the detection of other anti-phospholipid antibodies has been shown to be of clinical benefit. An anti-phosphatidyl serine ELISA is also available from Hycor Biomedical. To increase the efficiency of the clinical laboratory testing Hycor also offers an anti-cardiolipin screen and an anti-phosphatidyl serine screen ELISA.

PRINCIPLE OF THE PROCEDURE


The Autostat™II assay for detection of autoantibodies is a solid phase immunosorbent assay (ELISA) in which the analyte is indicated by a colour reaction of an enzyme and substrate. The Autostat™II wells are coated with purified antigen (1).

On adding diluted serum to the wells the antibodies (2) present bind to the antigen. After incubating at room temperature and washing away unbound material, horseradish peroxidase conjugated anti-IgG, anti-IgM or anti-IgA (3) is added, which binds to the immobilised antibodies.

Following further incubation and washing, tetra-methyl benzidine substrate (TMB) (4) is added to each well. The presence of the antigen-antibody-conjugate complex turns the substrate to a dark blue colour. Addition of the stop solution turns the colour to yellow.

The colour intensity is proportional to the amount of autoantibodies present in the original serum sample.

FIGURE 1

MATERIALS PROVIDED

- Anti-IgG, anti-IgM and anti-IgA conjugate: One vial of each containing 15ml of ready-to-use HRP conjugate. Conjugate diluent contains 0.05% Proclin 300. Conjugates are color coded pink.

- Microplate: One microplate is supplied which contains 12 strips of 8 breakapart wells. The wells are coated with purified cardiolipin from bovine heart.

- Standards: Three vials containing 1 ml of standard. The IgM and IgG standards are calibrated to International preparation of the Phospholipid Standardisation Laboratory, University of Louisville, USA., are in MPL U/ml and GPL U/ml respectively. The IgA standard is calibrated to arbitrary units and is in u/ml. Each of the standards contains human antisera and 0.09% sodium azide as a preservative. The concentrations allocated to the standards are:

IgG Standard	100 GPL U/ml
IgM Standard	60 MPL U/ml
IgA Standard	60 u/ml

- Positive control: Three vials containing 0.45 ml of concentrated positive control which contains human antisera and 0.09% sodium azide as a preservative.

- Negative control: One vial containing 0.45 ml of concentrated negative control which contains normal human serum and 0.09% sodium azide as a preservative.

- TMB substrate solution: One vial containing 15ml of ready-to-use tetra-methylbenzidine (TMB) substrate.

- Sample diluent: One bottle containing 100ml of ready-to-use sample diluent buffer. The buffer includes 0.09% Sodium azide. Sample diluent buffer is colour coded blue.

- Wash buffer: One bottle containing 50ml of ACA wash buffer concentrate. Wash buffer concentrate contains 0.06% Proclin 300.

- Stop solution: One bottle containing 20ml of 0.25M H₂SO₄ stop solution.

WARNINGS OR PRECAUTIONS

Warning - Potentially Hazardous and Biohazardous Materials

Sera used in the preparation of the standards and controls have been tested for the presence of antibodies to Human Immunodeficiency Virus (HIV 1 and 2), as well as for Hepatitis B Surface Antigen (HBsAg) and HCV and found to be negative. All material is tested with FDA approved assays. Because no test method can offer complete assurance that HIV, HBsAg or other infectious agents are absent, it is recommended that human serum based products be handled with the same precautions used for patient specimens.

Dispose of reagent solutions containing sodium azide and thimerosal as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.

The performance data represented here were obtained using specific reagents listed in the package insert. Do not use reagents from other manufacturers in the kits. Do not dilute or adulterate the kit reagents, unless directed by the kit protocol. Do not use the substrate solution if it has begun to turn blue.

Do not use heat-inactivated serum.

Reagents contain preservatives which may be toxic if ingested. Do not pipette by mouth. Avoid contact of reagents or patient samples with skin or mucous membranes. If contact occurs, immediately flush with large quantities of water. Avoid splashing or creation of aerosols. Reusable glassware must be thoroughly washed and rinsed so that it is free of all detergents.

Microplate washing is important. Improperly washed wells will give erroneous results. Do not allow the wells to dry between incubations. Do not vary reagents and incubation temperatures above or below room temperature (18 - 25°C).

MATERIALS REQUIRED BUT NOT SUPPLIED

Support reagents are available from your supplier.

1. Distilled or deionized water.
2. Wash bottle, automated or semi-automated microwell plate washing system.
3. Rack for sample dilution.
4. Micropipettes including multichannels capable of accurately delivering 5-1000µl (less than 3% cv).
5. Reagent reservoirs for multichannel pipettes.
6. One-liter graduated cylinder.
7. Disposal basins and 0.5% sodium hypochlorite (50 ml bleach in 950ml water).
8. Single or dual wavelength microplate reader with 450nm filter. If dual wavelength is used, set the reference filter to 600 - 650 nm.
9. Paper towels, pipette tips and timer

STORAGE AND SHELF LIFE

Store kit components at 2-8°C and do not use after the expiry date on the box outer label. Before use all components should be allowed to warm up to ambient temperature (18-25°C). After use, the plate should be resealed, the bottle caps replaced and tightened and the kit stored at 2-8°C. The opened kit should be used within three months.

INTERFERENCES

Grossly haemolysed, lipaemic or microbiologically contaminated samples should not be used. Samples with abnormally elevated levels of haemoglobin, bilirubin and especially EDTA may interfere with assay performance and accuracy.

SPECIMEN COLLECTION AND PREPARATION

Use serum in this procedure. It is most important to preserve the chemical integrity of a blood specimen from the moment it is collected until it is assayed. Obtain patient samples by non-traumatic venipuncture, using a vacuum tube or sterile syringe. If a syringe is used, transfer the blood immediately to a vacuum tube (plain red-top or serum separator).

Allow samples to clot at room temperature (18-25°C) for at least 20-30 minutes, until the clot just begins to retract. Spin the sample in a centrifuge. Following the centrifugation, transfer the cell-free serum to a tightly stoppered storage tube.

Do not use sera samples showing signs of haemolysis. If it is necessary to store a sample prior to analysis, it is recommended that, for a period of up to 72 hours, store the sample in a sealed container at 2-8°C. Freeze samples at -20°C if longer storage is required. Avoid repeat freeze-thawing.

LIMITATIONS OF USE

A negative result should not be used as a sole criterion to rule out anti-phospholipid syndrome or other autoimmune disease, but must be taken in relation to other clinical observations and diagnostic tests. While the precision of the Autostat™II kit is sufficient to allow samples to be measured in single determinations, this is done at the clinical laboratory's discretion. It is advised that duplicate determinations should be used to enable identification of potential pipetting error or to allow for confirmation in the equivocal range.

It should be noted that anti-cardiolipin antibodies occur at low levels in other autoimmune and non-autoimmune conditions. Therefore all other clinical observations and diagnostic tests should be taken into account for clinical diagnosis.

REAGENT PREPARATION

- Wash buffer: Measure and dilute 50 ml of the 20 X wash buffer to one liter with distilled or deionized water. Mix well before use. Store this solution at 2-8°C if it is not to be used immediately. The diluted wash buffer is stable at 2-8°C for one week.

Important note: Do not use wash buffer from other Autostat™II kits as this will affect the assay and produce unreliable results.

- Positive and negative controls: These are provided in a concentrated form and should be diluted 1/100 with sample diluent buffer before use. Prepare fresh control dilutions before each assay run. Vortex all samples and controls before testing.

The conjugate, substrate and stop solution in Autostat™II kits are supplied in a ready-to-use format.

MANUAL ASSAY PROTOCOL

1. Bring all reagents to room temperature (18-25°C).
2. Select sufficient microwells for the test. Remove protective covering and select sufficient wells to accommodate the patient samples, standards and assay controls. Each sample is recommended to be tested in duplicate.
3. Dilute all serum samples and assay controls 1/100 in sample diluent by adding 5µl to 495µl sample diluent. Standards do not require dilution.
4. Pipette 100µl of the standards, diluted control or diluted patient sample into the wells. To achieve blanking on the plate reader add a 'no serum' control of 100µl of sample diluent to the first two wells. This will act as the zero point for the curve fit.
5. Incubate the wells at room temperature (18-25°C) for 30 minutes.
6. Wash the wells three times with diluted wash buffer. This can be done manually with a multichannel pipette or on an automatic plate washer. Empty the wells, invert and tap dry on paper towel.
7. Add 100µl of ready-to-use conjugate to each well.
8. Incubate the wells at room temperature (18-25°C) for 15 minutes.
9. Repeat washing as in section 6. above.
10. Add 100µl of ready-to-use TMB substrate to each well.
11. Incubate the wells at room temperature (18-25°C) for 15 minutes.
12. Add 50µl of stop solution to each well. Tap gently to ensure uniform color distribution and read within 15 minutes.
13. To read the plate, ensure the base is free from moisture and no air bubbles are in the wells. Read the absorbance of the well contents at 450nm on a suitable plate reader. On readers equipped with a dual wavelength facility set the reference filter to 600 - 650nm.
14. Subtract the blank (or mean of blanks) from the optical densities of the standard, controls and patient samples. If the assay was performed in duplicate, the mean of the wells should be taken.

CALCULATION OF RESULTS

Semi-Quantitative Calculation:

Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit standard. Calculate the mean, blank corrected, absorbance value (OD) for duplicates of the kit controls and patient samples. Using the following algorithm, calculate the concentration of each of the samples:

$$\frac{\text{Concentration of standard}}{\text{OD of standard}} \times \text{OD of sample or control}$$

The standard concentrations are 100 GPL U/ml, 60 MPL U/ml or 60 arb. u/ml for IgG, IgM and IgA classes respectively.

TYPICAL RESULTS

The data below are for example only and should not be used for calculation of results. The results were calculated using the formula described in the 'Calculation of Results' section.

		OD	OD	Mean OD	Result
IgG ACA	Reagent Blank	0.049	0.051	0.050	
	Standard	2.040	1.970	2.005	
	Sample 1	0.195	0.187	0.191	7.2
	Sample 2	1.581	1.529	1.555	77.0
IgM ACA	Reagent Blank	0.070	0.068	0.069	
	Standard	1.273	1.311	1.292	
	Sample 1	0.304	0.269	0.287	10.7
	Sample 2	0.714	0.686	0.700	31.0
IgA ACA	Reagent Blank	0.122	0.116	0.119	
	Standard	1.136	1.334	1.235	
	Sample 1	0.631	0.606	0.618	26.8
	Sample 2	0.139	0.151	0.145	1.4

INTERPRETATION OF RESULTS

The assay cut-offs were determined by combining data from a panel of a minimum of 29 asymptomatic normal samples and a panel of ACA and other autoimmune positives. These are as follows:

Parameter	U/ml	Negative	Low Pos	Mod Pos	High Pos
IgG	GPL	< 10	10 - 20	20 - 40	> 40
IgM	MPL	< 10	10 - 20	20 - 30	> 30
IgA	arb. u/ml	< 10	10 - 20	20 - 30	> 30

MANUAL PERFORMANCE CHARACTERISTICS

Analytical Specificity and Sensitivity

Of 41 'normal' samples tested in the IgG assay, 3 samples were positive; 2 designated as 'weak' and one as 'moderate'. Of 37 samples supplied as potentially ACA positive, 24 (64.9%) were positive. 37 samples from people with a suspected autoimmune disease were tested. 12 of these were determined as weak positives and 3 as moderate positives.

Of 38 'normal' samples tested in the IgM assay, all gave a negative result. Of 25 samples supplied as potentially ACA positive, 16 (64%) were positive. 37 samples from people with a suspected autoimmune disease were tested. 2 of these samples were determined as strong positives. These were from suspected SLE and RA patients, in whom elevated levels of ACA are not unexpected. 11 of these samples were determined as weak positives.

Of 40 'normal' samples tested in the IgA assay, all gave a negative result. Of 22 samples supplied as potentially ACA positive, 10 (45.4%) were positive. Of 37 samples from people with a suspected autoimmune disease, 33 (89.2%) were negative. Of the 4 that were positive, 2 were from SLE patients, in whom elevated levels of ACA are not unexpected.

The sensitivity of the assay was established by calculation of the mean plus two standard deviations of a minimum of 22 replicates of the zero standard which gave the following values:

IgG ACA: 0.38 GPL U/ml

IgM ACA: 0.49 MPL U/ml

IgA ACA: 3.6 arb. u/ml

Relative Specificity and Sensitivity

The assay was compared with another commercially available test. The results for each isotype are shown below. Equivocal results have been omitted from the calculations.

IgG	Autostat II		Competitor	Relative sensitivity	Relative specificity	Overall agreement
	Positive	Negative				
	15	0		100%	97.3%	98.1%
	1	36				
IgM	Autostat II		Competitor	Relative sensitivity	Relative specificity	Overall agreement
	Positive	Negative				
	31	2		93.9%	90.9%	92.4%
	3	30				
IgA	Autostat II		Competitor	Relative sensitivity	Relative specificity	Overall agreement
	Positive	Negative				
	6	2		75%	100%	95.6%
	0	37				

Reproducibility

Inter-assay variation		A	B	C
IgG ACA	Mean of duplicates	8.2	20.4	49.7
	%CV	8.8	9.1	7.7
IgM ACA	Mean of duplicates	22.2	29.9	59.6
	%CV	10.2	12.4	10.2
IgA ACA	Mean of duplicates	12.8	56.5	66.3
	%CV	15.2	11.1	9.8

Intra-assay variation		A	B	C
IgG ACA	Mean of duplicates	7.4	18.6	81.1
	%CV	8.6	7.4	5.9
IgM ACA	Mean of duplicates	22.2	29.9	59.6
	%CV	3.1	5.2	4.3
IgA ACA	Mean of duplicates	12.7	35.7	52.2
	%CV	6.7	10.1	6.8

QUALITY CONTROL

Good laboratory practice indicates that with each assay run, one or more quality control samples of known antibody level should be analyzed as though they were clinical samples. Positive and negative control samples are supplied with each kit, which may be assayed with each run. The results of these quality control samples should fall within the limits indicated on the Certificate of Analysis.

Should the results fall outwith this range repeat the assay using freshly prepared controls. Should the results continue to fall outside the specific range, and after equipment, adherence to the protocol and laboratory procedure have been verified, seek assistance from the supplier. Do not report patient results if the control results fall outwith the acceptable ranges.

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

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BIBLIOGRAPHY

- Harris EN, Gharavi AE, Patel SP, Hughes GRV (1987) Evaluation of anti-cardiolipin antibody test. Report of an International Workshop, held 4 April 1986. Clin. Exp. Immunol., **68**, 215-222.
- Gharavi AE, Harris EN, Asherson RA, Hughes GRV (1987) Anti-cardiolipin antibodies: isotype distribution and phospholipid specificity. Ann. Rheum. Dis., **46**, 1-6.
- Escalante A, Brey RL, Mitchell Jr BD, Dreiner U (1995) Accuracy of anti-cardiolipin antibodies in identifying a history of thrombosis among patients with Systemic Lupus Erythematosus. Am. J. of Medicine, **98**, 559-565.
- Paggi A, Caccavo D, Ferri GM, et al. (1994) Anti-cardiolipin antibodies in autoimmune thyroid diseases. Clin. Endocrinol., **40**, 329-333.
- Cronin ME, Biswas RM, Van der Straeton C, Fleisher TA, Klippel JH. (1988) IgG and IgM anti-cardiolipin antibodies in patients with lupus with anti-cardiolipin antibody associated syndromes. J. Rheumatol., **15**, 795-798.
- Rote NS, Dostal-Johnson D, Ware Branch MS, Ware Branch D. (1990) Anti-phospholipid antibodies and recurrent pregnancy loss: correlation between the activated partial thromboplastin time and antibodies against phosphatidyl serine and cardiolipin. Am. J. Obstet. Gynecol., **163**, 575-584.

Hycor Biomedical
7272 Chapman Ave
Garden Grove
California 92841
U.S.A.
+1 800 382 2527 (Tel)
+1 714 901 1264 (Fax)
www.hycorbiomedical.com

Advena Ltd.
Pure Offices
Plato Close
Warwick CV34 6WE
U.K.
+44(0) 1926 800153 (Tel)

Performance Data/Performances/Typische Ergebnisse/Caratteristiche/Datos de actuación.

	Manual
Analytical sensitivity/Limite de détection/Analytische Empfindlichkeit/Sensibilità analitica/Sensibilidad analítica	IgG 0.38 GPL/U/ml IgM 0.49 MPL U/ml IgA 3.6 u/ml
Intra-assay variation/ la reproductibilité intra-lot/Intra Testvariationen/Riproducibilità intralotto/Variación intraanálisis	IgG: 8.6% 7.4% 5.9% IgM: 3.1% 5.2% 4.3% IgA: 6.7% 10.1% 6.8%
Inter-assay variation/ la reproductibilité inter-lot/ Intra Testvariationen /Riproducibilità interlotto/ Variación interanálisis	IgG: 8.8% 9.1% 7.7% IgM: 10.2% 12.4% 10.2% IgA: 15.2% 11.1% 9.8%

	English	Français	Deutsch	Italiano	Español
LOT	Batch / Lot code	Code du lot	Chargenbezeichnung	Codice del lotto	Código de lote
Expiry / Use by	Expiry / Use by	Date du expiration	Verwendbar bis	Utilizzare entro	Fecha de caducidad
Manufactured by	Manufactured by	Fabriqué par	Hergestellt von	Fabbricato da	Fabricado por
Store at	Store at	Conserver à	Lagerung bei	Conservare a	Conservar a
See instructions for use	See instructions for use	Consulter les instructions d'utilisation	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso
IVD	For in vitro diagnostic use only	Pour utilisation in vitro uniquement	Nur zur in vitro Verwendung	Per uso diagnostico in vitro	Sólo para el uso diagnóstico in vitro
CONTROL -	Negative control	Contrôle négatif	Negativ Kontrolle	Controllo negativo	Control negativo
CONTROL +	Positive control	Contrôle positif	Positiv Kontrolle	Controllo positivo	Control positivo
REF	Catalogue number	Code produit	Bestellnummer	Numero di catalogo	Número de catálogo
ECIREP	Authorized representative	Représentant Autorisé	Autorisierter Repräsentant	Rappresentante Autorizzato	Representante Autorizado
Σ	Contains sufficient for <n> tests	Contenu suffisant pour <n> tests	Enthält ausreichend Reagenzien für <n> Tests	Contenuto sufficiente per <n> test	Contiene lo necesario para <n> test
CAL	Calibrator	Calibrateur	Kalibrator	Calibratore	Calibrador
CONJ	Conjugate	Conjugué	Konjugat	Conjugate	Conjugar
MT PLATE	Microtiter plate	Plaque de microtitration	Mikrotiterplatte	Piatto Microtiter	Plato Microtiter
SAMPLEDIL	Sample diluent	Diluant	Proben verdüner	Diluire prova	Diluir muestra
BUF WASH 20X	Wash buffer Concentrate	Solution de lavage Concentré	Waschpuffer Konzentrat	lavari di tampone concentrare	lavado de defensa concentrar
SUBS	Substrate	Substrat	Substrat	Substrato	Substrato
SOLN STOP	Stop solution	Solution d'arrêt	Stop Lösung	terminare la soluzione	Parar la solución