



# Direct Enzymatic HbA1c Reagents and Lysis Buffer Safety Data Sheet 1107

according to Regulation (EC) No. 453/2010

Date of issue: 03/06/2015

Revision date: 04/06/2015

Version: B

## SECTION 1 Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form	Mixture
Product name	Direct Enzymatic HbA1c Reagent R1A, Reagent R1B, Reagent R2 and Lysis Buffer
Product code	DZ168A DZ168A-R1A, DZ168A-R1B, DZ168A-R2, DZ168A-L
Product group	Blend

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

No additional information available

#### 1.2.2. Uses advised against

No additional information available

### 1.3. Details of the supplier of the safety data sheet

Supplier: Diazyme Laboratories  
12889 Gregg Court  
92064 Poway  
T 858-455-4768

Authorized Representative: MDSS GmbH  
Schiffgraben 41  
Germany  
(+49) 511-6262-8630

### 1.4. Emergency telephone number

Emergency number Contact your local health authority or poison control center in an emergency. Manufacturer contact number for the US is as follows: (858) 455-4768

## SECTION 2 Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

### 2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

### 2.3. Other hazards

Contains PBT/vPvB substances  $\geq$  0.1% assessed in accordance with REACH Annex XIII

## SECTION 3 Composition/information on ingredients

### 3.1. Substance

Not applicable

### 3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
CHES	(CAS No) 103-47-9 (EC no) 203-115-6	2 (Lysis Buffer only)	Eye Irrit. 2, H319 STOT SE 3, H335 Skin Irrit. 2, H315

# Direct Enzymatic HbA1c Reagents and Lysis Buffer Safety Data Sheet 1107

according to Regulation (EC) No. 453/2010

Triton X-100 substance listed as REACH Candidate (4-(1,1,3,3-tetramethylbutyl)phenol, ethoxylated [covering well-defined substances and UVCB substances, polymers and homologues]) non-classified PBT substance non-classified vPvB substance	(CAS No) 9002-93-1	0.5 (R1A and R2 only)	Acute Tox. 4 (Oral), H302 Eye Dam. 1, H318
--	--------------------	--------------------------	---

Full text of H-phrases: see section 16

## SECTION 4 First aid measures

### 4.1. Description of first aid measures

First-aid measures after inhalation	Remove person to fresh air and keep comfortable for breathing.
First-aid measures after skin contact	Wash skin with plenty of water.
First-aid measures after eye contact	Rinse eyes with water as a precaution.
First-aid measures after ingestion	Call a poison center or a doctor if you feel unwell.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5 Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media Water spray. Dry powder. Foam. Carbon dioxide.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire Toxic fumes may be released.

### 5.3. Advice for firefighters

Protection during firefighting Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

## SECTION 6 Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Emergency procedures Ventilate spillage area.

#### 6.1.2. For emergency responders

Protective equipment Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

### 6.2. Environmental precautions

Avoid release to the environment.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up Take up liquid spill into absorbent material.  
Other information Dispose of materials or solid residues at an authorized site.

### 6.4. Reference to other sections

For further information refer to section 13.

## SECTION 7 Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling Ensure good ventilation of the work station. Wear personal protective equipment.  
Hygiene measures Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions Store in a well-ventilated place. Keep cool.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8 Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

# Direct Enzymatic HbA1c Reagents and Lysis Buffer Safety Data Sheet 1107

according to Regulation (EC) No. 453/2010

## 8.2. Exposure controls

Appropriate engineering controls	Ensure good ventilation of the work station.
Hand protection	Protective gloves
Eye protection	Safety glasses
Skin and body protection	Wear suitable protective clothing
Respiratory protection	In case of insufficient ventilation, wear suitable respiratory equipment
Environmental exposure controls	Avoid release to the environment.

## SECTION 9 Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state	Liquid
Colour	No data available
Odour	No data available
Odour threshold	No data available
pH	No data available
Relative evaporation rate (butylacetate=1)	No data available
Melting point	Not applicable
Freezing point	No data available
Boiling point	No data available
Flash point	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Flammability (solid, gas)	Not applicable
Vapour pressure	No data available
Relative vapour density at 20 °C	No data available
Relative density	No data available
Solubility	No data available
Log Pow	No data available
Viscosity, kinematic	No data available
Viscosity, dynamic	No data available
Explosive properties	No data available
Oxidising properties	No data available
Explosive limits	No data available

### 9.2. Other information

No additional information available

## SECTION 10 Stability and reactivity

### 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

### 10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

## SECTION 11 Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	Not classified
----------------	----------------



# Direct Enzymatic HbA1c Reagents and Lysis Buffer Safety Data Sheet 1107

according to Regulation (EC) No. 453/2010

Proper Shipping Name (IATA)	Not applicable
Proper Shipping Name (ADN)	Not applicable
Proper Shipping Name (RID)	Not applicable

## 14.3. Transport hazard class(es)

### ADR

Transport hazard class(es) (ADR)	Not applicable
----------------------------------	----------------

### IMDG

Transport hazard class(es) (IMDG)	Not applicable
-----------------------------------	----------------

### IATA

Transport hazard class(es) (IATA)	Not applicable
-----------------------------------	----------------

### ADN

Transport hazard class(es) (ADN)	Not applicable
----------------------------------	----------------

### RID

Transport hazard class(es) (RID)	Not applicable
----------------------------------	----------------

## 14.4. Packing group

Packing group (ADR)	Not applicable
Packing group (IMDG)	Not applicable
Packing group (IATA)	Not applicable
Packing group (ADN)	Not applicable
Packing group (RID)	Not applicable

## 14.5. Environmental hazards

Dangerous for the environment	No
Marine pollutant	No
Other information	No supplementary information available

## 14.6. Special precautions for user

### - Overland transport

No data available

### - Transport by sea

No data available

### - Air transport

No data available

### - Inland waterway transport

Carriage prohibited (ADN)	No
Not subject to ADN	No

### - Rail transport

Carriage prohibited (RID)	No
---------------------------	----

## 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15 Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions

Contains a substance on the REACH candidate list in concentration  $\geq 0.1\%$  or with a lower specific limit: 4-(1,1,3,3-tetramethylbutyl)phenol, ethoxylated [covering well-defined substances and UVCB substances, polymers and homologues] (CAS 9002-93-1)

Contains no REACH Annex XIV substances

# Direct Enzymatic HbA1c Reagents and Lysis Buffer Safety Data Sheet 1107

according to Regulation (EC) No. 453/2010

## 15.1.2. National regulations

### Germany

VwVwS Annex reference Water hazard class (WGK) 2, hazard to waters (Classification according to VwVwS, Annex 4)  
12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

### Netherlands

SZW-lijst van kankerverwekkende stoffen None of the components are listed  
SZW-lijst van mutagene stoffen None of the components are listed  
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding None of the components are listed  
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid None of the components are listed  
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling None of the components are listed

## 15.2. Chemical safety assessment

No additional information available

## SECTION 16 Other information

Full text of H- and EUH-statements

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H302	Harmful if swallowed
H315	Causes skin irritation
H318	Causes serious eye damage
H319	Causes serious eye irritation
H335	May cause respiratory irritation

SDS EU Diazyme with CLP

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product*