

GROWTH HORMONE AND GROWTH FACTOR RELATED PRODUCTS

IMMUNOASSAYS

Cat. No.	Name	Sample	Calibration range Limit of detection
RMEE35R	ALS Human ELISA	Serum, Plasma (EDTA, heparin, citrate)	0 – 40 mU/ml
	Sandwich ELISA, HRP-labelled antibody	10 µl	0.23 mU/ml
RCD017R	Growth Hormone Human ELISA	Serum	1 – 50 ng/ml
	Sandwich ELISA, HRP-labelled antibody	25 µl	0.2 ng/ml
RMEE022	Growth Hormone Sensitive Human ELISA	Serum, Plasma (EDTA, heparin, citrate)	50 – 1000 pg/ml
	Sandwich ELISA, Biotin-labelled antibody	100 µl	1.66 pg/ml
RSHAKRGH-010R	Growth Hormone Rat ELISA	Serum, Plasma (EDTA, heparin, citrate)	31.3 – 2000 pg/ml
	Sandwich ELISA, Biotin-labelled antibody	5 µl	16 pg/ml
RA19003R	Growth Hormone (GH) Rat ELISA	Cell culture medium, Plasma	
	Sandwich ELISA		0.5 ng/ml
RMEE20	IGF-1 Human ELISA	Serum, Plasma (EDTA, heparin, citrate)	1 – 50 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	20 µl	0.09 ng/ml
RMEE25R	IGF-1 Mouse/Rat ELISA	Serum	0.5 – 18 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	50 µl	0.029 ng/ml
RMEE30	IGF-2 Human ELISA	Serum, Plasma (EDTA, heparin, citrate)	0.02 – 3600 ng/ml
	Sandwich ELISA, Biotin-labelled antibody		0.02 ng/ml
RCD018R	IGFBP-1 Human ELISA	Serum	1 – 250 µg/l
	Sandwich ELISA, HRP-labelled antibody	25 µl	0.5 µg/l
RMEE01	IGFBP-1 Human ELISA	Serum, Cell culture medium, Urine, Amniotic fluid, Milk, Saliva	0 – 8 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	20 µl	0.02 ng/ml

GROWTH HORMONE AND GROWTH FACTOR RELATED PRODUCTS

Cat. No.	Name	Sample	Calibration range Limit of detection
RMEE05R	IGFBP-2 ELISA	Amniotic fluid, Cerebrospinal fluid, Milk, Plasma (EDTA), Saliva, Serum, Urine	2 – 80 ng/ml
	Sandwich ELISA, HRP-labelled antibody	15 µl	0.2 ng/ml
RMEE08R	IGFBP-2 Mouse/Rat ELISA	Serum	0.031 – 2 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	10 µl	0.01 ng/ml
RMEE04R	IGFBP-3 functional Human ELISA	Serum, Plasma (EDTA, heparin)	0.4 – 30 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	10 µl	0.18 ng/ml
RMEE03A	IGFBP-3 total Human ELISA	Cerebrospinal fluid, Plasma (Heparin), Serum	0.4 – 30 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	10 µl	0.1 ng/ml
RMEE031R	IGFBP-3 Mouse/Rat ELISA	Serum, Plasma, Cell culture supernatant	0.078 – 5 ng/ml
	Sandwich ELISA, Biotin-labelled antibody	5 µl	0.018 ng/ml

› Contact Information



BioVendor – Laboratorní medicína a.s.
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185, Fax: +420 549 211 460
E-mail: info@biovendor.com

BioVendor GmbH
Otto-Hahn-Straße 16, 34123 Kassel, Germany
Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
Phone: +43 1 890 9025, Fax: +43 1 890 5163
E-mail: infoAustria@biovendor.com

BioVendor, LLC
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
Phone: +1-800-404-7807, Phone: +1-828-575-9250
Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com