

PULMONARY DISEASE

IMMUNOASSAYS

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RSHAKRAL020SR	Albumin Urinary (S-type) Rat Assay	Immunturbidimetry/nephelometry	Serum, Urine, Plasma	6.17–500 µg/ml
RAG017R	Angiotensin-Like Protein 7 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Cell culture supernatant	0.063–4 ng/ml
RAQP4/96/2R	Aquaporin-4 Autoantibody Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	1.5–80 units/ml
RAG016R	Caspase-1 Mouse ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Cell culture supernatant	0.156–10 ng/ml
RD191022200	Club Cell Protein (CC16) Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	1.57 – 50 ng/ml
RD191180200R	CTRP9 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	25–800 pg/ml
RSHAKRIN-011RUR	Insulin (RTU) Mouse ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Culture medium	100–12000 pg/ml
RAG018R	Irisin ELISA	Competitive ELISA, Immobilized antigen	Serum, Plasma, Cell culture supernatant	0.001–5 µg/ml
RSCYK243882R	KL-6 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum	201 – 4020 U/ml
RD191259200R	Pancreatic-Derived Factor Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	0.0625–4 ng/ml
RAG019R	Periostin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell culture supernatant	78–5000 pg/ml
RAG020R	Periostin Mouse ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell culture supernatant	3–2000 pg/ml
RD191021100	PMN Elastase Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma (EDTA, citrate), CSF, BALF, Seminal plasma, Exudate	15.6 – 1000 ng/ml
RD191493200CS	S100A8 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	0.05 – 20 ng/ml
RD194218200CS	S100A9 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	1.563 – 50 ng/ml
RD191221200R	S100A12 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.063–4 ng/ml
RD191221200R	S100A12 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Bronchoalveolar lavage, Amniotic fluid	0.063 – 4 ng/ml
RD191225200R	S100G Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	25–800 pg/ml
RD191373200R	Secretoglobulin 3A2 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Amniotic fluid, Bronchoalveolar lavage, Plasma (EDTA, citrate, heparin), Serum	31.3 – 2000 pg/ml
RD191139200R	Surfactant Protein A Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (heparin, citrate), BALF, Amniotic fluid	1 – 50 ng/ml
RD194059101	Surfactant Protein D Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, BALF, Amniotic fluid	1.56 – 100 ng/ml
RD191400200CS	T-Cell Immunoglobulin and Mucin Domain Containing 4 (TIM4) Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	0.1 – 3.2 ng/ml
RD191158100R	Trefoil Factor 1 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin)	0.125 – 4 ng/ml
RD191159100R	Trefoil Factor 2 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma	0.25 – 10 ng/ml
RD191160200R	Trefoil Factor 3 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Bronchoalveolar lavage, Plasma	0.075 – 2.4 ng/ml
RD194442200CS	VEGF-D Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	31.25 – 2000 pg/ml
RZnT8/96R	Zinc Transporter 8 Autoantibody Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (citrate, heparin)	10–2000 units/ml
RD191229100R	Zymogen Granule Membrane Protein 16 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	1.5–30 ng/ml

PULMONARY DISEASE

PROTEINS

Cat. No.	Description	Species	Type	Source
RD172022100	Club Cell Protein (CC16)	Human	Recombinant	<i>E. coli</i>
RD172153100	Complement C1q Tumor Necrosis Factor-Related Protein 1	Human	Recombinant	<i>E. coli</i>
RD172203100	Interleukin-33	Human	Recombinant	<i>E. coli</i>
RD172208100	S100A1	Human	Recombinant	<i>E. coli</i>
RD172211100	S100A4	Human	Recombinant	<i>E. coli</i>
RD172219100	S100A10	Human	Recombinant	<i>E. coli</i>
RD172220100	S100A11	Human	Recombinant	<i>E. coli</i>
RD172221100	S100A12	Human	Recombinant	<i>E. coli</i>
RD172222100	S100A13	Human	Recombinant	<i>E. coli</i>
RD172225100	S100G	Human	Recombinant	<i>E. coli</i>
RD172373100	Secretoglobin Family 3A Member 2	Human	Recombinant	<i>E. coli</i>
RD172158100	Trefoil Factor 1	Human	Recombinant	<i>E. coli</i>
RD172159100	Trefoil Factor 2	Human	Recombinant	<i>E. coli</i>
RD172160100	Trefoil Factor 3	Human	Recombinant	<i>E. coli</i>
RD372160100	Trefoil Factor 3	Rat	Recombinant	<i>E. coli</i>

ANTIBODIES

Cat. No.	Description	Species	Type	Antigen
RD181022220-01	Club Cell Protein (CC16)	Human	Rabbit PAb	Human urine
RD181203100	Interleukin-33	Human	Rabbit PAb	<i>E. coli</i>
RD184203100	Interleukin-33	Human	Sheep PAb	<i>E. coli</i>
RD181158100	Trefoil Factor 1	Human	Rabbit PAb	<i>E. coli</i>
RD184158100	Trefoil Factor 1	Human	Sheep PAb	<i>E. coli</i>
RD181159100	Trefoil Factor 2	Human	Rabbit PAb	<i>E. coli</i>
RD184159100	Trefoil Factor 2	Human	Sheep PAb	<i>E. coli</i>
RD181160100	Trefoil Factor 3	Human	Rabbit PAb	<i>E. coli</i>
RD184160100	Trefoil Factor 3	Human	Sheep PAb	<i>E. coli</i>

› Contact Information



BioVendor – Laboratorní medicína a.s.
 Karasek 1767/1, 621 00 Brno, Czech Republic
 Phone: +420 549 124 185, Fax: +420 549 211 460
 E-mail: info@biovendor.com

BioVendor GmbH
 Otto-Hahn-Straße 16, 34123 Kassel, Germany
 Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
 E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
 Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
 Phone: +43 1 890 9025, Fax: +43 1 890 5163
 E-mail: infoAustria@biovendor.com

BioVendor, LLC
 128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
 Phone: +1-800-404-7807, Phone: +1-828-575-9250
 Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com