

RENAL DISEASE

IMMUNOASSAYS

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RSHAKRBS018R	Albumin Bovine ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant	0.78 - 50 ng/ml
RSHAKRAL121R	Albumin Mouse ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Cell culture medium, Urine, Plasma	50 - 1000 ng/ml
RSHAKRAL120R	Albumin Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant, Urine	50 - 1000 ng/ml
RA19001R	Albumin Rat ELISA High Sensitivity	Competitive ELISA	Urine	8 - 1000 ng/ml
RSHAKRAL020SR	Albumin Urinary (S-type) Rat Assay	Immunoturbidimetry/nephelometry	Serum, Urine, Plasma	6.17 - 500 µg/ml
RSHAKRAL021SR	Albumin Urinary Mouse (S-type) Assay	Immunoturbidimetry/nephelometry	Serum, Urine	6.17 - 500 µg/ml
RD192406200R	Aminoacylase-1 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (citrate, heparin), Urine	0.94 - 30 ng/ml
RA05880R	Angiotensin II Human/Rat ELISA	Sandwich ELISA, AChE-labelled antibody	Plasma (EDTA)	0.98 - 125 pg/ml
RD191328200R	Apoptosis Inhibitor of Macrophage Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Plasma (heparin), Serum	0.16 - 10 ng/ml
RCD003R	Beta2-Microglobulin Human ELISA	Competitive ELISA, Immobilized antigen	Serum	0.2 - 10 mg/l
RSCPRH111R	Beta2-Microglobulin Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Urine	0.31 - 10 ng/ml
RD491034200R	Clusterin Canine ELISA (ApoJ)	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	4 - 128 ng/ml
RD194034200R	Clusterin Human ELISA (ApoJ)	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine, CSF, Cell lysate	5 - 160 ng/ml
RD391034200CS	Clusterin Rat ELISA (ApoJ)	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	4 - 128 ng/ml
RD491009100R	Cystatin C Canine ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Urine	0.31 - 10 ng/ml
RD191009100	Cystatin C Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Urine, CSF	200 - 10000 ng/ml

RENAL DISEASE

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RD120009	Cystatin C Human PENTASET	ELISA development kit	Cell culture medium, ELISA, Others, Plasma, Serum	0.78 - 25 ng/ml
RD291009200R	Cystatin C Mouse ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.156 - 10 ng/ml
RD220009	Cystatin C Mouse PENTASET	ELISA development kit	Cell culture medium, ELISA, Others, Plasma, Serum	0.31 - 10 ng/ml
RD391009200R	Cystatin C Rat ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.078 - 25 ng/ml
RD20009	Cystatin C Rat PENTASET	ELISA development kit	Cell culture medium, ELISA, Others, Plasma, Serum	0.78 - 25 ng/ml
RAF013R	Erythropoietin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant	1.6 - 100 mIU/mL
RD191037100	Fetuin-A Human ELISA (AHSG)	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma	2 - 100 ng/ml
RD191407100R	Haptoglobin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Cerebrospinal fluid, Plasma (citrate, EDTA, heparin), Serum, Urine	3.13 - 200 ng/ml
RHBA001R	KIM-1 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Urine	0 - 10 ng/ml
RRBA001R	KIM-1 Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Urine	0 - 10 ng/ml
RD191102200R	Lipocalin 2/ NGAL Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA), Urine, Stool	0.3 - 10 ng/ml
RD191446200R	Osteopontin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Cerebrospinal fluid, Milk, Serum, Plasma, Urine	0.25 - 16 ng/ml
RD191046100R	Proguanylin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma	0.3 - 10 ng/ml
RMIO01R	Prorenin Human ELISA	Direct ELISA, HRP-labelled antibody	Plasma	0 - 10 ng/ml
RD191167200R	Total Prorenin/Renin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Amniotic fluid, Plasma (EDTA, citrate, heparin) Serum, Urine	25 - 1600 pg/ml
RD191113100R	Prostaglandin D Synthase (Lipocalin-type) Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Urine, CSF	1 - 40 ng/ml
RD191069200R	Prouroguanylin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate)	0.8 - 40 ng/ml
RD194169200R	Active Renin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Amniotic fluid, Plasma (EDTA, citrate, heparin), Serum, Urine	4.5 - 144 µIU/ml
RDS004R	Renin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma (EDTA)	0 - 270 pg/ml
RD191116200R	sRAGE Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	50 - 3200 pg/ml
RD191408200R	suPAR Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine	7.8 - 500 pg/ml
RD191160200R	Trefoil Factor 3 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine, Bronchoalveolar lavage	0.075 - 2.4 ng/ml
RD391160200R	Trefoil Factor 3 Rat ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	62.5 - 4000 pg/ml
RD191343100R	Urinary Trypsin Inhibitor Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Urine	0.63 - 20 ng/ml
RD491163200R	Uromodulin Canine ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	0.16 - 10 ng/ml
RD191163200R	Uromodulin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine	0.5 - 32 ng/ml

RENAL DISEASE

PROTEINS

Cat. No.	Description	Species	Type	Source
RD172112025	Adiponectin Globular	Human	Recombinant	<i>E. coli</i>
RD472034100	Clusterin (Apo J)	Canine	Recombinant	<i>E. coli</i>
RD472034100-HEK	Clusterin (Apo J)	Canine	Recombinant	HEK293
RD172034100	Clusterin (Apo J)	Human	Recombinant	HEK293
RD172034050-S	Clusterin (Apo J)	Human	Native	Plasma
RD372034100	Clusterin (Apo J)	Rat	Recombinant	<i>E. coli</i>
RD472009100	Cystatin C	Canine	Recombinant	<i>E. coli</i>
RD172009100	Cystatin C	Human	Native	Urine
RD172009100-H	Cystatin C	Human	Recombinant	<i>E. coli</i>
RD272009100	Cystatin C	Mouse	Recombinant	<i>E. coli</i>
RD372009100	Cystatin C	Rat	Recombinant	<i>E. coli</i>
RD172037100	Fetuin-A (AHSG)	Human	Native	Plasma
RD172241100	Fibroblast Growth Factor 23 C-terminal peptide	Human	Recombinant	<i>E. coli</i>
RD172150100	KIM-1	Human	Recombinant	HEK293
RD172099100	Klotho	Human	Recombinant	CHO
RD172370100	Leukocyte Cell-Derived Chemotaxin-2	Human	Recombinant	<i>E. coli</i>
RD172102050	Lipocalin-2/NGAL	Human	Recombinant	<i>E. coli</i>
RD172240100	Liver FABP	Human	Recombinant	<i>E. coli</i>
RD172240101	Liver FABP, Tagless	Human	Recombinant	<i>E. coli</i>
RD172170100	Osteoactivin/GPNMB	Human	Recombinant	HEK293
RD172446100	Osteopontin	Human	Recombinant	HEK293
RD172330100	Podocin	Human	Recombinant	<i>E. coli</i>
RD172046100	Proguanylin	Human	Recombinant	<i>E. coli</i>
RD172113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Recombinant	<i>E. coli</i>
RD172069100	Prouroguanylin	Human	Recombinant	<i>E. coli</i>
RD172116100	sRAGE	Human	Recombinant	<i>E. coli</i>
RD172116100-HEK	sRAGE	Human	Recombinant	HEK293
RD172167100-HEK	Renin	Human	Recombinant	HEK293
RD172264100	Renin Receptor	Human	Recombinant	<i>E. coli</i>
RD172116025-HEK	Soluble Advanced Glycation End Product-Specific Receptor	Human	Recombinant	HEK293
RD172160100	Trefoil Factor 3	Human	Recombinant	<i>E. coli</i>
RD372160100	Trefoil Factor 3	Rat	Recombinant	<i>E. coli</i>
RD463163100	Uromodulin	Canine	Native	Canine urine
RD663163050	Uromodulin	Feline	Native	Feline urine
RD172163100	Uromodulin	Human	Native	Human urine
RD563163100	Uromodulin	Porcine	Native	Porcine urine

RENAL DISEASE

ANTIBODIES

Cat. No.	Description	Spec. React.	Type	Antigen
RD184328100	Apoptosis Inhibitor of Macrophage	Human	Sheep PAb	<i>E. coli</i>
RD481034100	Clusterin (Apo J)	Canine	Rabbit PAb	<i>E. coli</i>
RD484034100	Clusterin (Apo J)	Canine	Sheep PAb	<i>E. coli</i>
RD182034110-H3	Clusterin (Apo J)	Human	Mouse MAb, Clone Hs-3	Human sperm extract
RD381034100	Clusterin (Apo J)	Rat	Rabbit PAb	<i>E. coli</i>
RD481009100	Cystatin C	Canine	Rabbit PAb	<i>E. coli</i>
RD484009100	Cystatin C	Canine	Sheep PAb	<i>E. coli</i>
RD181009220-01	Cystatin C	Human	Rabbit PAb	Human urine
RD181009100	Cystatin C	Human	Rabbit PAb	<i>E. coli</i>
RD184009100	Cystatin C	Human	Sheep PAb	<i>E. coli</i>
RD184009220	Cystatin C	Human	Sheep PAb	Human urine
RD281009100	Cystatin C	Mouse	Rabbit PAb	<i>E. coli</i>
RD284009100	Cystatin C	Mouse	Sheep PAb	<i>E. coli</i>
RD384009100	Cystatin C	Rat	Sheep PAb	<i>E. coli</i>
RD184036100	Fetuin-A (AHSG)	Human	Goat PAb	Human serum
RD181036100	Fetuin-A (AHSG)	Human	Rabbit PAb	Human serum
RD181109100	Fibroblast Growth Factor 23	Human	Rabbit PAb	<i>E. coli</i>
RD181170100	Osteoactivin/GPNMB	Human	Rabbit PAb	HEK293
RD184170100	Osteoactivin/GPNMB	Human	Sheep PAb	HEK293
RD181113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Rabbit PAb	<i>E. coli</i>
RD184113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Sheep PAb	<i>E. coli</i>
RD184116100	sRAGE	Human	Sheep PAb	<i>E. coli</i>
RD181221100	S100A12	Human	Rabbit PAb	<i>E. coli</i>
RD184221100	S100A12	Human	Sheep PAb	<i>E. coli</i>
RD181160100	Trefoil Factor 3	Human	Rabbit PAb	<i>E. coli</i>
RD184160100	Trefoil Factor 3	Human	Sheep PAb	<i>E. coli</i>
RD181343100	Urinary Trypsin Inhibitor	Human	Rabbit PAb	Human Urine
RD181163100	Uromodulin	Human	Rabbit PAb	<i>E. coli</i>
RD184163100	Uromodulin	Human	Sheep PAb	<i>E. coli</i>

› Contact Information



BioVendor – Laboratorní medicína a.s.
 Karasek 1767/1, 621 00 Brno, Czech Republic
 Phone: +420 549 124 185, Fax: +420 549 211 460
 E-mail: info@biovendor.com

BioVendor GmbH
 Otto-Hahn-Straße 16, 34123 Kassel, Germany
 Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
 E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
 Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
 Phone: +43 1 890 9025, Fax: +43 1 890 5163
 E-mail: infoAustria@biovendor.com

BioVendor, LLC
 128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
 Phone: +1-800-404-7807, Phone: +1-828-575-9250
 Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com