

# STEROID HORMONES

## IMMUNOASSAYS

Cat. No.	Name	Sample	Calibration range	Limit of detection
<b>RCD007R</b>	<b>5<math>\alpha</math>-Androstane-3<math>\alpha</math>,17<math>\beta</math>-diol Glucuronide ELISA</b>	Serum	0.25 – 50 ng/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 0.1 ng/ml
<b>RCD021R</b>	<b>17<math>\alpha</math>-Hydroxyprogesterone ELISA</b>	Serum	0.15 – 20 ng/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 0.11 ng/ml
<b>RTC010R</b>	<b>17-OH-Progesterone free in saliva ELISA</b>	Saliva	10 – 1000 pg/ml	
	Competitive ELISA, Immobilized antigen			25 $\mu$ l 2.5 pg/ml
<b>RCD002R</b>	<b>Aldosterone ELISA</b>	Serum, Plasma, Urine	15 – 1000 pg/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 10 pg/ml
<b>RCD001R</b>	<b>Androstenedione ELISA</b>	Serum	0.1 – 10 ng/ml	
	Competitive ELISA, Immobilized antigen			25 $\mu$ l 0.05 ng/ml
<b>RTC011R</b>	<b>Androstenedione ELISA</b>	Saliva	20 – 1000 pg/ml	
	Competitive ELISA			50 $\mu$ l 5 pg/ml
<b>RTC011R</b>	<b>Androstenedione ELISA</b>	Saliva	20 – 1000 pg/ml	
	Competitive ELISA			50 $\mu$ l 5 pg/ml
<b>RTC011R</b>	<b>Androstenedione free in saliva ELISA</b>	Saliva	20 – 1000 pg/ml	
	Competitive ELISA			50 $\mu$ l 5 pg/ml
<b>RD192234200R</b>	<b>Corticosteroid-Binding Globulin Human ELISA</b>	Serum, Plasma (EDTA, citrate, heparin)	3.13 – 200 ng/ml	
	Sandwich ELISA, Biotin-labelled antibody			5 $\mu$ l 0.01 ng/ml
<b>RSCYK240R</b>	<b>Corticosterone Human ELISA (Multispecies specificity)</b>	Serum, Urine, Plasma, Cell culture supernatant	0.21 – 50 ng/ml	
	Competitive ELISA, Immobilized antibody			
<b>RTC002R</b>	<b>Corticosterone Mouse/Rat ELISA</b>	Serum, Plasma	15 – 2250 ng/ml	
	Competitive ELISA, Immobilized antigen			10 $\mu$ l 4.1 ng/ml
<b>RCD004R</b>	<b>Cortisol ELISA</b>	Serum	0.5-60 $\mu$ g/dl	
	Competitive ELISA, Immobilized antigen			20 $\mu$ l 0.4 $\mu$ g/dl
<b>RCD005R</b>	<b>Cortisol (Saliva) Human ELISA</b>	Saliva	1 – 100 ng/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 1.0 ng/ml
<b>RCD008R</b>	<b>Dehydroepiandrosterone Sulfate ELISA</b>	Serum	0.005 – 10 $\mu$ g/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 0.005 $\mu$ g/ml
<b>RTC012R</b>	<b>DHEA free in saliva ELISA</b>	Saliva	10 – 1440 pg/ml	
	Competitive ELISA			50 $\mu$ l 2.2 pg/ml
<b>RCD009R</b>	<b>Dihydrotestosterone ELISA</b>	Serum	25 – 2500 pg/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 6 pg/ml
<b>RCD011R</b>	<b>Estradiol ELISA</b>	Serum	20 – 3200 pg/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 10 pg/ml
<b>RTC013R</b>	<b>Estradiol free in saliva ELISA</b>	Saliva	1 – 100 pg/ml	
	Competitive ELISA			100 $\mu$ l 0.4 pg/ml
<b>RTC014R</b>	<b>Estriol free in saliva ELISA</b>	Saliva	2.5 – 4000 pg/ml	
	Competitive ELISA			50 $\mu$ l 1.4 pg/ml
<b>RTC009R</b>	<b>Estradiol Rat ELISA</b>	Plasma	2.5 – 1280 pg/ml	
	Competitive ELISA, Immobilized antigen			
<b>RTC014R</b>	<b>Estriol free in saliva ELISA</b>	Saliva	2.5 – 4000 pg/ml	
	Competitive ELISA			50 $\mu$ l 1.4 pg/ml
<b>RCD010R</b>	<b>Estrone ELISA</b>	Serum	15 – 2000 pg/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 10 pg/ml
<b>RTC015R</b>	<b>Estrone free in saliva ELISA</b>	Saliva	10 – 300 pg/ml	
	Competitive ELISA			100 $\mu$ l 4.7 pg/ml
<b>RIS009R</b>	<b>Free Estriol Human ELISA</b>	Serum	0.3 – 40 ng/ml	
	Sandwich ELISA, HRP-labelled antibody			10 $\mu$ l 0.075 ng/ml
<b>RCD022R</b>	<b>Pregnenolone ELISA</b>	Serum	0.1 – 25 ng/ml	
	Competitive ELISA, Immobilized antigen			50 $\mu$ l 0.05 ng/ml

## STEROID HORMONES

Cat. No.	Name	Sample	Calibration range	Limit of detection
<b>RCD020R</b>	<b>Progesterone ELISA</b> Competitive ELISA, Immobilized antigen	Serum	0.3 – 60 ng/ml	0.1 ng/ml
<b>RTC016R</b>	<b>Progesterone free in saliva ELISA</b> Competitive ELISA	Saliva	10 – 5000 pg/ml	5 pg/ml
<b>RTC008R</b>	<b>Progesterone Mouse/Rat ELISA</b> Sandwich ELISA, HRP-labelled antibody	Serum, Plasma	0.4 – 100 ng/ml	0.04 ng/ml
<b>RTC900R</b>	<b>Rat Quality Control Set (Corticosterone, Prolactin, Testosterone, TSH)</b>			
<b>RCD027R</b>	<b>Testosterone ELISA</b> Competitive ELISA, Immobilized antigen	Serum	0.08-16 ng/ml	0.022 ng/ml
<b>RTC017R</b>	<b>Testosterone free in saliva ELISA</b> Competitive ELISA	Saliva	10 – 1000 pg/ml	2.2 pg/ml
<b>RCD024R</b>	<b>Sex Hormone Binding Globulin Human ELISA</b> Sandwich ELISA, HRP-labelled antibody	Serum	3.3 – 295 nmol/l	0.1 nmol/l
<b>RCD016R</b>	<b>Testosterone, Free Human ELISA</b> Competitive ELISA, Immobilized antigen	Serum	0.25 – 125 pg/ml	0.17 pg/ml
<b>RTC001R</b>	<b>Testosterone Mouse/Rat ELISA</b> Competitive ELISA, Immobilized antigen	Serum, Plasma	0.1 – 25 ng/ml	0.066 ng/ml

## PROTEINS

Cat. No.	Description	Species	Type	Source
RD272464100	Corticosteroid-Binding Globulin	Mouse	Recombinant	HEK293
RD172151100	Corticotropin-Releasing Factor-Binding Protein	Human	Recombinant	<i>E. coli</i>
RD172099100	Klotho Human CHO	Human	Recombinant	CHO
RGB001	Sex Hormone-Binding Globulin	Human	Recombinant	CHO

## ANTIBODIES

Cat. No.	Description	Spec. React.	Type	Antigen
RD184151100	Corticotropin-Releasing Factor-Binding Protein	Human	Sheep PAb	<i>E. coli</i>

### › Contact Information



**BioVendor – Laboratorní medicína a.s.**  
 Karasek 1767/1, 621 00 Brno, Czech Republic  
 Phone: +420 549 124 185, Fax: +420 549 211 460  
 E-mail: info@biovendor.com

**BioVendor GmbH**  
 Otto-Hahn-Straße 16, 34123 Kassel, Germany  
 Phone: +49 6221 4339 100, Fax: +49 6221 4339 111  
 E-mail: infoEU@biovendor.com

› [www.biovendor.com](http://www.biovendor.com)

**BioVendor GesmbH**  
 Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria  
 Phone: +43 1 890 9025, Fax: +43 1 890 5163  
 E-mail: infoAustria@biovendor.com

**BioVendor, LLC**  
 128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA  
 Phone: +1-800-404-7807, Phone: +1-828-575-9250  
 Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com