

CARDIOVASCULAR DISEASE

IMMUNOASSAYS

Cat. No.	Description	Assay Format	Sample Type	Assay Range
REA201/96	ADMA ELISA	Competitive ELISA, Immobilized antigen	Serum, Plasma (EDTA)	0.1 – 5 µmol/l
REA209/96R	ADMA ELISA, High Sensitivity	Competitive ELISA, Immobilized antigen	Serum, Plasma (EDTA, heparin), Cell culture medium	
RD191204200R	Allograft Inflammatory Factor 1 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, Citrate)	31.3 – 1000 pg/ml
RSHAKHB48R	Apo B-48 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (heparin, EDTA)	2.5 – 160 ng/ml
RD191236100R	Apolipoprotein H/Beta2-GP1 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Urine	1.25 – 80 ng/ml
RD191328200R	Apoptosis Inhibitor of Macrophage Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (heparin)	0.16 – 10 ng/ml
RA19002R	Atriopeptin Rat ELISA		Plasma, Cardia Tissue	16 – 2000 pg/ml
RD192026200CS	Cardiotrophin-1 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma (EDTA)	5 pg/ml – 1 ng/ml
RAF006R	Caspase-8 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Cell culture supernatant, Cell lysate	0.16 – 10 ng/ml
RA05482.96R	CGRP Rat ELISA kit	Sandwich ELISA, AChE-labelled antibody	Serum, Cerebrospinal fluid, Plasma, Cell culture supernatant, Nervous tissue	4 – 500 pg/ml
RD191394200R	Collagen Triple Helix Repeat-Containing Protein 1	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Tissue extract	78 – 2500 pg/ml
RA019482R	CRGP Affinity Sorbent			
RD191180200R	CTRP9 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	25 – 800 pg/ml

CARDIOVASCULAR DISEASE

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RD191329200R	Cyclophilin A Human ELISA	Direct ELISA, Biotin-labelled antibody	Serum, Plasma	0.39 - 25 ng/ml
RD193324200R	Decorin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	0.125 - 4 ng/ml
RD191060200R	Epidermal FABP Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin), Urine, Milk, Tissue extract, Tissue culture medium	0.5 - 40 ng/ml
RD191037100	Fetuin-A (AHSG) Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	2 - 100 ng/ml
RD191135200R	GDF-15/MIC-1 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	70 - 4480 pg/ml
RD191230200R	Glucose-Regulated Protein 78 Human ELISA	Direct ELISA, Biotin-labelled antibody	Serum, Plasma	0.78 - 50 ng/ml
RD191407100R	Haptoglobin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Cerebrospinal fluid, Plasma (citrate, EDTA, heparin), Serum, Urine	3.13 - 200 ng/ml
REA205/96R	Homoarginine ELISA	Competitive ELISA	Serum, Plasma (EDTA, citrate, heparin)	31.25 - 2000 pg/ml
RAF064R	Interleukin-33 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (citrate, heparin), Cell culture supernatant	4.0 - 250 pg/ml
RD191122200R	Lecithin-Cholesterol Acyltransferase Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate), CSF	2.5 - 160 ng/ml
RHK324R	MPO Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Plasma, Cell culture medium, Sputum	0.4 - 100 ng/ml
RSHAKDNPA-012R	NT-pro ANP Dog ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA)	12.5 - 400 pg/ml
RD191279200R	Paraoxonase 1 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (citrate)	1.25 - 80 ng/ml
RSCPP-PD03-E0R	Pentraxin 3 Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma	0.1 - 20 ng/ml
RAG019R	Periostin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell culture supernatant	78 - 5000 pg/ml
RAG020R	Periostin Mouse ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell culture supernatant	31 - 2000 pg/ml
RDS004R	Renin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma (EDTA)	0 - 270 pg/ml
RD191220100R	S100A11 Human ELISA (Multispecies specificity)	Sandwich ELISA, HRP-labelled antibody	Serum, Bronchoalveolar lavage, Amniotic fluid	0.313 - 100 ng/ml
RD191225200R	S100G Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	25 - 800 pg/ml
RD191048200R	S100P Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate)	0.0625 - 2.0 ng/ml
RAF093R	sCD40L Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell culture supernatant (only high levels)	0.16 - 10 ng/ml
RAF094R	sCD40L Mouse ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Cell culture supernatant	0.31 - 20 ng/ml
RLF-EK0128R	Troponin I Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Cell lysate	0 - 25 ng/ml
RD194442200CS	VEGF-D Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin)	31.25 - 2000 pg/ml

CARDIOVASCULAR DISEASE

PROTEINS

Cat. No.	Description	Species	Type	Source
RD172410100	6-phosphogluconate dehydrogenase, decarboxylating	Human	Recombinant	<i>E. coli</i>
RD172204100	Allograft Inflammatory Factor 1	Human	Recombinant	<i>E. coli</i>
RD172406100	Aminoacylase-1	Human	Recombinant	<i>E. coli</i>
RD172026100	Cardiotrophin-1	Human	Recombinant	<i>E. coli</i>
RD172180100	CTRP9A	Human	Recombinant	<i>E. coli</i>
RD172329100	Cyclophilin A	Human	Recombinant	<i>E. coli</i>
RD272455100	Cyclophilin A	Mouse	Recombinant	<i>E. coli</i>
RD372454100	Cyclophilin A	Rat	Recombinant	<i>E. coli</i>
RD172173100	Endothelial Lipase	Human	Recombinant	HEK293
RD172060100	Epidermal Fatty Acid Binding Protein	Human	Recombinant	<i>E. coli</i>
RD172037100	Fetuin-A (AHSG) NATIVE	Human	Native	Human plasma
RD172393100	FGF-18	Human	Recombinant	<i>E. coli</i>
RD172203100	Interleukin-33	Human	Recombinant	<i>E. coli</i>
RD172122100	Lecithin-Cholesterol Acyltransferase	Human	Recombinant	<i>E. coli</i>
RD172122100-HEK	Lecithin-Cholesterol Acyltransferase	Human	Recombinant	HEK293
RD172045100	Osteoblast Specific Factor 2	Human	Recombinant	<i>E. coli</i>
RD172045025-HEK	Osteoblast Specific Factor 2	Human	Recombinant	HEK293
RD172045025	Osteoblast Specific Factor 2	Human	Recombinant	<i>E. coli</i>
RD172279100	Paraoxonase 1	Human	Recombinant	<i>E. coli</i>
RD172316100	Peptidoglycan Recognition Protein 1	Human	Recombinant	<i>E. coli</i>
RD172167100-HEK	Renin	Human	Recombinant	HEK293
RD172264100	Renin Receptor	Human	Recombinant	<i>E. coli</i>
RD172208025	S100A1	Human	Recombinant	<i>E. coli</i>
RD172208100	S100A1	Human	Recombinant	<i>E. coli</i>
RD172211100	S100A4	Human	Recombinant	<i>E. coli</i>
RD172219100	S100A10	Human	Recombinant	<i>E. coli</i>
RD172220100	S100A11	Human	Recombinant	<i>E. coli</i>
RD172222100	S100A13	Human	Recombinant	<i>E. coli</i>
RD172225100	S100G	Human	Recombinant	<i>E. coli</i>
RD172416100	Ubiquitin Carboxyl-Terminal Hydrolase Isozyme L1	Human	Recombinant	<i>E. coli</i>

CARDIOVASCULAR DISEASE

ANTIBODIES

Cat. No.	Description	Spec. React.	Type	Antigen
RD184204100	Allograft Inflammatory Factor 1	Human	Sheep PAb	<i>E. coli</i>
RD184328100	Apoptosis Inhibitor of Macrophage	Human	Sheep PAb	<i>E. coli</i>
RD1820261001F3	Cardiotrophin-1	Human	Mouse MAb, Clone 1F3	<i>E. coli</i>
RD182026100-01	Cardiotrophin-1	Human	Mouse MAb, Clone 3G6D9	<i>E. coli</i>
RD181026100	Cardiotrophin-1	Human	Rabbit PAb	<i>E. coli</i>
RD181060100	Epidermal Fatty Acid Binding Protein	Human	Rabbit PAb	<i>E. coli</i>
RD184060100	Epidermal Fatty Acid Binding Protein	Human	Sheep PAb	<i>E. coli</i>
RD184036100	Fetuin-A (AHSG)	Human	Goat PAb	Human serum
RD181036100	Fetuin-A (AHSG)	Human	Rabbit PAb	Human serum
RD181203100	Interleukin-33	Human	Rabbit PAb	<i>E. coli</i>
RD184203100	Interleukin-33	Human	Sheep PAb	<i>E. coli</i>
RD181057050	Myosin (Cardiac) Light Chain	Human	Rabbit PAb	<i>E. coli</i>
RD181045050	Osteoblast Specific Factor 2	Human	Rabbit PAb	<i>E. coli</i>
RD184045100	Osteoblast Specific Factor 2	Human	Sheep PAb	<i>E. coli</i>
RD184219100	S100A10	Human	Sheep PAb	<i>E. coli</i>
RD184220100	S100A11	Human	Sheep PAb	<i>E. coli</i>
RD184222100	S100A13	Human	Sheep PAb	<i>E. coli</i>
RD184225100	S100G	Human	Sheep PAb	<i>E. coli</i>

› Contact Information



BioVendor – Laboratorní medicína a.s.
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185, Fax: +420 549 211 460
E-mail: info@biovendor.com

BioVendor GmbH
Otto-Hahn-Straße 16, 34123 Kassel, Germany
Phone: +49 6221 4339 100, Fax: +49 6221 4339 111
E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
Nußdorfer Straße 20/10, 1090 Vienna, Austria
Phone: +43 1 890 9025, Fax: +43 1 890 9025-15
E-mail: infoAustria@biovendor.com

BioVendor, LLC
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
Phone: +1-800-404-7807, Phone: +1-828-575-9250
Fax: +1-828-575-9251, E-mail: infoUSA@biovendor.com



BioMedical Specimen Bank
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185, Fax: +420 549 211 460
E-mail: info@biomedicalspecimenbank.com

› www.biomedicalspecimenbank.com